

## Drying Racks



188050000

189330023



HS23243A

### Drying Racks, Benchtop, Bel-Art

Single-sided. **188050000** .....\$170.00

Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 43 x 17 x 48 cm(H).

Double-sided

2-Tier. **189330022** .....\$391.00

3-Tier. **189330023** .....\$539.00

**NEW**

### Drying Rack, Benchtop/Wallmount, Heathrow

Single-sided. **HS23243A** .....\$146.00

Dry up to 90 pieces of labware with this self-standing or wall-mounted rack. Constructed of vinyl-coated wire hairpin loops with a polystyrene tray and rack. Dimensions: 46.2 x 18.2 x 52.5 cm(H).

## Electrophoresis



Electrophoresis System, MyGel™ Mini. E1101 .....\$607.00

## Eyewear

Absorbs  
99.9% UV



Glasses, Safety, integral side shield. **43574002** .....\$6.00



Flexible frame and adjustable elastic strap provides a tight seal.

Goggles, Safety, vented (chemical splash). **97-2001** .....\$4.00

## Filter & Blotting Papers



Ahlstrom - Munksjo  
Filter Paper

### Qualitative, Standard Grades (for gravity filtration)

#### Slow flow (4 mL/min), 1.5 µm retention, Grade 610

(Whatman 5; S&S 602)

9 cm, Pk/100. **6100-0900** .....\$12.20

12.5 cm, Pk/100. **6100-1250** .....\$15.25

15 cm, Pk/100. **6100-1500** .....\$18.40

20.5 cm, Pk/100. **6100-2050** .....\$29.05

24 cm, Pk/100. **6100-2400** .....\$40.45

#### Medium flow (30 mL/min), 2 µm retention, Grade 642

(Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)

9 cm, Pk/100. **6420-0900** .....\$9.50

15 cm, Pk/100. **6420-1500** .....\$17.05

#### Medium flow (55 mL/min), 2.5 µm retention, Grade 601

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

4.25 cm, Pk/100. **6010-0425** .....\$7.25

9 cm, Pk/100. **6010-0900** .....\$9.80

12.5 cm, Pk/100. **6010-1250** .....\$14.10

15 cm, Pk/100. **6010-1500** .....\$17.15

18.5 cm, Pk/100. **6010-1850** .....\$24.85

20.5 cm, Pk/100. **6010-2050** .....\$26.80

33 cm, Pk/100. **6010-3300** .....\$72.40



#### Medium flow (60 mL/min), 6 µm retention, Grade 613

(Whatman Student; S&S 595, 0903)

9 cm, Pk/100. **6130-0900** .....\$8.35

15 cm, Pk/100. **6130-1500** .....\$14.50

20.5 cm, Pk/100. **6130-2050** .....\$21.80

33 cm, Pk/100. **6130-3300** .....\$48.45

#### Medium flow (200 mL/min), 10 µm retention, Grade 631

(Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)

20.5 cm, Pk/100. **6310-2050** .....\$27.00

33 cm, Pk/100. **6310-3300** .....\$73.15

#### Fast flow (235 mL/min), 25 µm retention, creped, Grade 615

9 cm, Pk/100. **6150-0900** .....\$8.60

12.5 cm, Pk/100. **6150-1250** .....\$12.50

15 cm, Pk/100. **6150-1500** .....\$15.00

20.5 cm, Pk/100. **6150-2050** .....\$22.50

33 cm, Pk/100. **6150-3300** .....\$49.90

#### Very Fast flow (360 mL/min), 35 µm retention, creped, Grade 617

9 cm, Pk/50. **6170-0900** .....\$7.85

12.5 cm, Pk/50. **6170-1250** .....\$10.45

15 cm, Pk/50. **6170-1500** .....\$12.15

20.5 cm, Pk/50. **6170-2050** .....\$18.55

33 cm, Pk/50. **6170-3300** .....\$39.80

### Qualitative, Pleated (Fluted)

#### Slow flow (30 mL/min), 2 µm retention, Grade 542

(Whatman 2V)

12.5 cm, Pk/100. **5420-1250** .....\$78.15

18.5 cm, Pk/100. **5420-1850** .....\$88.45

Continued

**Qualitative, Wet-Strengthened Grades** (for vacuum filtration)

- Medium flow** (65 mL/min), 3 µm retention, **Grade 909**  
 (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)  
 32 cm, Pk/100. **9090-3200** .....\$78.85
- Fast flow** (245 mL/min), 27 µm retention, creped, **Grade 962**  
 (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)  
 11 cm, Pk/100. **9620-1100** .....\$11.70

**Quantitative, Ashless**

(Designed for use in analytical and gravimetric analysis)

- Slow flow** (4 mL/min), 1.5 µm retention, **Grade 94**  
 (Whatman 42; S&S 589/3)  
 9 cm, Pk/100. **0940-0900** .....\$33.00  
 15 cm, Pk/100. **0940-1500** .....\$56.85
- Medium flow** (30 mL/min), 2 µm retention, **Grade 74**  
 (Whatman 40; S&S 589/2, 589/4)  
 9 cm, Pk/100. **0740-0900** .....\$33.00  
 15 cm, Pk/100. **0740-1500** .....\$56.85
- Fast flow** (220 mL/min), 15 µm retention, **Grade 55**  
 (Whatman 54, 541; S&S 1505)  
 9 cm, Pk/100. **0550-0900** .....\$37.90  
 15 cm, Pk/100. **0550-1500** .....\$59.95

**Blotting/Chromatography/Electrophoresis Papers**

- 0.19 mm thick**, 125 mm/30 min absorption, **Grade 601**  
 (Whatman 1, 1CHR; S&S 593, 597, 2043A)  
 46 x 57 cm, Pk/100. **6018-4657** .....\$147.00
- 0.34 mm thick**, 140 mm/30 min absorption, **Grade 238**  
 (Whatman 3MM, S&S 598,231<sup>o</sup>)  
 20 x 20 cm, Pk/100. **2388-2020** .....\$75.25  
 46 x 57 cm, Pk/100. **2388-4657** .....\$286.00
- 0.42 mm thick**, 150 mm/30 min absorption, **Grade 237**  
 (Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703;  
 Albet 2040B, FP591; Munktell 3S/H, FN7a)  
 46 x 57 cm, Pk/100. **2378-4657** .....\$329.00
- 0.95 mm thick**, 180 mm/30 min absorption, **Grade 243**  
 (Whatman 17CHR; S&S 2668)  
 11 x 14 cm, Pk/100. **2438-1114** .....\$31.75  
 20 x 20 cm, Pk/100. **2438-2020** .....\$62.70
- 2.48 mm thick**, 240 mm/30 min absorption, **Grade 320**  
 (S&S 2727)  
 11 x 14 cm, Pk/50. **3208-1114** .....\$60.10  
 20 x 20 cm, Pk/50. **3208-2020** .....\$87.45

**NEW**



Smooth tips & polypropylene grips

- Forceps, Filter, 11.5 cm. 300-025** .....\$4.40

**Freezer Racks**



**Freezer Racks, stainless steel, USA**

175V-13-2L

Restocking fee of 20% on all returned freezer racks listed below.  
 No returns on custom racks.

**Horizontal** (fits upright freezers)

- Horizontal, Holds 2" Boxes**
- |   |            |            |
|---|------------|------------|
|   | <u>1-9</u> | <u>≥10</u> |
| 9-Place (3L x 3H). <b>175HS-9-2</b> .....   | \$59.00    | 53.00      |
| Dimensions: 16.31(L) x 5.5(W) x 6.5"(H)     |            |            |
| 12-Place (3L x 4H). <b>175HS-12-2</b> ..... | \$78.00    | 71.00      |
| Dimensions: 16.31(L) x 5.5(W) x 8.75"(H)    |            |            |
| 15-Place (3L x 5H). <b>175HS-15-2</b> ..... | \$98.00    | 89.00      |
| Dimensions: 16.31(L) x 5.5(W) x 11"(H)      |            |            |
| 16-Place (4L x 4H). <b>175HL-16-2</b> ..... | \$105.00   | 95.00      |
| Dimensions: 22(L) x 5.5(W) x 8.75"(H)       |            |            |
| 20-Place (4L x 5H). <b>175HL-20-2</b> ..... | \$132.00   | 119.00     |
| Dimensions: 22(L) x 5.5(W) x 11"(H)         |            |            |
- Horizontal, Holds 3" Boxes**
- |   |         |       |
|---|---------|-------|
| 9-Place (3L x 3H). <b>175HS-9-3</b> .....   | \$61.00 | 55.00 |
| Dimensions: 16.31(L) x 5.5(W) x 9.5"(H)     |         |       |
| 12-Place (4L x 3H). <b>175HL-12-3</b> ..... | \$80.00 | 73.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H)        |         |       |



**Horizontal Drawer** (fits upright freezers)

- Horizontal Drawer, Holds 2" Boxes**
- |  |            |            |
|--|------------|------------|
|  | <u>1-9</u> | <u>≥10</u> |
| 9-Place (3L x 3H). <b>175HD-332</b> .....  | \$85.00    | 77.00      |
| Dimensions: 16.5(L) x 5.5(W) x 7.125"(H)   |            |            |
| 12-Place (3L x 4H). <b>175HD-342</b> ..... | \$111.00   | 101.00     |
| Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)     |            |            |
| 15-Place (3L x 5H). <b>175HD-352</b> ..... | \$141.00   | 128.00     |
| Dimensions: 16.5(L) x 5.5(W) x 11.812"(H)  |            |            |
| 16-Place (4L x 4H). <b>175HD-442</b> ..... | \$151.00   | 137.00     |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H)       |            |            |
| 20-Place (5L x 4H). <b>175HD-542</b> ..... | \$188.00   | 171.00     |
| Dimensions: 27.25(L) x 5.5(W) x 9.5"(H)    |            |            |
- Horizontal Drawer, Holds 5" Boxes**
- |   |            |            |
|---|------------|------------|
|   | <u>1-9</u> | <u>≥10</u> |
| 6-Place (3L x 2H). <b>175RD-325</b> ..... | \$74.00    | 67.00      |
| Dimensions: 18.125(L) x 6(W) x 10.75"(H)  |            |            |
| 8-Place (4L x 2H). <b>175RD-425</b> ..... | \$99.00    | 90.00      |
| Dimensions: 24(L) x 6(W) x 10.75"(H)      |            |            |
- Horizontal Drawer, Holds 15 mL tubes**
- |   |          |        |
|---|----------|--------|
| 160-Place (80L x 2H). <b>175HD-32ML-15</b> .....  | \$99.00  | 90.00  |
| Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)            |          |        |
| 220-Place (110L x 2H). <b>175HD-42ML-15</b> ..... | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H)              |          |        |
- Horizontal Drawer, Holds 50 mL tubes**
- |   |          |        |
|---|----------|--------|
| 84-Place (42L x 2H). <b>175HD-42ML-50</b> ..... | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H)            |          |        |

Continued