

Qualitative, Wet-Strengthened Grades (for vacuum filtration)

- Medium flow** (65 mL/min), 3 µm retention, **Grade 909**
 (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)
 32 cm, Pk/100. **9090-3200**\$78.85
- Fast flow** (245 mL/min), 27 µm retention, creped, **Grade 962**
 (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)
 11 cm, Pk/100. **9620-1100**\$11.70

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

- Slow flow** (4 mL/min), 1.5 µm retention, **Grade 94**
 (Whatman 42; S&S 589/3)
 9 cm, Pk/100. **0940-0900**\$33.00
 15 cm, Pk/100. **0940-1500**\$56.85
- Medium flow** (30 mL/min), 2 µm retention, **Grade 74**
 (Whatman 40; S&S 589/2, 589/4)
 9 cm, Pk/100. **0740-0900**\$33.00
 15 cm, Pk/100. **0740-1500**\$56.85
- Fast flow** (220 mL/min), 15 µm retention, **Grade 55**
 (Whatman 54, 541; S&S 1505)
 9 cm, Pk/100. **0550-0900**\$37.90
 15 cm, Pk/100. **0550-1500**\$59.95

Blotting/Chromatography/Electrophoresis Papers

- 0.19 mm thick**, 125 mm/30 min absorption, **Grade 601**
 (Whatman 1, 1CHR; S&S 593, 597, 2043A)
 46 x 57 cm, Pk/100. **6018-4657**\$147.00
- 0.34 mm thick**, 140 mm/30 min absorption, **Grade 238**
 (Whatman 3MM, S&S 598,231)
 20 x 20 cm, Pk/100. **2388-2020**\$75.25
 46 x 57 cm, Pk/100. **2388-4657**\$286.00
- 0.42 mm thick**, 150 mm/30 min absorption, **Grade 237**
 (Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703;
 Albet 2040B, FP591; Munktell 3S/H, FN7a)
 46 x 57 cm, Pk/100. **2378-4657**\$329.00
- 0.95 mm thick**, 180 mm/30 min absorption, **Grade 243**
 (Whatman 17CHR; S&S 2668)
 11 x 14 cm, Pk/100. **2438-1114**\$31.75
 20 x 20 cm, Pk/100. **2438-2020**\$62.70
- 2.48 mm thick**, 240 mm/30 min absorption, **Grade 320**
 (S&S 2727)
 11 x 14 cm, Pk/50. **3208-1114**\$60.10
 20 x 20 cm, Pk/50. **3208-2020**\$87.45

NEW



Smooth tips & polypropylene grips

- Forceps, Filter, 11.5 cm. 300-025**\$4.40

Freezer Racks



Freezer Racks, stainless steel, USA

175V-13-2L

Restocking fee of 20% on all returned freezer racks listed below.
 No returns on custom racks.

Horizontal (fits upright freezers)

- Horizontal, Holds 2" Boxes**
- | | | |
|---|------------|------------|
| | <u>1-9</u> | <u>≥10</u> |
| 9-Place (3L x 3H). 175HS-9-2 | \$59.00 | 53.00 |
| Dimensions: 16.31(L) x 5.5(W) x 6.5"(H) | | |
| 12-Place (3L x 4H). 175HS-12-2 | \$78.00 | 71.00 |
| Dimensions: 16.31(L) x 5.5(W) x 8.75"(H) | | |
| 15-Place (3L x 5H). 175HS-15-2 | \$98.00 | 89.00 |
| Dimensions: 16.31(L) x 5.5(W) x 11"(H) | | |
| 16-Place (4L x 4H). 175HL-16-2 | \$105.00 | 95.00 |
| Dimensions: 22(L) x 5.5(W) x 8.75"(H) | | |
| 20-Place (4L x 5H). 175HL-20-2 | \$132.00 | 119.00 |
| Dimensions: 22(L) x 5.5(W) x 11"(H) | | |
- Horizontal, Holds 3" Boxes**
- | | | |
|---|---------|-------|
| 9-Place (3L x 3H). 175HS-9-3 | \$61.00 | 55.00 |
| Dimensions: 16.31(L) x 5.5(W) x 9.5"(H) | | |
| 12-Place (4L x 3H). 175HL-12-3 | \$80.00 | 73.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H) | | |



Horizontal Drawer (fits upright freezers)

- Horizontal Drawer, Holds 2" Boxes**
- | | | |
|--|------------|------------|
| | <u>1-9</u> | <u>≥10</u> |
| 9-Place (3L x 3H). 175HD-332 | \$85.00 | 77.00 |
| Dimensions: 16.5(L) x 5.5(W) x 7.125"(H) | | |
| 12-Place (3L x 4H). 175HD-342 | \$111.00 | 101.00 |
| Dimensions: 16.5(L) x 5.5(W) x 9.5"(H) | | |
| 15-Place (3L x 5H). 175HD-352 | \$141.00 | 128.00 |
| Dimensions: 16.5(L) x 5.5(W) x 11.812"(H) | | |
| 16-Place (4L x 4H). 175HD-442 | \$151.00 | 137.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H) | | |
| 20-Place (5L x 4H). 175HD-542 | \$188.00 | 171.00 |
| Dimensions: 27.25(L) x 5.5(W) x 9.5"(H) | | |
- Horizontal Drawer, Holds 5" Boxes**
- | | | |
|---|------------|------------|
| | <u>1-9</u> | <u>≥10</u> |
| 6-Place (3L x 2H). 175RD-325 | \$74.00 | 67.00 |
| Dimensions: 18.125(L) x 6(W) x 10.75"(H) | | |
| 8-Place (4L x 2H). 175RD-425 | \$99.00 | 90.00 |
| Dimensions: 24(L) x 6(W) x 10.75"(H) | | |
- Horizontal Drawer, Holds 15 mL tubes**
- | | | |
|---|----------|--------|
| 160-Place (80L x 2H). 175HD-32ML-15 | \$99.00 | 90.00 |
| Dimensions: 16.5(L) x 5.5(W) x 9.5"(H) | | |
| 220-Place (110L x 2H). 175HD-42ML-15 | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H) | | |
- Horizontal Drawer, Holds 50 mL tubes**
- | | | |
|---|----------|--------|
| 84-Place (42L x 2H). 175HD-42ML-50 | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H) | | |

Continued

Horizontal Drawer-Tray (fits upright freezers)


Drawer-Tray, Holds 2" Boxes , aluminum trays	1-9	≥10
12-Place (3L x 4H). 175HDT-342	\$96.00	86.00
Dimensions: 16.63(L) x 5.5(W) x 8.5"(H)		
18-Place (3L x 6H). 175HDT-362	\$143.00	129.00
Dimensions: 16.63(L) x 5.5(W) x 12.75"(H)		
24-Place (4L x 6H). 175HDT-462	\$189.00	172.00
Dimensions: 22(L) x 5.5(W) x 12.75"(H)		

Vertical (fits chest freezers)

Vertical, Holds 2" Boxes		
13-Place. 175V-13-2L	\$88.00	79.00
Dimensions: 5.5(L) x 5.625(W) x 28.562"(H)		

Tube Storage Boxes (cryogenic) – see in GENERAL LAB SUPPLIES & EQUIPMENT section

Funnels



Funnels

Disposable, polypropylene, short stems, Simport		
★ 55 mm Top ID, takes 11 cm filter, Cs/100. F490-3	\$40.25	
★ 65 mm Top ID, takes 12.5 cm filter, Cs/100. F490-4	\$40.25	
Glass (borosilicate), long stem (150 mm)		
50 mm Top ID, Pk/6. 600GF614050	\$45.60	
65 mm Top ID, Pk/6. 600GF614065	\$54.90	
75 mm Top ID, Pk/6. 600GF614075	\$56.75	
Heavy-Duty, Molded-In Handle, polypropylene		
250 mm Top ID, 2.5 L capacity, each. 60057232	\$25.70	
Micro, polypropylene, Bel-Art		
24 mm Top ID, 2.6 mL capacity, Pk/12. 1146850024	\$41.05	
28 mm Top ID, 5.2 mL capacity, Pk/12. 146850028	\$33.60	
43 mm Top ID, 15.7 mL capacity, Pk/12. 146850043	\$70.80	
Powder, polypropylene		
Parallel Stem		
78 mm Top ID, Pk/12. 60057222	\$16.45	
100 mm Top ID, Pk/6. 146600100	\$15.95	
145 mm Top ID, Pk/6. 60057224	\$20.70	
Tapered Stem		
152 mm Top ID, each. 146820000	\$34.80	

Continued

Funnels (cont)

Standard Stem, Premium, polypropylene, Bel-Art	
Exclusive filter-step internal design provides speedier filtration; external ribs prevent air-lock.	
55 mm Top ID, takes 9 cm filter, Pk/12. 146760055	\$47.75
65 mm Top ID, takes 11 cm filter, Pk/12. 146760065	\$49.60
75 mm Top ID, takes 12.5 cm filter, Pk/6. 146760075	\$25.65
100 mm Top ID, takes 18.5 cm filter, Pk/4. 146760100	\$32.60
Standard Stem, Economy, polypropylene	
50 mm Top ID, each. 310FUNPP50	\$0.80
62 mm Top ID, each. 310FUNPP62	\$0.95
75 mm Top ID, each. 310FUNPP75	\$1.15
100 mm Top ID, each. 310FUNPP100	\$1.65

Gauze Pads

These 4 x 4" nonsterile gauze sponges are 12-ply.



Gauze Pads , 10 pk/200, Cs/2000. 44122000	\$88.45
---	---------

Germicidal Wipes

Sani-Cloth Plus	
Germicidal Wipes , Tub/160 440Q89072	\$10.00



Gloves & Hand Protection

Cotton Gloves	
Protect optics and photographs from fingerprints. Can also be used as glove liners.	
Small, Pk/12. GLV-190.10	\$5.25
Large, Pk/12. GLV-190.20	\$5.75
Purple Nitrile, powder-free, 5.5 mil, textured fingers, Halyard	
Nitrile gloves are hypoallergenic and provide greater puncture and chemical resistance than latex or vinyl gloves.	
Small, bx/100, Cs/10 Bx. 55081	\$139.00
Medium, bx/100, Cs/10 Bx. 55082	\$139.00
Large, bx/100, Cs/10 Bx. 55083	\$139.00
X-Large, bx/90, Cs/10 Bx. 55084	\$139.00

Nitrile Gloves



Gloves

Cotton

Continued