

Horizontal Drawer-Tray (fits upright freezers)


Drawer-Tray, Holds 2" Boxes , aluminum trays	1-9	≥10
12-Place (3L x 4H). 175HDT-342	\$96.00	86.00
Dimensions: 16.63(L) x 5.5(W) x 8.5"(H)		
18-Place (3L x 6H). 175HDT-362	\$143.00	129.00
Dimensions: 16.63(L) x 5.5(W) x 12.75"(H)		
24-Place (4L x 6H). 175HDT-462	\$189.00	172.00
Dimensions: 22(L) x 5.5(W) x 12.75"(H)		

Vertical (fits chest freezers)

Vertical, Holds 2" Boxes		
13-Place. 175V-13-2L	\$88.00	79.00
Dimensions: 5.5(L) x 5.625(W) x 28.562"(H)		

Tube Storage Boxes (cryogenic) – see in GENERAL LAB SUPPLIES & EQUIPMENT section

Funnels



Funnels

Disposable, polypropylene, short stems, Simport		
★ 55 mm Top ID, takes 11 cm filter, Cs/100. F490-3	\$40.25	
★ 65 mm Top ID, takes 12.5 cm filter, Cs/100. F490-4	\$40.25	
Glass (borosilicate), long stem (150 mm)		
50 mm Top ID, Pk/6. 600GF614050	\$45.60	
65 mm Top ID, Pk/6. 600GF614065	\$54.90	
75 mm Top ID, Pk/6. 600GF614075	\$56.75	
Heavy-Duty, Molded-In Handle, polypropylene		
250 mm Top ID, 2.5 L capacity, each. 60057232	\$25.70	
Micro, polypropylene, Bel-Art		
24 mm Top ID, 2.6 mL capacity, Pk/12. 1146850024	\$41.05	
28 mm Top ID, 5.2 mL capacity, Pk/12. 146850028	\$33.60	
43 mm Top ID, 15.7 mL capacity, Pk/12. 146850043	\$70.80	
Powder, polypropylene		
Parallel Stem		
78 mm Top ID, Pk/12. 60057222	\$16.45	
100 mm Top ID, Pk/6. 146600100	\$15.95	
145 mm Top ID, Pk/6. 60057224	\$20.70	
Tapered Stem		
152 mm Top ID, each. 146820000	\$34.80	

Continued

Funnels (cont)

Standard Stem, Premium, polypropylene, Bel-Art	
Exclusive filter-step internal design provides speedier filtration; external ribs prevent air-lock.	
55 mm Top ID, takes 9 cm filter, Pk/12. 146760055	\$47.75
65 mm Top ID, takes 11 cm filter, Pk/12. 146760065	\$49.60
75 mm Top ID, takes 12.5 cm filter, Pk/6. 146760075	\$25.65
100 mm Top ID, takes 18.5 cm filter, Pk/4. 146760100	\$32.60
Standard Stem, Economy, polypropylene	
50 mm Top ID, each. 310FUNPP50	\$0.80
62 mm Top ID, each. 310FUNPP62	\$0.95
75 mm Top ID, each. 310FUNPP75	\$1.15
100 mm Top ID, each. 310FUNPP100	\$1.65

Gauze Pads

These 4 x 4" nonsterile gauze sponges are 12-ply.



Gauze Pads , 10 pk/200, Cs/2000. 44122000	\$88.45
---	---------

Germicidal Wipes

Sani-Cloth Plus	
Germicidal Wipes , Tub/160 440Q89072	\$10.00



Gloves & Hand Protection

Cotton Gloves	
Protect optics and photographs from fingerprints. Can also be used as glove liners.	
Small, Pk/12. GLV-190.10	\$5.25
Large, Pk/12. GLV-190.20	\$5.75
Purple Nitrile, powder-free, 5.5 mil, textured fingers, Halyard	
Nitrile gloves are hypoallergenic and provide greater puncture and chemical resistance than latex or vinyl gloves.	
Small, bx/100, Cs/10 Bx. 55081	\$139.00
Medium, bx/100, Cs/10 Bx. 55082	\$139.00
Large, bx/100, Cs/10 Bx. 55083	\$139.00
X-Large, bx/90, Cs/10 Bx. 55084	\$139.00



Continued

Gloves (cont)

Purple Nitrile-XTRA, 12" Long, powder-free, 5.9 mil, Halyard
Offers additional protection against chemotherapy drugs and common chemicals.

- Small, bx/50, Cs/10 Bx. **50601**\$213.00
- Medium, bx/50, Cs/10 Bx. **50602**\$213.00
- Large, bx/50, Cs/10 Bx. **50603**\$213.00
- X-Large, bx/50, Cs/10 Bx. **50604**\$213.00


NEW



Safely grip hot or cold labware with this silicone rubber hand protector partially lined with mini suction cups. Withstands temperatures from 260°C (500°F) to -57°C (-70°F).

- Hand Protector, red. 600HNDPTC**\$19.70

Homogenizers



Rapidly homogenize, emulsify, suspend and disrupt 0.1 mL to 250 mL samples. Includes a 5 x 50 mm flat-bottom probe & 7 x 50 mm saw tooth probe for use with micro tubes.

- Homogenizer, Handheld. D1000**.....\$1,238.00



Homogenize a variety of sample types, from soft tissue to bone, in 35 seconds or less.
Benchmark BeadBlaster™ 24

- Homogenizer, Benchtop, Microtube. D2400**.....\$8,042.00

- Holds (24) 2.0 mL tubes, or (6) 5 mL tubes w/optional tube holder
- Speed adjusts from 4.00 to 7.00 m/s in 0.05 m/s increments
- Stores 50 programs in memory, each with a maximum of 10 cycles
- 2-year warranty

Homogenizer Prefilled Tubes

2.0 mL screw-capped tubes are prefilled with your choice of the following beads:

- Glass, acid washed
 - 0.1 mm, Pk/50. **D1031-01**\$132.00
 - 0.5 mm, Pk/50. **D1031-05**\$132.00
- Zirconium Beads, acid washed/heat treated/nuclease-free
 - 1.5 mm, Pk/50. **D1032-15**\$184.00
 - 3.0 mm, Pk/50. **D1032-30**\$184.00
- Stainless Steel Beads, acid washed
 - 2.8 mm, Pk/50. **D1033-28**\$184.00

Hot Plates, Magnetic Stirrers & Hot Plate Stirrer



- Hot Plates, analog**
- 6" Diameter Cast Iron Burner, 1000W. **310HOTPLT**\$31.60
- 7.5 x 7.5" Ceramic Platform, 80-380°C. **H4000-H**\$293.00



- Magnetic Stirrers, analog, Benchmark**
- 3.4 x 3.4" Platform, 100-1500 rpm, MiniMag™. **S1005**\$136.00 (1-year warranty)
- 7.1 x 7.1" Ceramic Platform, 60-1500 rpm. **H4000-S**\$293.00 (2-year warranty)

- Hot Plate Stirrer, 80-380°C, Benchmark. H4000-HS**\$438.00

Ice Pans



Temperature Range: -196 to 100°C

Ice Pans w/o lids, expanded PVC foam, Bel-Art, USA

- 1 L Capacity (Mini), 15.5 x 15.5 x 8.5 cm(H)
 - Blue. **188481101**\$40.65
 - Purple. **188481106**\$40.65
- 4 L Capacity (Midi), 33 x 24.5 x 12.5 cm(H)
 - Blue. **188484101**\$64.80
 - Purple. **188484106**\$64.80
- 9 L Capacity (Maxi), 44 x 35 x 13 cm(H)
 - Black. **188489102**\$76.20
 - Blue. **188489101**\$76.20