

**Pipets, Transfer, LDPE (low-density polyethylene)**

Graduated, Nonsterile, 15 cm long

- ★ 5 mL capacity, 1 mL graduated, 25 drops/mL
  - Bx/500. **PTP-02**.....\$21.30
  - Cs/10 Bx. **PTP-02**.....\$196.00
- 5 mL capacity, thick-walled
  - Cs/10 Bx. **188BP1200**.....\$164.00
- 7 mL capacity, 3 mL graduated, 14 drops/mL
  - Bx/500. **PTP-01**.....\$21.30
  - Cs/10 Bx. **PTP-01**.....\$196.00
- Graduated, Sterile Pk/1 (individually wrapped), 15 cm long
  - 5 mL capacity, 1 mL graduated, Simport
    - Bx/400. **P200-521S**.....\$49.85




These mini, air displacement pipettors provide an economical and accurate method for pipetting fixed volumes. Pipettors feature a two-step plunger and are fully autoclavable.

**Pipettors, Single Channel**

**Fixed Volume**, pocket-size, Diamond® Jr.

- 5 µL, red. **3322-5**.....\$27.75
- 10 µL, pink. **3322-10**.....\$27.75
- 20 µL, green. **3322-20**.....\$27.75
- 50 µL, purple. **3322-50**.....\$27.75
- 100 µL, orange. **3322-100**.....\$27.75

Single channel, adjustable volume pipettors are fully autoclavable.



**Adjustable Volume**, Diamond® Apex™ (3-year warranty)


- 0.5-10 µL, 0.02 µL increments, red. **3351-10**.....\$198.00
- 5-50 µL, 0.1 µL increments, yellow. **3351-50**.....\$198.00
- 10-100 µL, 0.1 µL increments, blue. **3351-100**.....\$198.00
- 100-1000 µL, 1 µL increments, orange. **3351-1000**.....\$198.00



**Pipettor Holders, Carousel, Flip & Grip™, Heathrow**

- Blue. **HS120026**.....\$92.45
- Lime. **HS120056**.....\$92.45

**Pitchers**



Polypropylene




Stainless Steel

**Pitchers**

- Polypropylene, graduated (mL), tall form, tapered
  - 1000 mL, 50 mL graduations, each. **289910000**.....\$7.00
  - 2000 mL, 100 mL graduations, each. **289920000**.....\$8.05
  - 3000 mL, 100 mL graduations, each. **289930000**.....\$13.80
- Stainless Steel, graduated (mL/oz), tapered, Polar Ware
  - 1000 mL/32 oz, each. **T1063**.....\$73.20

**Rockers, Shakers & Mixers**



**Rockers, Fixed Speed, Benchmark (2-year warranty)**

- 2-D (seesaw), 12 x 12" platform w/flat mat
  - BenchBlotter™. **BR1000**.....\$543.00
    - Fixed speed of 12 rocks/minute
    - Platform angle is ±8°
    - Cold room & incubator safe (4-45°C)
    - Load capacity is 2.3 kg (5 lbs)
- 3-D (nutating)
  - Small, 10.5 x 7.5" platform w/dimpled mat
    - MiniMixer™. **B3D1020**.....\$449.00
  - Large, 12 x 12" platform w/dimpled mat
    - BioMixer™. **B3D1320**.....\$558.00
      - Fixed speed of 24 rocks/minute
      - Platform angle is ±20°
      - Cold room & incubator safe (4-45°C)
      - Load capacity is 1.8 kg (4 lbs)



### Rockers, Variable Speed


Benchmark (2-year warranty)

- 2-D (seesaw), 14 x 12" platform w/flat mat  
 BenchRocker™ 2D. **BR2000** .....\$830.00  
 Stacking Platform. **BR2000-STACK**.....\$155.00  
 14 x 12" platform w/flat mat has 3" of clearance.  
 Stacking Platform Extenders, Pk/4. **BR2000-SP** .....\$55.00  
 Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 3.6 kg (8 lbs)

- 3-D (nutating), 12 x 12" platform w/dimpled and flat mats  
 BenchRocker™ 3D. **B3D2300** .....\$1,029.00  
 Stacking Platform. **BR1000-STACK**.....\$142.00  
 12 x 12" platform w/flat mat has 3" of clearance.  
 Stacking Platform Extenders, Pk/4. **BR2000-SP** .....\$55.00  
 Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)



### Rockers, Variable Speed

Reliable Scientific (5-year warranty), **USA**

- 2-D (seesaw)  
 Standard, 11 x 14" platform w/rubber mat  
 Single Platform. **55S1114** .....\$603.00  
 Double Platform, 3¼" apart. **55D1114**.....\$665.00  
 Large, 12 x 16" platform w/rubber mat  
 Single Platform. **55S1216**.....\$647.00  
 Double Platform, 3¼" apart. **55D1216**.....\$704.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



**Rocker, Tube**, 12 x 4.5" grooved mat, Benchmark. **M2100** .....\$405.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is ±30°
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs )
- 2-year warranty

### Orbital Shakers



Aerate and mix with these low to high capacity orbital shakers.

### Orbital Shakers, Benchmark (2-year warranty)

- Low Speed, 14 x 12" platform w/flat rubber mat  
 Orbi-Blotter™. **BT30** .....\$933.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

- Small, 9.5 x 8" platform w/flat rubber mat  
 Orbi-Shaker™ JR. **BT302** .....\$1,126.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min



BT3001

- Large, 13 x 12" platform w/flat rubber mat  
 Orbi-Shaker™. **BT3001** .....\$1,730.00

- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L, or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

*Optional Magic Clamp™ platform & accessories sold separately - specify configuration*

- Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR**....\$327.00

- Clamps, magnetic  
 250 mL Flask. **H1000-MR-250** .....\$46.10  
 500 mL Flask. **H1000-MR-500** .....\$46.10  
 1000 mL Flask. **H1000-MR-1000** .....\$70.50