


Tubes & Vials



- Tubes, Centrifuge, Conical** bottom, sterile, white writing patch
RNase-/DNase-/Pyrogen-free
- 15 mL, PP (polypropylene) tube/HDPE cap, natural
Rack/25, 20 racks, Cs/500. **BP1501**\$116.00
Maximum RCF is 12,000 x g. Frosted writing patch.
 - 15 mL, PS (polystyrene)
Rack/50, 10 racks, Cs/500. **6274**\$246.00
Maximum RCF is 3,000 x g.
 - 50 mL, natural, polypropylene (PP) tube/HDPE cap
Bag/25, 20 bags, Cs/500. **BP5002**\$128.00
Maximum RCF is 12,000 x g. Frosted writing patch.
 - Bag/25, 20 bags, Cs/500. **6297**\$200.00
Maximum RCF is 20,000 x g
 - Rack/25, 20 racks, Cs/500. **BP5001**\$135.00
Maximum RCF is 12,000 x g. Frosted writing patch.

Self-Standing



- Tubes, Centrifuge, Self-Standing** (skirted)
- 30 mL
Nonsterile, PP tube/PP cap, Cs/500. **109220**\$171.00
Dimensions: 23(OD) x 90 mm
 - Sterile, PP tube/HDPE cap
Cs/500, 20 bags/25. **C572-1**\$156.00
Molded-in graduations from 5 to 30 mL in 5 mL increments. Dimensions: 25.3(OD) x 111 mm

Plastic culture tubes with dual-position caps.



Tubes, Culture,

- 12 x 75 mm (5 mL)
Polystyrene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110428**\$89.50
- Polystyrene, Non sterile, Cs/1000. **110410**\$62.90
- Polypropylene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110438**\$89.70
- 13 x 100 mm (10 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1510**\$97.20
- 17 x 100 mm (16 mL)
Polypropylene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110178**\$143.00
- 16 x 125 mm (19 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1515**\$157.00
- 16 x 150 mm (23 mL)
★ #1 Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1517**\$161.00

Culture Tube Caps, polypropylene
Up to 16 mm Tubes




- Blue, Cs/1000. **118152B**\$99.65
- Green, Cs/1000. **118152G**\$99.65
- Red, Cs/1000. **118152R**\$99.65
- Yellow, Cs/1000. **118152Y**\$99.65



RNase-/DNase-free

- Tubes, Microcentrifuge, Economy**, polypropylene
- 0.5 mL, 0.1 mL graduations, natural, Pk/1000. **HS4422**\$38.45
 - 1.5 mL, 0.25 mL graduations, Heathrow
Natural, max RCF is 16,400 x g
Pk/500. **HS4323**\$18.70
Pk/1000. **HS4325**\$35.35
 - Black, max RCF is 18,000 x g, Pk/500. **HS4323K**\$20.80
 - 5 mL, natural, max RCF is 25,000 x g, Benchmark
Nonsterile, 4 bags/50, Pk/200. **C1005T5**\$65.40
Sterile, 10 bags/20, Pk/200. **C1005T5ST**\$71.45

Continued



Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable. Free of RNase/DNase.

Tubes, Microcentrifuge, Premium, PP, G-Tube™, Bio Plas, USA

Amber color protects light sensitive specimens.

- 0.6 mL, Pk/500.....\$27.75
Amber. **4048** Lavender. **4047** Natural. **4040**
Mixed (blue, green, orange, pink, & yellow). **4049**.....\$33.45
- 1.5 mL, 0.5 mL graduations, frosted area on side, Pk/500\$27.75
Amber. **4038** Blue. **4033** Lavender. **4037** Natural. **4030**
Mixed (blue, green, orange, pink, & yellow). **4039**.....\$33.45
- 2.0 mL, 0.5 mL graduations, frosted area on side
Natural, Pk/500. **4050**.....\$29.80



Tubes, PCR, 0.2 mL thin-wall, polypropylene, Simport

- Individual, flat caps, natural, Pk/1000. **T325-2VN**.....\$50.20



Vials, Cryogenic, round bottom, internal thread

- O-Ring Seal, 4 mL, Cs/1000. **T301-4**.....\$297.00
- Washer Seal, 4 mL, Cs/1000. **T311-4**.....\$297.00



Vinyl grip, **USA**

Vial Decapping Pliers

- 13 mm Crimp Seals. **D-13**.....\$72.00
- 20 mm Crimp Seals. **D-20**.....\$72.00

Tube Storage

Cryogenic Storage Boxes



SB2CC-64 Cryogenic Tube Storage Box

Tube Storage Boxes, Cardboard, vented

- 2" Box Only, 5¼ x 5¼ x 2", each. **SB2CC**.....\$1.80
- 2" Box w/dividers, white
 - 100-Cell (10 x 10), Pk/12. **SB2CC-100**.....\$28.45
 - ★ 81-Cell (9 x 9), Pk/12. **SB2CC-81**.....\$28.20
 - 64-Cell (8 x 8), Pk/12. **SB2CC-64**.....\$28.20
- 2" Mini box w/dividers, 3 x 3 x 2"
 - 25-Cell (5 x 5), Pk/12. **MB2C-25**.....\$30.65
 - Holds 13 mm vials and 4.0-5.0 mL cryovials.
- 3" Box Only, 5¼ x 5¼ x 3", each. **SB3C**.....\$1.80
- 3" Box w/dividers
 - 100-Cell (10 x 10), Pk/12. **SB3C-100**.....\$28.45
 - 81-Cell (9 x 9), Pk/12. **SB3C-81**.....\$28.20
 - 64-Cell (8 x 8), Pk/12. **SB3C-64**.....\$28.20
- Dividers only, Fit 2 & 3" Boxes
 - 100-Cell (10 x 10), Pk/12. **CD-100**.....\$11.85
 - Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.
 - 81-Cell (9 x 9), Pk/12. **CD-81**.....\$11.55
 - Holds 13 mm vials and 4.0-5.0 mL cryovials.
 - 64-Cell (8 x 8), Pk/12. **CD-64**.....\$11.55
 - Holds 14 mm vials.
 - 49-Cell (7 x 7), Pk/12. **CD-49**.....\$11.85
 - Holds 14 mm vials.
 - 25-Cell (5 x 5), 25 mm cell, Pk/12. **CD-25**.....\$11.55
 - 16-Cell (4 x 4), 31 mm cell, Pk/12. **CD-16**.....\$11.55
- 5" Box Only, 5¼ x 5¼ x 5", each. **SB4C**.....\$3.30
- 5" Box w/dividers
 - 64-Cell (8 x 8), Pk/12. **SB4C-64**.....\$43.50
 - 36-Cell (6 x 6), 15 mL tubes, Pk/12. **SB4C-36**.....\$43.50
 - 16-Cell (4 x 4), 50 mL tubes, Pk/12. **SB4C-16**.....\$43.50
- Dividers only, Fit 5" Box
 - 64-Cell (8 x 8), Pk/12. **CD4-64**.....\$11.55
 - 36-Cell (6 x 6), Pk/12. **CD4-36**.....\$11.55
 - 16-Cell (4 x 4), Pk/12. **CD4-16**.....\$11.55



Holds 1.0-2.0 mL tubes and cryovials. Features a locking hinged lid, and base with molded alphanumeric coordinates. Dimensions: 12.8 x 12.8 x 5.4 cm.

Tube Storage Boxes, Polypropylene

- 81-Well (9x9), 2 blue/2 yellow boxes, Pk/4. **60066501**.....\$33.30