

Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable. Free of RNase/DNase.

Tubes, Microcentrifuge, Premium, PP, G-Tube™, Bio Plas, USA

Amber color protects light sensitive specimens.

- 0.6 mL, Pk/500\$27.75
Amber. **4048** Lavender. **4047** Natural. **4040**
Mixed (blue, green, orange, pink, & yellow). **4049**\$33.45
- 1.5 mL, 0.5 mL graduations, frosted area on side, Pk/500\$27.75
Amber. **4038** Blue. **4033** Lavender. **4037** Natural. **4030**
Mixed (blue, green, orange, pink, & yellow). **4039**\$33.45
- 2.0 mL, 0.5 mL graduations, frosted area on side
Natural, Pk/500. **4050**\$29.80



Tubes, PCR, 0.2 mL thin-wall, polypropylene, Simport

- Individual, flat caps, natural, Pk/1000. **T325-2VN**\$50.20



Vials, Cryogenic, round bottom, internal thread

- O-Ring Seal, 4 mL, Cs/1000. **T301-4**\$297.00
- Washer Seal, 4 mL, Cs/1000. **T311-4**\$297.00



Vinyl grip, **USA**

Vial Decapping Pliers

- 13 mm Crimp Seals. **D-13**\$72.00
- 20 mm Crimp Seals. **D-20**\$72.00

Tube Storage

Cryogenic Storage Boxes



SB2CC-64 Cryogenic Tube Storage Box

Tube Storage Boxes, Cardboard, vented

- 2" Box Only, 5¼ x 5¼ x 2", each. **SB2CC**\$1.80
- 2" Box w/dividers, white
 - 100-Cell (10 x 10), Pk/12. **SB2CC-100**\$28.45
 - ★ 81-Cell (9 x 9), Pk/12. **SB2CC-81**\$28.20
 - 64-Cell (8 x 8), Pk/12. **SB2CC-64**\$28.20
- 2" Mini box w/dividers, 3 x 3 x 2"
 - 25-Cell (5 x 5), Pk/12. **MB2C-25**\$30.65
 - Holds 13 mm vials and 4.0-5.0 mL cryovials.
- 3" Box Only, 5¼ x 5¼ x 3", each. **SB3C**\$1.80
- 3" Box w/dividers
 - 100-Cell (10 x 10), Pk/12. **SB3C-100**\$28.45
 - 81-Cell (9 x 9), Pk/12. **SB3C-81**\$28.20
 - 64-Cell (8 x 8), Pk/12. **SB3C-64**\$28.20
- Dividers only, Fit 2 & 3" Boxes
 - 100-Cell (10 x 10), Pk/12. **CD-100**\$11.85
 - Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.
 - 81-Cell (9 x 9), Pk/12. **CD-81**\$11.55
 - Holds 13 mm vials and 4.0-5.0 mL cryovials.
 - 64-Cell (8 x 8), Pk/12. **CD-64**\$11.55
 - Holds 14 mm vials.
 - 49-Cell (7 x 7), Pk/12. **CD-49**\$11.85
 - Holds 14 mm vials.
 - 25-Cell (5 x 5), 25 mm cell, Pk/12. **CD-25**\$11.55
 - 16-Cell (4 x 4), 31 mm cell, Pk/12. **CD-16**\$11.55
- 5" Box Only, 5¾ x 5¾ x 5", each. **SB4C**\$3.30
- 5" Box w/dividers
 - 64-Cell (8 x 8), Pk/12. **SB4C-64**\$43.50
 - 36-Cell (6 x 6), 15 mL tubes, Pk/12. **SB4C-36**\$43.50
 - 16-Cell (4 x 4), 50 mL tubes, Pk/12. **SB4C-16**\$43.50
- Dividers only, Fit 5" Box
 - 64-Cell (8 x 8), Pk/12. **CD4-64**\$11.55
 - 36-Cell (6 x 6), Pk/12. **CD4-36**\$11.55
 - 16-Cell (4 x 4), Pk/12. **CD4-16**\$11.55



Holds 1.0-2.0 mL tubes and cryovials. Features a locking hinged lid, and base with molded alphanumeric coordinates. Dimensions: 12.8 x 12.8 x 5.4 cm.

Tube Storage Boxes, Polypropylene

- 81-Well (9x9), 2 blue/2 yellow boxes, Pk/4. **60066501**\$33.30

Tube Racks



Tube Racks, Grip, polypropylene, No-Wire™, Bel-Art
18-20 mm Tubes, 40-place, each. **HS290002** \$35.30



4-Way Rack

HS29022G




4-Way Rotatable Rack

HS29040A

Tube Racks, Interlocking, polypropylene, Heathrow
4-Way, mixed colors, Pk/5. **HS29022G** \$46.15
Holds (4) 50 mL tubes, (12) 15 mL tubes,
(32) 1.5 mL tubes, or (32) 0.5 mL tubes.
4-Way Rotatable, Rota-Rack™, four connected modules
Large Tubes, 22.4 x 11 x 11 cm, each. **HS29040A** \$47.80
Each module holds (2) 50 mL tubes, (8) 15 mL tubes,
(6) 20 mm tubes, or (10) 10-13 mm tubes. Mixed colors.




HS24322A/E




HS120112

Tube Racks, Mega Rack, polypropylene, fully-assembled
10-13 mm Tubes
216-Place, 33.7 x 26.7 x 9.1 cm(H)
Blue, each. **HS24322A** \$50.30
White, each. **HS24322E** \$50.30
432-Place, 67.4 x 26.7 x 9.1 cm(H)
Blue, each. **HS120113** \$87.70
White, each. **HS120112** \$87.70



80-Well




96-Well Reversible

Tube Racks, Microtube, polypropylene, Heathrow
80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm
Holds 1.5-2.0 mL tubes in conical wells.
Blue, Pk/5. **HS29025B** \$27.90
Pink, Pk/5. **HS29025D** \$27.90
Mixed, Pk/5. **HS29025G** \$27.90
96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.
Pink, Pk/5. **HS2345E** \$56.85

Available in blue,
green, orange,
purple & yellow.




Tube Rack, PCR, polypropylene, mixed colors, Heathrow
Reversible, Pk/5. **HS2344A** \$60.60
Holds (168) 0.2, (40) 0.5 and (12) 1.5 mL tubes.



Tube Racks, Peg, 14-17 mm Tubes, 50-place, polypropylene
White, Pk/2. **HS24311A** \$30.15



Half-Size



Full-Size

Tube Racks, Teardrop, OneRack™, POM, Heathrow
Half-Size
16 mm tubes, 36-place, 12.7 x 12.7 x 7 cm(H)
Blue, Pk/8. **HS27502B** \$115.00
Full-Size
13 mm tubes, 72-place, 20.2 x 10.4 x 5.9 cm(H)
Blue, Pk/8. **HS27511B** \$151.00
Magenta, Pk/8. **HS27511D** \$151.00
16 mm tubes, 72-place, 25 x 12.7 x 7 cm(H)
Blue, Pk/8. **HS27512B** \$151.00
Magenta, Pk/8. **HS27512D** \$151.00