

GENERAL LAB SUPPLIES & EQUIPMENT

Applicators




Applicators, 6" long

Cotton Tipped, wooden shaft

Nonsterile, 10 bg/bx, Bx/1000. 95-8702	\$14.05
Sterile Pk/1, Bx/100. 44024-106-1S	\$7.85
Sterile Pk/2, Bx/200. 44024-106-2S	\$9.80
Plain, wood, Bg/1000. 95-8701	\$8.60

Aprons



Protect yourself with this disposable 0.8-mil thick apron.

Aprons, disposable, polyethylene, 46"(L) x 28"(W)

0.8 mil Thick (0.02 mm), Bx/100. 105DA2846	\$16.65
1.75 mil Thick (0.044 mm), Bx/100. 59010485	\$38.00

Autoclaves

B4000-28-120

B4000-M



B4000-18

Autoclaves, benchtop, Benchmark (1-year warranty)

8 L, BioClave™ Mini. B4000-M	\$4,399.00↓
Chamber Dimensions: 6.7" ID x 12.5" (D)	
18 L, BioClave™ 18. B4000-18	\$6,019.00↓
Chamber Dimensions: 9.8" ID x 13.75" (D)	
28 L, BioClave™ 28. B4000-28-120	\$11,080.00
Chamber Dimensions: 12.6" ID x 16.4" (D)	

Bags



Collect biohazardous waste in this space-saving, autoclavable bag. Printed with a biohazard symbol and sterilization indicator patch. Measures 20.3 x 30.5 cm (8 x 12") and is 1.57 mil (0.04 mm) thick polypropylene.

HS10044H



Bags, Biohazard, Benchtop, 8 x 12", Bx/500. HS10320	\$102.00
Biohazard Bag Holder, Benchtop. HS10044H	\$27.00
Epoxy-coated wire holder measures 9.1 x 12.8 x 21.4 cm.	



Bags display a biohazard symbol and precautionary procedures. Sterilization indicator patch darkens when autoclaved. Polypropylene bags are bright red. Durable, extra heavy-duty bags provide maximum protection against dangerous tears and leaks.

Bags, Biohazard, Standard, 1.5 mil (0.038 mm) thick

14 x 19" (36 x 48 cm), Bx/200. HS10321	\$101.00
19 x 23" (48 x 58 cm), Bx/200. 131641923	\$144.00
25 x 35" (64 x 89 cm), Bx/200. 131642535	\$309.00

Bags, Biohazard, Heavy-Duty, 2.2 mil (0.055 mm) thick

12 x 18"(L), 9 L, Bx/200. 600UNU40105	\$62.00
12 x 24"(L), 12 L, Bx/200. 600UNU40106	\$74.25
18 x 22"(L), 16 L, Bx/200. 600UNU40103	\$89.00
25 x 35"(L), 37 L, Bx/200. 600UNU40104	\$155.00



Bag, Waste Disposal, HDPE, Cleanware™

8 x 10", Pk/50. 131741008	\$29.35
----------------------------------	---------



Seal polyethylene or polypropylene bags, sheets, or flat tubing from 0.01 to 0.8 mm thick. Sliding blade cuts excess material 19 mm from the seal.

Bag Sealer, electric, 12" (30.5 cm) Wide. HS28272	\$248.00
Bag Sealer Replacement Part Kit. HS28280	\$37.80

Balances

Portable



Model	DWP-1001
Capacity (g)	1000
Readability (g)	0.1
Modes	g, oz, ct, gn
Pan	14.6 cm (5.75") stainless steel
Power	4 AA/AC adapter
Price	\$72.00

Warranty: 1 year

Toploading

Benchmark
W3200-1200



Model	W3200-500	W3200-1200
Capacity (g)	500	1200
Readability (g)	0.001	0.01
Modes	g, oz, ct	g, oz, ct
Pan	11.5 cm (4.5")	16 cm (6.3")
Draft Shield	Yes, 3-door glass	No
Price	\$1,299.00	\$1,053.00

Warranty: 2 years

Analytical

Benchmark




Model	★ W3100-120	W3100A-120
Calibration	External (100 g weight provided)	Internal (Quick-Cal™ software)
Capacity (g)	120	120
Readability (g)	0.0001	0.0001
Modes	g, oz, lb, ct	g, oz, lb, ct
Pan	9 cm (3.5")	9 cm (3.5")
Price	\$1,948.00 ↓	\$2,067.00 ↓

Warranty: 2 years

Baths

BSH1002

Aluminum
Blocks




Dry Baths, digital, Benchmark (2-year warranty)

1-Block Capacity. BSH1001	\$527.00
2-Block Capacity. BSH1002	\$649.00
4-Block Capacity. BSH1004	\$1,029.00

Blocks sold separately – please specify configuration

Blocks, aluminum

Holds (48) 0.2 mL Tubes. BSW02	\$152.00
Holds (24) 0.5 mL Tubes. BSW05	\$152.00
Holds (24) 1.5 mL Tubes. BSW1500	\$152.00
Holds (24) 1.5/2.0 mL Tubes. BSW1520	\$152.00
Holds (20) 12/13 mL Tubes. BSW13	\$152.00
Holds (12) 15 mL Tubes. BSW15	\$152.00
Holds (5) 50 mL Tubes. BSW50	\$152.00



Racks are only included with model B2000-4.

Water Baths, digital, MyBath™, Benchmark (2-year warranty)

4 L, 3 tube racks. B2000-4	\$1,009.00
8 L. B2000-8	\$1,115.00
12 L. B2000-12	\$1,329.00

Beakers

Polypropylene
Beakers



Beakers

Glass (borosilicate), low form, double scale

Economy

10 mL, no graduations, Pk/12 600BG10	\$32.40
50 mL, Pk/12. 600BG50	\$39.10
100 mL, Pk/12. 600BG100	\$39.45
150 mL, Pk/12. 600BG150	\$42.25
250 mL, Pk/12. 600BG250	\$45.45
400 mL, Pk/12. 600BG400	\$57.25
600 mL, Pk/6. 600BG600	\$41.90
800 mL, Pk/6. 600BG800	\$75.00
1000 mL, Pk/6. 600BG1000	\$67.50
2000 mL, each. 600BG2000	\$21.20
4000 mL, each. 600BG4000	\$63.75

Heavy-Duty

50 mL, Pk/12. 600BG100350	\$40.60
100 mL, Pk/12. 600BG1003100	\$73.00
150 mL, Pk/12. 600BG1003150	\$81.00
250 mL, Pk/12. 600BG1003250	\$93.50
400 mL, Pk/12. 600BG1003400	\$109.00
600 mL, Pk/6. 600BG1003600	\$81.00
1000 mL, each. 600BG100310000	\$19.65

Polypropylene, premium, low form, single scale

100 mL, 25 mL graduations, each. 262110000	\$3.95
250 mL, 50 mL graduations, each. 262130000	\$6.55
600 mL, 100 mL graduations, each. 262150000	\$10.20
1000 mL, 200 mL graduations, each. 262160000	\$10.35
2000 mL, 500 mL graduations, each. 262170000	\$31.35
4000 mL, 1000 mL graduations, each. 262190000	\$52.25

NEW

Bench Protector



Protect your bench surface using this 40 mil absorbent liner backed with 0.5 mil polyethylene film.

Bench Liner, Roll, 50' x 20", white, Labmat™. **246750000** \$89.50

Bottles




Deliver controlled drops from standard 0.030" tip opening with these LDPE bottles. Sealer cap snaps on easily.

Bottles, Dispensing, squeeze bottles

LDPE (low-density polyethylene), Yorker cap

30 mL (1 oz), Bg/12. 116370001	\$27.75
185 mL (6 oz), Bg/12. 116370006	\$49.70

A long nozzle, textured sides, and slim profile make this chemically inert, HDPE container ideal for dispensing immersion oil.



HDPE (high-density polyethylene)

20 mL, each. 18145	\$8.10
---------------------------	--------

Amber glass protects light-sensitive liquids. Screw-cap dropper features a black rubber bulb and glass pipet.



Boston Round

Bottles, Dropping, Amber Glass (flint)

30 mL (1 oz), Bx/12. 3104053	\$18.35
Bottle Dropper Assembly (replacement)	
20/400, 7 x 76 mm Pipet, Pk/12. 3104053D	\$10.75
60 mL (2 oz), Bx/12. 3104054	\$26.85
Bottle Dropper Assembly (replacement)	
20/400, 7 x 89 mm Pipet, Pk/12. 3104054D	\$11.00

French Square



Bottles, Dropping, Clear Glass (flint)

15 mL (1/2 oz), Bx/12. 3104062	\$33.80
Bottle Dropper Assembly (replacement)	
20/400, 7 x 60 mm Pipet, Pk/12. 3104074	\$11.55
30 mL (1 oz), Bx/12. 3104063	\$38.50
Bottle Dropper Assembly (replacement)	
24/400, 7 x 75 mm Pipet, Pk/12. 3104075	\$14.70

Continued

Borosilicate glass bottle;
blue PP cap (GL45)



Bottles, Media, graduated, autoclavable

100 mL, Cs/10. 600BM0100	\$90.90
250 mL, Cs/10. 600BM0250	\$103.00
500 mL, Cs/10. 600BM0500	\$118.00
1000 mL, Cs/10. 600BM1000	\$168.00
Replacement Items	
Bottle Caps, GL45, blue, Pk/10. B3000-CAP	\$44.15
Bottle Drip-Free Rings, GL45, Pk/10. B3000-RIN	\$23.90



Globally Harmonized System
(GHS) wash bottles are vented
to prevent solvent drips.
Polyethylene bottles (LDPE)
with polypropylene (PP) caps.

Bottles, Wash, GHS Labeled, 500 mL (empty)

Acetone, Pk/6. 600370051	\$39.25
Bleach, Pk/6. 600370059	\$39.25
Deionized Water, Pk/4. 124160003	\$39.25
Ethanol, Pk/4. 600370052	\$39.25
Methanol, Pk/6. 600370053	\$39.25
Water, Pk/6. 600370061	\$39.25

Bottles, Wash, Safety Labeled, 500 mL (empty)

4-Color diamond, LDPE bottle/PP cap	
Soap, Pk/4. 117160014	\$39.25
Universal (label your own bottle), Pk/4. 117160009	\$39.25↓



Bottles, Wash, Integral Spout, LDPE bottle/PP cap

250 mL (8 oz), Cs/10. 601633	\$70.70
500 mL (16 oz), Cs/10. 601634	\$77.90↓



Bottles, Narrow-Mouth, HDPE bottle/polypropylene cap

125 mL, Pk/12. 106200015	\$24.50
250 mL, Pk/12. 106200016	\$35.20
500 mL, 28 mm cap, Pk/12. 106200007	\$40.95
1000 mL, Pk/6. 106200008	\$32.30



Bottles, Wide-Mouth

HDPE bottle/polypropylene cap	
125 mL, 38 mm cap, Pk/12. 60033407	\$24.50
250 mL, Pk/12. 106250006	\$35.20
500 mL, 53 mm cap, Pk/12. 60033409	\$40.95
1000 mL, Pk/6. 60033410	\$34.55
Polypropylene bottle and cap, autoclavable	
125 mL, 38 mm cap, Pk/12. 60033307	\$24.20
250 mL, 43 mm cap, Pk/12. 60033308	\$35.20
500 mL, 53 mm cap, Pk/12. 106320007	\$40.70
1000 mL, 63 mm cap, Pk/6. 60033310	\$35.95



Bottle Top Dispensers

0.5-5 mL. 600BTDR2	\$293.00
5-60 mL. 600BTDR5	\$349.00

Boxes



Matchbox Style

Box, Matchbox Style (Pill Box), cardboard, white

2 1/8 x 1 1/8 x 5/8", Pk/72. 20047	\$35.45
---	---------

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. **8342**

Each cardboard box can store approximately 540 paraffin blocks.




Boxes, Glass Disposal, cardboard, polyethylene bag included

Bench Model, 8 x 8 x 10", Pk/6. 246530002	\$116.00
Floor Model, 12 x 12 x 27", Pk/6. 246530001	\$160.00
Glass Disposal Box Cover, reusable, PP. 132040001	\$52.65
Replacement item for the floor model. Nonrefundable. Polypropylene cover is autoclavable.	

Continued

Hinged
Polystyrene



Boxes, Plastic, (see following chart)

Plastic Boxes				
Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
H-1	1 x 1 x ¼"	-	Reticles, Pointers	\$0.40
H-1/11B	1 x 1 x ¼"	1	Reticles, Pointers	\$0.50
H-54	2 9/16 x 2 9/16 x ½"	-	Filters	\$1.20
H-54/2916G	2 9/16 x 2 9/16 x ½"	1	Filters	\$1.35
H-28/334B2	3 ¾ x 1 ¾ x 5/8"	2	Slides, Stage Micrometers	\$1.40
H-10	2 7/8 x 1 ¾ x 3/8"	-	Small Parts	\$0.60
H-21/512B	5 ½ x 1 ¾ x ½"	1	Forceps & Scissors	\$1.80

Brushes



Brushes

- Bench, horse hair, 10" long. **BRS-969.10**\$9.95
 Bristle head: 4¾"(L) x 1½"(W); bristle length: 1"
 Flat Bristles (paint style)
 Camel Hair, 1" wide, 6½" total length. **BR665B (LB-37)**\$6.65
 Horsehair, heavy-duty bristles, metal handle. **33761880**\$0.50
 Round Bristles, camel hair, 8½" long. **BR661 (LB-30)**\$11.25
 Swivel-Head, black nylon bristles. **BR565 (LB-14)**\$4.70
 Bristle head: 4¼"(L) x 1½"(W); total length: 12"

Remove debris, rust, or corrosion in hard to reach places.



Pen Style (scratch), retractable, 11.5 cm (4½") long, Germany

- Brass. **BRS-295.00**\$12.20
 Nylon. **BRS-293.00**\$9.70
 Steel. **BRS-292.00**\$12.80

Burners



- Burner, Alcohol**, 120 mL, glass body, metal cap. **97-5313**\$7.50
 Wicks (replacement), 10-14 cm x 8 mm, Pk/12. **97-5310**\$5.45



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing. Height is 15.6 cm (6¼")

H-5550 gridded top, 3.8 cm (1½") diameter, produces a uniform blue cone of light. Height is 21.6 cm (8½")



Burners, Bunsen, natural gas, Humboldt Mfg, USA

- Standard, flame stabilizer. **H-5870**\$32.65
 High Temperature, Meker-type. **H-5550**\$98.80
 Accessory
 Bunsen Burner Natural Gas Connector
 3' Long. **H-2087**\$39.15



- Burner Lighter**, replaceable single flint, round file. **SOL-290.00**\$3.95
 Flints (replacement), Pk/5. **SOL-290.01**\$2.85

Carboys



Carboys with Spigot, polypropylene, round

- 10 L (2.5 gal). **60034051**\$80.95
 20 L (5 gal). **60034052**\$112.00

Continued

Centrifuges

GCH-24



Centrifuge, Hematocrit, 24-place. GCH-24.....\$2,147.00
24-place capillary rotor and reading disk lid are included.



C1012



C1008-C

Centrifuges, Mini

5,500 RPM (2000 x g), 12-place, MyFuge™ 12. **C1012**\$439.00
(1-year warranty)

6,000 RPM (2000 x g), 8-place, MyFuge™ Mini. **C1008-__** * \$333.00
*Specify clear (-C), green (-G), or purple (-P) for lid color.
(1-year warranty)

Clamps & Supports

Rectangular
Steel Base



Triangular
Cast-Iron Base



Support (Ring) Stands, finished metal base/plated support rod

Rectangular Steel Base

6 x 4" base, 5/16 x 18" rod. **3104SSS66**\$15.15

8 x 5" base, 3/8 x 20" rod. **3104SSS67**\$18.60

9 x 6" base, 1/2 x 24" rod. **3104SSS68**\$22.75

Rectangular Cast-Iron Base

Provides greater stability than steel base.

6 x 4" base, 5/16 x 18" rod. **3104SSC14**\$17.00


8 x 5" base, 3/8 x 20" rod. **3104SSC15**\$19.70

9 x 6" base, 1/2 x 24" rod. **3104SSC16**\$27.45

Triangular Cast-Iron Base

5" Legs, 1/2 x 24" rod. **1257G2324**\$25.00

6" Legs, 1/2 x 36" rod. **1257G23**\$30.65



600COBR3FR




600THCL01



600CLBURDD



97-4023



1257G49

Support Clamps

Burette, plastic coated jaws, 360° rotatable. **600COBR3FR**\$7.95

Burette, Double, deluxe. **600CLBURDD**\$21.75

Plastic covered spring jaws on a metal frame securely grip
burettes up to 30 mm OD. Fits rods up to 3/4" (19 mm) OD.

Hook Collar, nickel plated. **1257G49**\$6.20

Fits rods up to 1/2" (13 mm) in diameter.

Right-Angle, aluminum. **97-4023**\$7.85

Fits rods up to 5/8" in diameter.

Thermometer, aluminum, rotates 360°. **600THCL01**\$10.15

Fits rods up to 1/2" (13 mm) in diameter.



Support Rings, plated cast-iron

3" OD (2 1/2" ID). **600SRC103**\$8.70

4" OD (3 1/2" ID). **600SRC104**\$11.35

5" OD (4 1/2" ID). **600SRC105**\$12.10



Plain



Ceramic
Center

Wire Gauze Squares, galvanized iron wire w/turned edges

Plain

4 x 4". **31040004**\$1.10

5 x 5". **31040005**\$1.30

6 x 6". **31040006**\$1.55


Ceramic Center

4 x 4". **31040001**\$1.55

5 x 5". **31040002**\$1.70


6 x 6". **31040003**\$2.20

Continued



Wire Triangles, clay pipe, galvanized iron wire

1½" Pipe Length. 600CPT015	\$1.65
2" Pipe Length. 600CPT020	\$1.95
2½" Pipe Length. 600CPT025	\$2.10



Test Tube Clamp, brass wire, 13 cm (5"). **600TTCL02**

\$2.50



Hoffman



Mohr

Tubing Clamps, nickel plated

Hoffman, swivel lower jaw, Pk/3. 97-1801	\$7.80
Mohr Pinchcock, each. 3104MC42	\$2.40

Counters



Counters

Colony, handheld electronic. HS120000	\$290.00
--	----------



Handheld Tally, 0-9999 range. 600A	\$8.80
---	--------

Crucibles



Crucibles, Tall Form with Lid, porcelain

5 mL Capacity, Pk/6. 600JCT005	\$25.10
10 mL Capacity, Pk/6. 600JCT010	\$28.70
15 mL Capacity, Pk/6. 600JCT015	\$32.15
30 mL Capacity, Pk/6. 600JCT030	\$40.05
50 mL Capacity, Pk/6. 600JCT050	\$46.50
100 mL Capacity, Pk/6. 600JCT100	\$65.35


Cylinders

Cylinders, graduated

Glass (borosilicate),

Class B (economy), Single Scale

10 mL, 0.2 mL grad., ea. 310CSS10	\$3.45
25 mL, 0.5 mL grad., ea. 310CSS25	\$4.00
50 mL, 1.0 mL grad., ea. 310CSS50	\$4.65
100 mL, 1.0 mL grad., ea. 310CSS100	\$5.85
250 mL, 2.0 mL grad., ea. 310CSS250	\$9.10
500 mL, 5.0 mL grad., ea. 310CSS500	\$14.35
1000 mL, 10.0 mL grad., ea. 310CSS1000	\$23.60




Class A, individually certified, double metric scale

50 mL, 1.0 mL grad., ea. 600CY302050	\$41.50
100 mL, 1.0 mL grad., ea. 600CY3020100	\$49.40
250 mL, 2.0 mL grad., ea. 600CY3020250	\$72.95
500 mL, 5.0 mL grad., ea. 600CY3020500	\$94.50
1000 mL, 10.0 mL grad., ea. 600CY30201000	\$179.00

Polypropylene, single scale

10 mL, 0.1 mL grad., ea. 284520000	\$8.05
50 mL, 1.0 mL grad., ea. 284540000	\$10.70
100 mL, 1.0 mL grad., ea. 601079	\$13.10
250 mL, 2.0 mL grad., ea. 284560000	\$14.55
500 mL, 5.0 mL grad., ea. 284570000	\$20.00
1000 mL, 10.0 mL grad., ea. 601082	\$29.90
2000 mL, 20.0 mL grad., ea. 601094	\$55.40

Desiccators



Non-Vacuum



Vacuum

Glazed porcelain plate has 5 mm perforations with a 25 mm center hole.

Desiccators

Non-Vacuum, borosilicate glass

8" (20 cm) ID. 600DSGL210	\$147.00
--	----------

Vacuum, polycarbonate top/polypropylene bottom

6" (15 cm), 13.8 cm plate diameter. 60055204	\$77.00
8" (20 cm), 19 cm plate diameter. 60055205	\$96.70

10" (25 cm), 23.8 cm plate diameter. 60055206	\$123.00
--	----------

Drying Racks



189330023



HS23243A

Drying Racks, Benchtop, Bel-Art

Single-sided. **188050000**\$170.00

Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 43 x 17 x 48 cm(H).

Double-sided

2-Tier. **189330022**\$391.00

3-Tier. **189330023**\$539.00

NEW

Drying Rack, Benchtop/Wallmount, Heathrow

Single-sided. **HS23243A**\$146.00

Dry up to 90 pieces of labware with this self-standing or wall-mounted rack. Constructed of vinyl-coated wire hairpin loops with a polystyrene tray and rack. Dimensions: 46.2 x 18.2 x 52.5 cm(H).

Electrophoresis



Electrophoresis System, MyGel™ Mini. E1101\$607.00

Eyewear

Absorbs
99.9% UV



Glasses, Safety, integral side shield. **43574002**\$6.00



Flexible frame and adjustable elastic strap provides a tight seal.

Goggles, Safety, vented (chemical splash). **97-2001**\$4.00

Filter & Blotting Papers



Ahlstrom - Munksjo
Filter Paper

Qualitative, Standard Grades (for gravity filtration)

Slow flow (4 mL/min), 1.5 µm retention, Grade 610

(Whatman 5; S&S 602)

9 cm, Pk/100. **6100-0900**\$12.20

12.5 cm, Pk/100. **6100-1250**\$15.25

15 cm, Pk/100. **6100-1500**\$18.40

20.5 cm, Pk/100. **6100-2050**\$29.05

24 cm, Pk/100. **6100-2400**\$40.45

Medium flow (30 mL/min), 2 µm retention, Grade 642

(Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)

9 cm, Pk/100. **6420-0900**\$9.50

15 cm, Pk/100. **6420-1500**\$17.05

Medium flow (55 mL/min), 2.5 µm retention, Grade 601

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

4.25 cm, Pk/100. **6010-0425**\$7.25

9 cm, Pk/100. **6010-0900**\$9.80

12.5 cm, Pk/100. **6010-1250**\$14.10

15 cm, Pk/100. **6010-1500**\$17.15

18.5 cm, Pk/100. **6010-1850**\$24.85

20.5 cm, Pk/100. **6010-2050**\$26.80

33 cm, Pk/100. **6010-3300**\$72.40

Medium flow (60 mL/min), 6 µm retention, Grade 613

(Whatman Student; S&S 595, 0903)

9 cm, Pk/100. **6130-0900**\$8.35

15 cm, Pk/100. **6130-1500**\$14.50

20.5 cm, Pk/100. **6130-2050**\$21.80

33 cm, Pk/100. **6130-3300**\$48.45

Medium flow (200 mL/min), 10 µm retention, Grade 631

(Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)

20.5 cm, Pk/100. **6310-2050**\$27.00

33 cm, Pk/100. **6310-3300**\$73.15

Fast flow (235 mL/min), 25 µm retention, creped, Grade 615

9 cm, Pk/100. **6150-0900**\$8.60

12.5 cm, Pk/100. **6150-1250**\$12.50

15 cm, Pk/100. **6150-1500**\$15.00

20.5 cm, Pk/100. **6150-2050**\$22.50

33 cm, Pk/100. **6150-3300**\$49.90

Very Fast flow (360 mL/min), 35 µm retention, creped, Grade 617

9 cm, Pk/50. **6170-0900**\$7.85

12.5 cm, Pk/50. **6170-1250**\$10.45

15 cm, Pk/50. **6170-1500**\$12.15

20.5 cm, Pk/50. **6170-2050**\$18.55

33 cm, Pk/50. **6170-3300**\$39.80

Qualitative, Pleated (Fluted)

Slow flow (30 mL/min), 2 µm retention, Grade 542

(Whatman 2V)

12.5 cm, Pk/100. **5420-1250**\$78.15

18.5 cm, Pk/100. **5420-1850**\$88.45

Continued

Qualitative, Wet-Strengthened Grades (for vacuum filtration)

Medium flow (65 mL/min), 3 µm retention, **Grade 909**

(Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)

32 cm, Pk/100. **9090-3200**\$78.85

Fast flow (245 mL/min), 27 µm retention, creped, **Grade 962**

(Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)

11 cm, Pk/100. **9620-1100**\$11.70

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

Slow flow (4 mL/min), 1.5 µm retention, **Grade 94**

(Whatman 42; S&S 589/3)

9 cm, Pk/100. **0940-0900**\$33.00

15 cm, Pk/100. **0940-1500**\$56.85

Medium flow (30 mL/min), 2 µm retention, **Grade 74**

(Whatman 40; S&S 589/2, 589/4)

9 cm, Pk/100. **0740-0900**\$33.00

15 cm, Pk/100. **0740-1500**\$56.85

Fast flow (220 mL/min), 15 µm retention, **Grade 55**

(Whatman 54, 541; S&S 1505)

9 cm, Pk/100. **0550-0900**\$37.90

15 cm, Pk/100. **0550-1500**\$59.95

Blotting/Chromatography/Electrophoresis Papers

0.19 mm thick, 125 mm/30 min absorption, **Grade 601**

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

46 x 57 cm, Pk/100. **6018-4657**\$147.00

0.34 mm thick, 140 mm/30 min absorption, **Grade 238**

(Whatman 3MM, S&S 598, 231⁺)

20 x 20 cm, Pk/100. **2388-2020**\$75.25

46 x 57 cm, Pk/100. **2388-4657**\$286.00

0.42 mm thick, 150 mm/30 min absorption, **Grade 237**

(Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703;

Albet 2040B, FP591; Munktell 3S/H, FN7a)

46 x 57 cm, Pk/100. **2378-4657**\$329.00

0.95 mm thick, 180 mm/30 min absorption, **Grade 243**

(Whatman 17CHR; S&S 2668)

11 x 14 cm, Pk/100. **2438-1114**\$31.75

20 x 20 cm, Pk/100. **2438-2020**\$62.70

2.48 mm thick, 240 mm/30 min absorption, **Grade 320**

(S&S 2727)

11 x 14 cm, Pk/50. **3208-1114**\$60.10


20 x 20 cm, Pk/50. **3208-2020**\$87.45



Smooth tips &
polypropylene grips

Forceps, Filter, 11.5 cm. 300-025\$4.40

Freezer Racks



Freezer Racks, stainless steel, USA

175V-13-2L

Restocking fee of 20% on all returned freezer racks listed below.

No returns on custom racks.

Horizontal (fits upright freezers)

Horizontal, Holds 2" Boxes

	1-9	≥10
9-Place (3L x 3H). 175HS-9-2	\$59.00	53.00
Dimensions: 16.31(L) x 5.5(W) x 6.5"(H)		
12-Place (3L x 4H). 175HS-12-2	\$78.00	71.00
Dimensions: 16.31(L) x 5.5(W) x 8.75"(H)		
15-Place (3L x 5H). 175HS-15-2	\$98.00	89.00
Dimensions: 16.31(L) x 5.5(W) x 11"(H)		
16-Place (4L x 4H). 175HL-16-2	\$105.00	95.00
Dimensions: 22(L) x 5.5(W) x 8.75"(H)		
20-Place (4L x 5H). 175HL-20-2	\$132.00	119.00
Dimensions: 22(L) x 5.5(W) x 11"(H)		

Horizontal, Holds 3" Boxes

9-Place (3L x 3H). 175HS-9-3	\$61.00	55.00
Dimensions: 16.31(L) x 5.5(W) x 9.5"(H)		
12-Place (4L x 3H). 175HL-12-3	\$80.00	73.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		



Horizontal Drawer (fits upright freezers)

Horizontal Drawer, Holds 2" Boxes

	1-9	≥10
9-Place (3L x 3H). 175HD-332	\$85.00	77.00
Dimensions: 16.5(L) x 5.5(W) x 7.125"(H)		
12-Place (3L x 4H). 175HD-342	\$111.00	101.00
Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)		
15-Place (3L x 5H). 175HD-352	\$141.00	128.00
Dimensions: 16.5(L) x 5.5(W) x 11.812"(H)		
16-Place (4L x 4H). 175HD-442	\$151.00	137.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		
20-Place (5L x 4H). 175HD-542	\$188.00	171.00
Dimensions: 27.25(L) x 5.5(W) x 9.5"(H)		

Horizontal Drawer, Holds 5" Boxes

	1-9	≥10
6-Place (3L x 2H). 175RD-325	\$74.00	67.00
Dimensions: 18.125(L) x 6(W) x 10.75"(H)		
8-Place (4L x 2H). 175RD-425	\$99.00	90.00
Dimensions: 24(L) x 6(W) x 10.75"(H)		

Horizontal Drawer, Holds 15 mL tubes

160-Place (80L x 2H). 175HD-32ML-15	\$99.00	90.00
Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)		
220-Place (110L x 2H). 175HD-42ML-15	\$129.00	118.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		

Horizontal Drawer, Holds 50 mL tubes

84-Place (42L x 2H). 175HD-42ML-50	\$129.00	118.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		

Continued

Horizontal Drawer-Tray (fits upright freezers)

Drawer-Tray, Holds 2" Boxes, aluminum trays	1-9	≥10
12-Place (3L x 4H). 175HDT-342	\$96.00	86.00
Dimensions: 16.63(L) x 5.5(W) x 8.5"(H)		
18-Place (3L x 6H). 175HDT-362	\$143.00	129.00
Dimensions: 16.63(L) x 5.5(W) x 12.75"(H)		
24-Place (4L x 6H). 175HDT-462	\$189.00	172.00
Dimensions: 22(L) x 5.5(W) x 12.75"(H)		

Vertical (fits chest freezers)

Vertical, Holds 2" Boxes

13-Place. 175V-13-2L	\$88.00	79.00
Dimensions: 5.5(L) x 5.625(W) x 28.562"(H)		

Tube Storage Boxes (cryogenic) – see in GENERAL LAB SUPPLIES & EQUIPMENT section

Funnels



Funnels

Disposable, polypropylene, short stems, Simport		
★ 55 mm Top ID, takes 11 cm filter, Cs/100. F490-3	\$40.25	
★ 65 mm Top ID, takes 12.5 cm filter, Cs/100. F490-4	\$40.25	
Glass (borosilicate), long stem (150 mm)		
50 mm Top ID, Pk/6. 600GF614050	\$45.60	
65 mm Top ID, Pk/6. 600GF614065	\$54.90	
75 mm Top ID, Pk/6. 600GF614075	\$56.75	
Heavy-Duty, Molded-In Handle, polypropylene		
250 mm Top ID, 2.5 L capacity, each. 60057232	\$25.70	
Micro, polypropylene, Bel-Art		
24 mm Top ID, 2.6 mL capacity, Pk/12. 1146850024	\$41.05	
28 mm Top ID, 5.2 mL capacity, Pk/12. 146850028	\$33.60	
43 mm Top ID, 15.7 mL capacity, Pk/12. 146850043	\$70.80	
Powder, polypropylene		
Parallel Stem		
78 mm Top ID, Pk/12. 60057222	\$16.45	
100 mm Top ID, Pk/6. 146600100	\$15.95	
145 mm Top ID, Pk/6. 60057224	\$20.70	
Tapered Stem		
152 mm Top ID, each. 146820000	\$34.80	

Funnels (cont)

Standard Stem, Premium, polypropylene, Bel-Art

Exclusive filter-step internal design provides speedier filtration; external ribs prevent air-lock.

55 mm Top ID, takes 9 cm filter, Pk/12. **146760055**

65 mm Top ID, takes 11 cm filter, Pk/12. **146760065**

75 mm Top ID, takes 12.5 cm filter, Pk/6. **146760075**

100 mm Top ID, takes 18.5 cm filter, Pk/4. **146760100**

Standard Stem, Economy, polypropylene

50 mm Top ID, each. **310FUNPP50**

62 mm Top ID, each. **310FUNPP62**

75 mm Top ID, each. **310FUNPP75**

100 mm Top ID, each. **310FUNPP100**

Gauze Pads

These 4 x 4" nonsterile gauze sponges are 12-ply.



Gauze Pads, 10 pk/200, Cs/2000. **44122000**

Germicidal Wipes

Sani-Cloth Plus



Germicidal Wipes, Tub/160 **440Q89072**

Gloves & Hand Protection

Cotton Gloves



Nitrile Gloves



Gloves

Cotton

Protect optics and photographs from fingerprints.

Can also be used as glove liners.

Small, Pk/12. **GLV-190.10**

Large, Pk/12. **GLV-190.20**

Purple Nitrile, powder-free, 5.5 mil, textured fingers, Halyard

Nitrile gloves are hypoallergenic and provide greater

puncture and chemical resistance than latex or vinyl gloves.

Small, bx/100, Cs/10 Bx. **55081**

Medium, bx/100, Cs/10 Bx. **55082**

Large, bx/100, Cs/10 Bx. **55083**

X-Large, bx/90, Cs/10 Bx. **55084**

Continued

Continued

Gloves (cont)

Purple Nitrile-XTRA, 12" Long, powder-free, 5.9 mil, Halyard
Offers additional protection against chemotherapy drugs and common chemicals.

Small, bx/50, Cs/10 Bx. 50601	\$213.00
Medium, bx/50, Cs/10 Bx. 50602	\$213.00
Large, bx/50, Cs/10 Bx. 50603	\$213.00
X-Large, bx/50, Cs/10 Bx. 50604	\$213.00



Safely grip hot or cold labware with this silicone rubber hand protector partially lined with mini suction cups. Withstands temperatures from 260°C (500°F) to -57°C (-70°F).

Hand Protector, red. 600HNDPTC\$19.70

Homogenizers



Rapidly homogenize, emulsify, suspend and disrupt 0.1 mL to 250 mL samples. Includes a 5 x 50 mm flat-bottom probe & 7 x 50 mm saw tooth probe for use with micro tubes.

Homogenizer, Handheld. D1000\$1,238.00



Homogenize a variety of sample types, from soft tissue to bone, in 35 seconds or less.
Benchmark BeadBlaster™ 24

Homogenizer, Benchtop, Microtube. D2400\$8,042.00

- Holds (24) 2.0 mL tubes, or (6) 5 mL tubes w/optional tube holder
- Speed adjusts from 4.00 to 7.00 m/s in 0.05 m/s increments
- Stores 50 programs in memory, each with a maximum of 10 cycles
- 2-year warranty

Homogenizer Prefilled Tubes

2.0 mL screw-capped tubes are prefilled with your choice of the following beads:

Glass, acid washed

0.1 mm, Pk/50. **D1031-01**\$132.00

0.5 mm, Pk/50. **D1031-05**\$132.00

Zirconium Beads, acid washed/heat treated/nuclease-free

1.5 mm, Pk/50. **D1032-15**\$184.00

3.0 mm, Pk/50. **D1032-30**\$184.00

Stainless Steel Beads, acid washed

2.8 mm, Pk/50. **D1033-28**\$184.00

Hot Plates, Magnetic Stirrers & Hot Plate Stirrer

310HOTPLT



H4000-H

Hot Plates, analog

6" Diameter Cast Iron Burner, 1000W. **310HOTPLT**\$31.60

7.5 x 7.5" Ceramic Platform, 80-380°C. **H4000-H**\$293.00



S1005

H4000-S



Magnetic Stirrers, analog, Benchmark

3.4 x 3.4" Platform, 100-1500 rpm, MiniMag™. **S1005**\$136.00
(1-year warranty)

7.1 x 7.1" Ceramic Platform, 60-1500 rpm. **H4000-S**\$293.00
(2-year warranty)

Hot Plate Stirrer, 80-380°C, Benchmark. H4000-HS\$438.00

Ice Pans



Temperature Range: -196 to 100°C

Ice Pans w/o lids, expanded PVC foam, Bel-Art, USA

1 L Capacity (Mini), 15.5 x 15.5 x 8.5 cm(H)

Blue. **188481101**\$40.65

Purple. **188481106**\$40.65

4 L Capacity (Midi), 33 x 24.5 x 12.5 cm(H)

Blue. **188484101**\$64.80

Purple. **188484106**\$64.80

9 L Capacity (Maxi), 44 x 35 x 13 cm(H)

Black. **188489102**\$76.20

Blue. **188489101**\$76.20

Incubators



Incubators, Mini Digital, MyTemp™, Benchmark (2-year warranty)

Heat Only, 1° above ambient to 60°C. **H2200-H**\$569.00

Heat/Cool, 15° below ambient to 60°C. **H2200-HC**\$857.00

Heat/Cool/CO₂, 10° below ambient to 55°C. **H2300-HC2** ..\$3,990.00 ↓

- 20 L capacity (0.8 cu ft)
- Digital temperature control with LED display
- Accepts bottles and flasks up to 2 L
- Two adjustable/removable shelves
- 9.3"(D) x 10.3"(W) x 12.8"(H) interior;
14.5"(D) x 13.2"(W) x 18.7"(H) exterior



Incubators, Dual-Convection, SureTemp™ (2-year warranty)

40 L. **H2505-40**\$1,854.00

13.7"(D) x 13.7"(W) x 13.7"(H) interior

70 L. **H2505-70**\$2,363.00

13.7"(D) x 17.7"(W) x 17.7"(H) interior

130 L. **H2505-130**\$3,116.00

17.7"(D) x 21.7"(W) x 21.7"(H) interior

- Precisely control temperature up to 75°C
- Record temperatures with SureCheck™ data logging
- Select either gravity or mechanical convection mode

NEW



Incubator, Shaking, Incu-Shaker™ Mini. H1001-M.....\$3,199.00↓
9.5" x 11.5" platform w/ rubber mat

- 5° above ambient to 70°C
- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Domed lid accepts flasks up to 1 L
- 16.2"(D) x 11"(W) x 13"(H) exterior

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

Steel Platform, Magic Clamp™, 9.5 x 11". **H1000-MR**\$302.00

Clamps, magnetic

250 mL Flask. **H1000-MR-250**\$46.10

500 mL Flask. **H1000-MR-500**\$46.10

1000 mL Flask. **H1000-MR-1000**\$70.50

Mortar and Pestle Sets

Mortar and Pestle Sets

Glass, **USA**

59 mL (2 oz) Capacity. **826-2**.....\$16.35

118 mL (4 oz) Capacity. **826-3**.....\$17.20

237 mL (8 oz) Capacity. **826-4**.....\$20.40

Porcelain, outside wall and rim glazed

30 mL Capacity, 63 mm top OD. **LPC-154**\$5.00

70 mL Capacity. **600JMD070**\$12.20

80 mL Capacity, 80 mm top OD. **LPC-154-1**\$6.90

150 mL Capacity. **600JMD150**\$16.05

320 mL Capacity, 135 mm top OD. **LPC-156**\$11.50



NEW

NEW

Parafilm



Parafilm, thermoplastic sealing film, **USA**

2" x 250' (5 cm x 76 m) Roll. **HS234526A**\$34.00

4" x 125' (10 cm x 38 m) Roll. **HS234526B**\$34.90

Holds one or two 2" rolls,
or one 4" roll.



Parafilm Dispenser, ABS plastic, blue. HS234525B\$47.70

Continued

Petri Dishes



Petri Dishes, polystyrene

35 x 10 mm, vented lid, stacking ring, Simport Cs/500, sterile sleeve/20. D210-15	\$158.00
60 x 15 mm, stackable, Kord-Valmark Cs/500, sterile sleeve/20. 2901	\$112.00
100 x 15 mm, Kord-Valmark Slippable (rimless lid for automation) Standard, Cs/500, sterile sleeve/25. 2917	\$119.00
★ Stackable Standard, Cs/500, sterile sleeve/25. 2900	\$119.00
100 x 25 mm, stackable, Kord-Valmark Cs/340, sterile sleeve/17. 2906	\$132.00

pH Paper



Measure pH with 1 pH unit accuracy from 1-14.
Package contains 100 paper strips, 2 x 1/4", in a plastic vial with color chart.

pH Paper, pH 1-14, Pk/100. 600PHP100\$3.90

Pipetting



Dropper Bulbs

Amber Latex

1 mL Capacity, 43 mm long, Pk/12. 20621A	\$7.35
1 mL Capacity, 43 mm long, Pk/72. HS20622A	\$38.65
2 mL Capacity, 52 mm long, Pk/12. 20621B	\$7.35
2 mL Capacity, 52 mm long, Pk/72. HS20622B	\$40.10

Natural Rubber


1 mL Capacity, 30 mm long, black, Pk/12. 200-1B	\$3.90
2 mL Capacity, 32 mm long, red, Pk/12. 200-2B	\$6.55

Natural Rubber w/Flange, blue

1 mL Capacity, 29 mm long, Pk/24. 200-1C	\$13.30
---	---------



Natural Rubber w/Flange



glass w/rubber bulb; 75 mm (3")


Medicine (Eye) Dropper, Pk/12. 97-4301\$4.00



3-Valve



Pipette Pump



Pipette Pump II

Pipet Fillers

3-Valve, rubber, black. 207B	\$16.50
Wheel-Type, Pipette Pump™ Fits Pipets Up To 2 mL, blue. 378970000	\$38.00
Fits Pipets Up To 10 mL, green. 378980000	\$40.75
Fits Pipets Up To 25 mL, red. 378990000	\$54.05
Wheel-Type, Pipette Pump™ II, quick-release Fits Pipets Up To 2 mL, blue. 379111002	\$46.90
Fits Pipets Up To 10 mL, green. 379111010	\$46.90
Fits Pipets Up To 25 mL, red. 379111025	\$62.20



Package includes: Diamond SeroLogic portable pipette controller; two rechargeable batteries; battery charger; three PTFE membrane filters – (1) 0.20 µm & (2) 0.45 µm (one installed); desktop stand and wall mounting bracket.

Pipet Fillers

Motorized, rechargeable

Diamond SeroLogic (1-year warranty). **3380**\$319.00



Dimensions: 29 x 8.8 x 40 cm

Pipet Rack, magnetized, ABS, blue, Heathrow. HS20612M\$86.70



Dispense fluids accurately and safely with these color-coded and plugged pipets. Negative graduations provide extra volume (not available on 25 mL). RNase-/DNase-/Pyrogen-free, and noncytotoxic. Calibrated "To Deliver" (TD).

Pipets, Serological, polystyrene, UniPlast™

Nonsterile


5 mL, blue band, pk/25, Cs/375. 1755	\$110.00
10 mL, orange band, pk/25, Cs/250. 1775	\$79.65
25 mL, red band, pk/10, Cs/100. 1784	\$74.15

Sterile Pk/1 (individually paper/plastic wrapped)

1 mL, yellow band, Cs/500. 1700	\$110.00
2 mL, green band, Cs/500. 1720	\$118.00
5 mL, blue band, reverse grad., Cs/250. 1740	\$88.00
10 mL, orange band, reverse grad., Cs/200. 1760	\$80.50
25 mL, red band, reverse grad., Cs/100. 1780	\$85.00
50 mL, black band, reverse grad., Cs/100. 1820	\$184.00

Sterile Bg/25

1 mL, yellow band, Cs/1000. 1710	\$178.00
10 mL, orange band, Cs/250. 1770	\$90.00
25 mL, red band, Cs/100. 1782	\$77.15



Pipets, Transfer, LDPE (low-density polyethylene)

Graduated, Nonsterile, 15 cm long

★ #1	5 mL capacity, 1 mL graduated, 25 drops/mL	
	Bx/500. PTP-02	\$21.30
	Cs/10 Bx. PTP-02	\$196.00
	5 mL capacity, thick-walled	
	Cs/10 Bx. 188BP1200	\$164.00
	7 mL capacity, 3 mL graduated, 14 drops/mL	
	Bx/500. PTP-01	\$21.30
	Cs/10 Bx. PTP-01	\$196.00
	Graduated, Sterile Pk/1 (individually wrapped), 15 cm long	
	5 mL capacity, 1 mL graduated, Simport	
	Bx/400. P200-521S	\$49.85



These mini, air displacement pipettors provide an economical and accurate method for pipetting fixed volumes. Pipettors feature a two-step plunger and are fully autoclavable.

Pipettors, Single Channel

Fixed Volume, pocket-size, Diamond® Jr.

5 µL, red. 3322-5	\$27.75
10 µL, pink. 3322-10	\$27.75
20 µL, green. 3322-20	\$27.75
50 µL, purple. 3322-50	\$27.75
100 µL, orange. 3322-100	\$27.75

Single channel, adjustable volume pipettors are fully autoclavable.



Adjustable Volume, Diamond® Apex™ (3-year warranty)


0.5-10 µL, 0.02 µL increments, red. 3351-10	\$198.00
5-50 µL, 0.1 µL increments, yellow. 3351-50	\$198.00
10-100 µL, 0.1 µL increments, blue. 3351-100	\$198.00
100-1000 µL, 1 µL increments, orange. 3351-1000	\$198.00



Pipettor Holders, Carousel, Flip & Grip™, Heathrow

Blue. HS120026	\$92.45
Lime. HS120056	\$92.45

Pitchers



Polypropylene




Stainless Steel

Pitchers

Polypropylene, graduated (mL), tall form, tapered

1000 mL, 50 mL graduations, each. 289910000	\$7.00
2000 mL, 100 mL graduations, each. 289920000	\$8.05
3000 mL, 100 mL graduations, each. 289930000	\$13.80
Stainless Steel, graduated (mL/oz), tapered, Polar Ware	
1000 mL/32 oz, each. T1063	\$73.20

Rockers, Shakers & Mixers



Rockers, Fixed Speed, Benchmark (2-year warranty)

2-D (seesaw), 12 x 12" platform w/flat mat

BenchBlotter™. **BR1000**.....\$543.00

- Fixed speed of 12 rocks/minute
- Platform angle is $\pm 8^\circ$
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)

3-D (nutating)

Small, 10.5 x 7.5" platform w/dimpled mat

MiniMixer™. **B3D1020**.....\$449.00

Large, 12 x 12" platform w/dimpled mat

BioMixer™. **B3D1320**.....\$558.00

- Fixed speed of 24 rocks/minute
- Platform angle is $\pm 20^\circ$
- Cold room & incubator safe (4-45°C)
- Load capacity is 1.8 kg (4 lbs)

Continued



Rockers, Variable Speed

Benchmark (2-year warranty)

2-D (seesaw), 14 x 12" platform w/flat mat

BenchRocker™ 2D. **BR2000**\$830.00

Stacking Platform. **BR2000-STACK**.....\$155.00

14 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP**\$55.00

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 3.6 kg (8 lbs)

3-D (nutating), 12 x 12" platform w/dimpled and flat mats

BenchRocker™ 3D. **B3D2300**\$1,029.00

Stacking Platform. **BR1000-STACK**.....\$142.00

12 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP**\$55.00

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)



Rockers, Variable Speed

Reliable Scientific (5-year warranty), **USA**

2-D (seesaw)

Standard, 11 x 14" platform w/rubber mat

Single Platform. **55S1114**\$603.00

Double Platform, 3¼" apart. **55D1114**.....\$665.00

Large, 12 x 16" platform w/rubber mat

Single Platform. **55S1216**\$647.00

Double Platform, 3¼" apart. **55D1216**.....\$704.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



Rocker, Tube, 12 x 4.5" grooved mat, Benchmark. **M2100**\$405.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is ±30°
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs)
- 2-year warranty

Orbital Shakers



BT30



BT302

Aerate and mix with these low to high capacity orbital shakers.

Orbital Shakers, Benchmark (2-year warranty)

Low Speed, 14 x 12" platform w/flat rubber mat

Orbi-Blotter™. **BT30**\$933.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

Small, 9.5 x 8" platform w/flat rubber mat

Orbi-Shaker™ JR. **BT302**\$1,126.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min



BT3001

Large, 13 x 12" platform w/flat rubber mat

Orbi-Shaker™. **BT3001**\$1,730.00

- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L, or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR**....\$327.00


Clamps, magnetic

250 mL Flask. **H1000-MR-250**\$46.10

500 mL Flask. **H1000-MR-500**\$46.10

1000 mL Flask. **H1000-MR-1000**\$70.50

Rotating Mixers



Tube Rotator, programmable, Benchmark (2-year warranty)

Roto-Bot Mini. **R4040**\$679.00
 Roto-Bot. **R4045**\$1,019.00

- Speed adjusts from 1-80 rpm
- Tube capacity (starter set)
 - R4040: 12 x 1.5 mL, 2 x 15 mL & 2 x 50 mL
 - R4045: 30 x 1.5 mL, 8 x 15 mL & 2 x 50 mL
- Timer is programmable from 1 sec – 99 hrs, 59 min / Cont.

Vortexers



Vortex Mixers, Benchmark (2-year warranty)

Vornado™, mini vortexer, fixed 2800 rpm. **BV101**_*\$229.00

*Specify Gray, Blue or Purple cup head

BenchMixer™, single cup head, 200-3200 rpm. **BV1000**\$296.00

Flat Head Attachment, 3". **BV1000-FLAT**\$44.75

Rulers, Meter Stick & Tape Measure

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. **4461**\$0.95

White. **4460**\$0.55

Stainless Steel, 12 mm wide, 0.5 mm thick, **USA. 605-065**\$6.25

Meter Stick, lacquered wood, plain ends. **605-050**\$4.90

One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

Tape Measure, 152.4 cm/60", reinforced fiberglass. **11036**\$3.70

One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.

Sharps Containers



4402269



440321391

Sharps Containers, red polypropylene

5.4 Quart, horizontal entry lid, Prevent™, each. **4402269**\$9.90
 Dimensions: 12" x 4.75" x 11"(H)

12 Gallon, hinged vertical entry lid, each. **440321391**\$30.75
 Dimensions: 18.25" x 12.75" x 18.75"(H)

Specimen Containers

Screw-Cap, 90 mL



Screw-Cap,
480 & 960 mL



Snap-Cap

Specimen Containers

Screw-Cap, polypropylene (PP)

90 mL, sterile, Cs/400. **840175**\$150.00

480 mL (16 oz), Cs/24. **59014510**\$71.95

83 mm (3 1/4") opening

960 mL (32 oz), Cs/12. **59014511**\$59.05

111 mm (4 3/8") opening

Snap-Cap, white HDPE bottom/white LDPE lid

480 mL (16 oz), Cs/100. **271016W**\$86.20

9.9 cm opening x 11.1 cm high

1000 mL (34 oz), Cs/100. **271032W**\$129.00

11.7 cm opening x 13.7 cm high

Spot Plate



Spot Plate, glazed white porcelain, 21 mm diameter wells

12-Well, 115 x 95 x 7 mm, each. **3104187701**\$4.45

Stirring Bars & Rods



Stirring Bars, Teflon®

Centered Sphere, Saturn Spinbar®

40 mm (1½") Long. 371160002	\$31.30
50 mm (2") Long. 371160003	\$58.80

Micro, Spinbar®

White

5 x 2 mm. 371190005	\$6.40
7 x 2 mm. 371190007	\$6.40
10 x 3 mm. 371190010	\$7.10
12.7 x 3 mm. 371190127	\$7.30
15 x 1.5 mm. 371190015	\$7.90

Red

7 x 2 mm. 371210013	\$6.80
8 x 1.5 mm. 371210001	\$6.05

Octagonal, molded-on pivot ring

White

12.7 x 8 mm (½ x 5/16"). 371101256	\$6.00
15.9 x 8 mm (5/8 x 5/16"). 371100058	\$6.00
25 x 8 mm (371100138). 600MSB25	\$7.00
40 x 8 mm (371101128). 600MSB40	\$10.70
48.5 x 10 mm (371100238). 600MSB48	\$12.35
63.5 x 9.5 mm (2½ x 3/8"). 371102128	\$14.60

Red, 12.7 x 3.2 mm (½ x 1/8"), no pivot ring. **371090034**\$7.25

Stirring Bar Assortment, octagonal, Pk/6. **371600000**\$45.80

Includes two of each size: 25.4 x 8 mm (1 x 5/16"),
38.1 x 8 mm (1½ x 5/16"), and 50.8 x 8 mm (2 x 5/16").

Polygonal, molded-on pivot ring, Spinbar®

12 x 4.5 mm (7/16 x 3/16"). **371220012**\$6.35

Triangular, vanes w/cross-mounted magnets, Spinvane®

Fits 1 mL Vials. **371340000**\$12.35

Fits 3/5 mL Vials. **371350000**\$12.35



Stirring Bar Retriever, magnetic

45.7 cm (18") Long. **377730000**\$21.65



Stirring Rods, glass, annealed ends

150 x 5 mm (6 x 3/16"), Pk/12. 600GSR006	\$5.05
Rubber Policeman, fits 5 mm rod, Pk/12. 201BB	\$6.65
254 x 6 mm (10 x 1/4"), Pk/12. 600GSR010	\$7.05

Surgical Drape

Non-fenestrated towels
are sterile, latex-free
and poly-lined.



Surgical Drape, 18 x 26", Bx/50. **44025-515**\$28.30

Tape

Kraft paper with
rubber resin adhesive.



Tape, Label, Write-On™, Bel-Art

1" Core, 15 yd (540") roll
¾" Wide

Green, Pk/4. 134820075	\$33.50
White, Pk/4. 134800075	\$31.70

3" Core, 40 yd roll


½" Wide

White, each. 134630005	\$13.20
White, Pk/4. 134830075	\$70.00

¾" Wide

Green, each. 134633075	\$17.50
White, Pk/4. 134830075	\$70.00
Yellow, each. 134632075	\$17.50

134730000



Label Tape Dispensers

Multi-Roll. 134730000	\$67.60
Single-Roll. 134610000	\$149.00

Thermal Cycler



Versatile multi-format block holds (96) 0.2 mL tubes/strips, (39) 0.5 mL tubes, or (1) nonskirted 96-well PCR plate. In-situ adapter sold separately. Touchscreen control and height adjustable heated lid. 2-year warranty.



0.2 mL Tubes



96-Well Plate



0.5 mL Tubes



In-Situ Adapter

Thermal Cycler, multi-format block, Benchmark. **T5000-96** ... \$5,156.00
In-Situ Adapter, holds 4 slides. **T5000-A-IS**\$558.00

Trays



Trays

Cafeteria Style, polypropylene, 12 x 16 1/4", NSF, **USA**

Orange. **105CT1216O**\$5.95
Red. **105CT1216R**\$5.95
Black. **105CT1216B**\$5.95



Contain spills or transport lab products with this rimmed, durable, large white tray. Holds 5 gallons (18.9 L). *Not autoclavable.*


HDPE (high-density polyethylene), **USA 1053506**\$43.75
Dimensions: 26 x 18 x 3 1/2" (66 x 45.7 x 8.9 cm)



Polypropylene, autoclavable

12 x 7 3/4 x 5 1/8", Bel-Art, each. **162600000**\$31.90↓
15 x 14 x 5", Pk/6. **60081721**\$148.00
20 x 17 x 5", Pk/6. **60081722**\$156.00

Tubes & Vials



- Tubes, Centrifuge, Conical** bottom, sterile, white writing patch
RNase-/DNase-/Pyrogen-free
- 15 mL, PP (polypropylene) tube/HDPE cap, natural
Rack/25, 20 racks, Cs/500. **BP1501**\$116.00
Maximum RCF is 12,000 x g. Frosted writing patch.
- 15 mL, PS (polystyrene)
Rack/50, 10 racks, Cs/500. **6274**\$246.00
Maximum RCF is 3,000 x g.
- 50 mL, natural, polypropylene (PP) tube/HDPE cap
Bag/25, 20 bags, Cs/500. **BP5002**\$128.00
Maximum RCF is 12,000 x g. Frosted writing patch.
- Bag/25, 20 bags, Cs/500. **6297**\$200.00
Maximum RCF is 20,000 x g
- Rack/25, 20 racks, Cs/500. **BP5001**\$135.00
Maximum RCF is 12,000 x g. Frosted writing patch.

Self-Standing



- Tubes, Centrifuge, Self-Standing** (skirted)
- 30 mL
Nonsterile, PP tube/PP cap, Cs/500. **109220**\$171.00
Dimensions: 23(OD) x 90 mm
- Sterile, PP tube/HDPE cap
Cs/500, 20 bags/25. **C572-1**\$156.00
Molded-in graduations from 5 to 30 mL in 5 mL increments. Dimensions: 25.3(OD) x 111 mm

Plastic culture tubes with dual-position caps.



Tubes, Culture,

- 12 x 75 mm (5 mL)
Polystyrene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110428**\$89.50
- Polystyrene, Non sterile, Cs/1000. **110410**\$62.90
- Polypropylene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110438**\$89.70
- 13 x 100 mm (10 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1510**\$97.20
- 17 x 100 mm (16 mL)
Polypropylene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110178**\$143.00
- 16 x 125 mm (19 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1515**\$157.00
- 16 x 150 mm (23 mL)
★ #1 Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1517**\$161.00

Culture Tube Caps, polypropylene Up to 16 mm Tubes



- Blue, Cs/1000. **118152B**\$99.65
- Green, Cs/1000. **118152G**\$99.65
- Red, Cs/1000. **118152R**\$99.65
- Yellow, Cs/1000. **118152Y**\$99.65



Black



RNase/DNase-free

Tubes, Microcentrifuge, Economy, polypropylene

- 0.5 mL, 0.1 mL graduations, natural, Pk/1000. **HS4422**\$38.45
- 1.5 mL, 0.25 mL graduations, Heathrow
Natural, max RCF is 16,400 x g
Pk/500. **HS4323**\$18.70
Pk/1000. **HS4325**\$35.35
- Black, max RCF is 18,000 x g, Pk/500. **HS4323K**\$20.80
- 5 mL, natural, max RCF is 25,000 x g, Benchmark
Nonsterile, 4 bags/50, Pk/200. **C1005T5**\$65.40
Sterile, 10 bags/20, Pk/200. **C1005T5ST**\$71.45

Continued



Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable. Free of RNase/DNase.

Tubes, Microcentrifuge, Premium, PP, G-Tube™, Bio Plas, **USA**

Amber color protects light sensitive specimens.

0.6 mL, Pk/500.....	\$27.75
Amber. 4048 Lavender. 4047 Natural. 4040 Mixed (blue, green, orange, pink, & yellow). 4049	\$33.45
1.5 mL, 0.5 mL graduations, frosted area on side, Pk/500	\$27.75
Amber. 4038 Blue. 4033 Lavender. 4037 Natural. 4030 Mixed (blue, green, orange, pink, & yellow). 4039	\$33.45
2.0 mL, 0.5 mL graduations, frosted area on side Natural, Pk/500. 4050	\$29.80



Tubes, PCR, 0.2 mL thin-wall, polypropylene, Simport

Individual, flat caps, natural, Pk/1000. **T325-2VN**\$50.20



Vials, Cryogenic, round bottom, internal thread

O-Ring Seal, 4 mL, Cs/1000. **T301-4**\$297.00
Washer Seal, 4 mL, Cs/1000. **T311-4**\$297.00



Vinyl grip, **USA**

Vial Decapping Pliers

13 mm Crimp Seals. **D-13**\$72.00
20 mm Crimp Seals. **D-20**\$72.00

Tube Storage

Cryogenic Storage Boxes



SB2CC-64 Cryogenic Tube Storage Box

Tube Storage Boxes, Cardboard, vented

2" Box Only, 5¼ x 5¼ x 2", each. SB2CC	\$1.80
2" Box w/dividers, white	
100-Cell (10 x 10), Pk/12. SB2CC-100	\$28.45
★ 81-Cell (9 x 9), Pk/12. SB2CC-81	\$28.20
64-Cell (8 x 8), Pk/12. SB2CC-64	\$28.20
2" Mini box w/dividers, 3 x 3 x 2"	
25-Cell (5 x 5), Pk/12. MB2C-25	\$30.65
Holds 13 mm vials and 4.0-5.0 mL cryovials.	
3" Box Only, 5¼ x 5¼ x 3", each. SB3C	\$1.80
3" Box w/dividers	
100-Cell (10 x 10), Pk/12. SB3C-100	\$28.45
81-Cell (9 x 9), Pk/12. SB3C-81	\$28.20
64-Cell (8 x 8), Pk/12. SB3C-64	\$28.20
Dividers only, Fit 2 & 3" Boxes	
100-Cell (10 x 10), Pk/12. CD-100	\$11.85
Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.	
81-Cell (9 x 9), Pk/12. CD-81	\$11.55
Holds 13 mm vials and 4.0-5.0 mL cryovials.	
64-Cell (8 x 8), Pk/12. CD-64	\$11.55
Holds 14 mm vials.	
49-Cell (7 x 7), Pk/12. CD-49	\$11.85
Holds 14 mm vials.	
25-Cell (5 x 5), 25 mm cell, Pk/12. CD-25	\$11.55
16-Cell (4 x 4), 31 mm cell, Pk/12. CD-16	\$11.55
5" Box Only, 5¼ x 5¼ x 5", each. SB4C	\$3.30
5" Box w/dividers	
64-Cell (8 x 8), Pk/12. SB4C-64	\$43.50
36-Cell (6 x 6), 15 mL tubes, Pk/12. SB4C-36	\$43.50
16-Cell (4 x 4), 50 mL tubes, Pk/12. SB4C-16	\$43.50
Dividers only, Fit 5" Box	
64-Cell (8 x 8), Pk/12. CD4-64	\$11.55
36-Cell (6 x 6), Pk/12. CD4-36	\$11.55
16-Cell (4 x 4), Pk/12. CD4-16	\$11.55



Holds 1.0-2.0 mL tubes and cryovials. Features a locking hinged lid, and base with molded alphanumeric coordinates. Dimensions: 12.8 x 12.8 x 5.4 cm.

Tube Storage Boxes, Polypropylene

81-Well (9x9), 2 blue/2 yellow boxes, Pk/4. **60066501**\$33.30


Tube Racks



Tube Racks, Grip, polypropylene, No-Wire™, Bel-Art
18-20 mm Tubes, 40-place, each. **HS290002**\$35.30



Tube Racks, Interlocking, polypropylene, Heathrow
4-Way, mixed colors, Pk/5. **HS29022G**\$46.15
Holds (4) 50 mL tubes, (12) 15 mL tubes,
(32) 1.5 mL tubes, or (32) 0.5 mL tubes.
4-Way Rotatable, Rota-Rack™, four connected modules
Large Tubes, 22.4 x 11 x 11 cm, each. **HS29040A**\$47.80
Each module holds (2) 50 mL tubes, (8) 15 mL tubes,
(6) 20 mm tubes, or (10) 10-13 mm tubes. Mixed colors.



Tube Racks, Mega Rack, polypropylene, fully-assembled
10-13 mm Tubes
216-Place, 33.7 x 26.7 x 9.1 cm(H)
Blue, each. **HS24322A**\$50.30
White, each. **HS24322E**\$50.30
432-Place, 67.4 x 26.7 x 9.1 cm(H)
Blue, each. **HS120113**\$87.70
White, each. **HS120112**\$87.70



80-Well



96-Well Reversible

Tube Racks, Microtube, polypropylene, Heathrow
80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm
Holds 1.5-2.0 mL tubes in conical wells.
Blue, Pk/5. **HS29025B**\$27.90
Pink, Pk/5. **HS29025D**\$27.90
Mixed, Pk/5. **HS29025G**\$27.90
96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.
Pink, Pk/5. **HS2345E**\$56.85


Available in blue,
green, orange,
purple & yellow.



Tube Rack, PCR, polypropylene, mixed colors, Heathrow
Reversible, Pk/5. **HS2344A**\$60.60
Holds (168) 0.2, (40) 0.5 and (12) 1.5 mL tubes.



Tube Racks, Peg, 14-17 mm Tubes, 50-place, polypropylene
White, Pk/2. **HS24311A**\$30.15




Half-Size



Full-Size

Tube Racks, Teardrop, OneRack™, POM, Heathrow
Half-Size
16 mm tubes, 36-place, 12.7 x 12.7 x 7 cm(H)
Blue, Pk/8. **HS27502B**\$115.00
Full-Size
13 mm tubes, 72-place, 20.2 x 10.4 x 5.9 cm(H)
Blue, Pk/8. **HS27511B**\$151.00
Magenta, Pk/8. **HS27511D**\$151.00
16 mm tubes, 72-place, 25 x 12.7 x 7 cm(H)
Blue, Pk/8. **HS27512B**\$151.00
Magenta, Pk/8. **HS27512D**\$151.00



10-13 mm Tubes, 108-place (9 x 12), each. **345120759**\$32.00
13-16 mm Tubes
108-Place (9 x 12), each. **345120763**\$32.00
15-16 mm Tubes
180-Place (9 x 20), each. **187860770**\$61.60↓
Dimensions: 38.6 x 18.3 x 6.4 cm(H)

Hexagonal



Standard

★ Small, 4.6 x 4.6 x 0.8 cm(D), Bx/500. HS1420A	\$24.55
Medium, 8.5 x 8.5 x 2.4 cm(D), Bx/500. HS1420B	\$39.90
Large, 14 x 14 x 2.2 cm(D), Bx/500. HS1420C	\$65.40

Small, 3.8 cm top dia., Bx/500. HS14251A	\$18.05
Medium, 6.4 cm top dia., Bx/500. HS14251B	\$31.15
Large, 12 x 14 cm top OD, Bx/500. HS14251C	\$81.50
(25433-104)	
X-Large, 14 x 16 cm top OD, Bx/500. HS14251D	\$104.00
(25433-106)	

Small, 5.8 x 4.3 x 1.3 cm(D), Bx/500. HS1419A	\$40.85
Medium, 13.2 x 8.3 x 2.6 cm(D), Bx/500. HS1419B	\$80.05
Large, 18.3 x 10.8 x 2.6 cm(D), Bx/500. HS1419C	\$116.00

Extra low-lint, nonabrasive, highly absorbent, 1-ply wipes are perfect for microscope maintenance.

4½ x 8½" Wipes. 34155	\$4.00	214.00
------------------------------------	--------	--------



WypAll® L40, limited use. 43005701	\$7.45	117.00
--	--------	--------

WypAll® X60, pk/100, cs/9 pk, Lot/9 cs. **43034770**\$2,610.00
Measures 11 x 23"