


Scissors, Operating (cont.)


- Curved, S/B (sharp/blunt) tips
- 11.5 cm (4½"), Economy. **310-101** \$7.65
 - 14 cm (5½"), Economy. **12-125** \$8.55
 - 14.5 cm (5¾"), Standard. **10-560** \$15.40 ↓
 - 14 cm (5½"), German. **100-560** \$37.90
 - 16.5 cm (6½"), Economy. **12-135** \$10.75
 - 16.5 cm (6½"), Standard. **10-610** \$25.25
 - 16.5 cm (6½"), 60 mm blades, Roboz. **RS-6849** \$81.20
- Curved, B/B (blunt/blunt) tips
- 11.5 cm (4½"), Economy. **310-102** \$7.85
 - 14 cm (5½"), Economy. **12-185** \$8.55
 - 14 cm (5½"), Standard. **10-568** \$21.80
 - 16.5 cm (6½"), Economy. **12-195** \$10.75
 - 16.5 cm (6½"), Standard. **10-618** \$25.25



- Scissors, Stevens Tenotomy, 11.5 cm (4½"), stainless steel**
 Straight, blunt tips, Standard. **13-520** \$24.75



- Scissors, Strabismus, stainless steel**
- Straight
- 10 cm (4"), Economy. **26-105** \$10.60
 - 11.5 cm (4½"), Standard. **13-524** \$23.65
- Curved
- 11 cm (4¼"), Economy. **26-110** \$10.60
 - 11.5 cm (4½"), Standard. **13-526** \$24.15



- Scissors, Bandage, stainless steel, economy**
- Knowles, straight, sharp probe tips, 14 cm (5½"). **1070** \$6.25
 - Lister, angled, blunt probe tips, 14 cm (5½"). **1050** \$5.00
- Can be used as enterotomy scissors.



Cut through everything from wire to bone with these durable scissors. Bottom blade is serrated and has a footplate for safety. Blades are stainless steel & handles are plastic. Autoclavable to 120° C.

- Utility Shears, 18 cm (7¼"), economy. 310-341** \$5.70

Scrapers

Remove paraffin, ice, etc., safely from your laboratory equipment with these plastic scrapers.



- Scrapers, plastic**
- 5 cm (2") Wide Blade. **311-245** \$3.70
 - 10 cm (4") Wide Blade. **311-250** \$4.70
 - 15 cm (6") Wide Blade. **311-255** \$5.70

Spatulas & Scoops



- Spatulas, Hardwood Handle, stainless steel blade**
- 10 cm (4") Blade. **14-365B** \$5.40
 - 15 cm (6") Blade. **14-365D** \$7.05
 - 20 cm (8") Blade. **20-497** \$6.05



- Spatulas, Heavy-Duty, stainless steel**
 Round End/Square End, both flat
- 30.5 cm (12"). **320-425** \$16.10
 - 46 cm (18"). **320-430** \$21.50



- Spatulas, Micro, stainless steel**
- Square End/Spoon End
- 19 cm (7½"). **320-411** \$4.50
 - 23 cm (9"). **320-414** \$4.65
- Round End/Tapered End, 18 cm (7") long. **320-370** \$5.00
 Round end: 8 x 51 mm; Tapered end: 6 x 51 mm
- V-Shaped End/Round End, 16 cm (6¼") long. **HS15906** \$7.45



- V-Shaped (Reagent Digger), 16.5 cm (6½"). **600LSCSS6** \$2.60
 Handle, wooden, 13 cm long, each. **600LSHND1** \$3.00