

Teasing Needles, Probes & Pin Vises

97-1302



Teasing Needles

Angular, wooden handle, 14.5 cm (5¾"), Pk/12. **97-1302**..... \$4.95
 Straight, plastic handle, 14.5 cm (5¾"), Pk/12. **97-1313** \$4.65



Stainless Steel, straight, 14 cm. **M28** \$5.00

Blunt tip is angled and tapered. Hexagonal handle is 6 mm wide.



Probes, Light Mall, 15.5 cm (6")

Nickel Plated Steel. **15-800** \$2.90
 (5800)
 Stainless Steel. **15-800SS** \$3.80




Probes, Bowman Lacrimal, double-ended

Economy, sterling silver, 12.7 cm (5") long
 Size 0000-000 (~0.4 – 0.5 mm). **24-706**..... \$11.65
 Size 00-0 (~0.6 – 0.7 mm). **24-708**..... \$11.65
 Size 1-2 (~0.9 – 1.0 mm). **24-710**..... \$11.65
 Size 3-4 (~1.1 – 1.2 mm). **24-712**..... \$11.65



Probe/Seekers, stainless steel

Double-Ended, 15.5 cm. **320-216A** \$4.90
 One end tapered to a point; one end curved with a rounded tip.
 Single End, 12.5 cm (5"). **320-217A**..... \$4.75
 Curved, fine point tip with hexagonal handle.



Probe with Bulb Tips, 25 cm (10"), stainless, German. 541-807 \$12.00 ↓



Probe with Eye & Bulb Tip, flexible, 15.5 cm (6"). 23-496..... \$2.25



PIN-220.00




PIN-212.00

Hold needles, probes, wires, drill bits, etc.

Pin Vise, Swivel-Head, 0-3.3 mm opening. PIN-220.00 \$8.35
 Two double-ended collets are provided; one in screw chuck, one inside handle.

Pin Vise, Mini, 0-0.8 mm opening. PIN-212.00 \$5.20
 blackened tool steel

Trocar



Trocar, 4-in-1 set, 13 cm, stainless steel. 1596..... \$17.50

Accessories/Maintenance/Storage



Cutting Boards, white high-density polyethylene
 10 x 6 x ½" (25.4 x 15.2 x 1.3 cm). **1054319** \$6.95
 18 x 12 x ½" (45.7 x 30.5 x 1.3 cm). **105CB1218** \$16.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1½ x 3" opening. DEM-345.00..... \$49.95



Dissecting Cases

6-Loop with snap fastener, black leatherette. **1257370** \$1.30
 12-Loop with snap fastener, black leatherette. **1257371** \$1.40
 13-Loop with zipper, brown leatherette. **36900** \$7.80
 Features red velveteen interior, protector flap, and ID pouch.