

Dissecting Kits

*Inquire about customizing a kit for your surgical procedure.
Kit components are discounted 10%.*

Basic kit. 38501.....\$5.30

Contents: Dissecting Case, 6-loop (1257370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

Mouse kit. 38551.....\$17.85

Contents: Dissecting Case, 6-loop (1257370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

Dissecting Pans, stainless steel (without wax)

20⁷/₈ x 12¹³/₁₆ x 2¹/₂". **AM20029**.....\$21.90
25³/₄ x 17³/₄ x 2¹/₄". **AM68357**.....\$64.90



Insect Pins, 38 mm (1¹/₂") long, nylon head, Czech Republic

Black Enameled, Pk/100.....\$7.00

#000, 0.25 mm dia. **355425** #1, 0.40 mm dia. **355440**
#00, 0.30 mm dia. **355430** #2, 0.45 mm dia. **355445**
#0, 0.35 mm dia. **355435** #3, 0.50 mm dia. **355450**

Stainless Steel, Pk/100.....\$8.75

#000, 0.25 mm dia. **355525** #1, 0.40 mm dia. **355540**
#00, 0.30 mm dia. **355530** #2, 0.45 mm dia. **355545**
#0, 0.35 mm dia. **355535** #3, 0.50 mm dia. **355550**



Instrument Cup, silicone, 3¹/₂" opening, 2¹/₂" deep. **150-1005**.....\$57.00



Instrument Storage Boxes, 10 x 4 x 2", 18-8 stainless steel

Small, 7 x 3 x 1¹/₂". **35-406**.....\$26.60
Large, 10 x 4 x 2". **35-408**.....\$35.00

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. **4461**.....\$0.95

White. **4460**.....\$0.55

Stainless Steel, 12 mm wide, **USA. 605-065**.....\$6.25



Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

Sterilizer, Hot Bead, Germinator™ 500, USA. 5-1450.....\$426.00



Glass Beads (replacement), 1.5 mm, Bg/12 oz. **5-1455**.....\$30.00



Safely incinerate microorganisms in 5-7 seconds at 815°C (1500°F) within this ceramic, stainless steel-encased chamber. Compact and flameless operation allows this unit to be used in anaerobic and aerobic chambers. Ideal for sterilizing loops and needles, glass tube/pipette mouths, and various glass and metal instruments. 1-year warranty.

Sterilizer, Infrared Heat, BactiZapper™. B1000.....\$476.00

6100A



Sterilizing Trays, plastic, autoclavable, PST, **USA**

B/L/M (base, lid & mat) included.

6 x 2¹/₂ x 3/4" (15.2 x 6.4 x 1.9 cm). **2600A**.....\$23.00
(500252)



6¹/₂ x 4 x 3/4" (16.2 x 10.2 x 1.9 cm). **4600A**.....\$33.00

10 x 6 x 3/4" (25.4 x 15.2 x 1.9 cm). **6100A**.....\$67.00
(500254)



Instrument Trays

Polypropylene, deep w/flat cover, autoclavable

8¹/₂ x 6 x 2³/₄" (22 x 15 x 7 cm), Pk/2. **60081736**.....\$25.15

Continued