

SLIDE & CASSETTE PROCESSING

Microscope Slides

Slides, general, 75 x 25 x 1 mm

Plain

Globe, ground edges, 90° corners
Cs/20 Bx, bx/72. **1301** \$152.00

Frosted (etched)

Globe, beveled edges, clipped corners
Cs/20 Bx, bx/72. **1334** \$204.00

White End


Globe, beveled edges, clipped corners
Cs/20 Bx, bx/72. **1334W** \$238.00

Superfrost®, "SPECIMEN". **4951-001T** \$327.00
(12-550-143, 48311-600)

Superfrost® Plus



Concavity, 2-Well



Slides, specialty

Charged (silanated), 75 x 25 x 1 mm
Recommended for adhering tissue sections
and cell preparations.

Diamond™, white end, ground edges, 90° corners
Cs/20 Bx, bx/72. **1358W** \$475.00

Superfrost® Plus, white end, Erie (EpreDia)
Cs/20 Bx, bx/72. **4951PLUS-001** \$682.00 ↓
Preferred slide for Ventana instrumentation.
(12-550-15, 48311-703)

Concavity

75 x 25 x 1.3 mm, 16 mm dia., 0.8 mm deep, economy
1-Well, Bx/12. **600CS3X11** \$7.75

2-Well, Bx/12. **600CS3X12** \$10.75

76 x 26 x 3.2 mm, 18 mm dia., 1.75 mm deep, Erie (EpreDia)
1-Well, Bx/12. **1526** \$97.35

Fused Quartz, NSG Type OZ®, 3" x 1" x 1 mm, **USA**
QMS3X1X1, each \$20.00

Fused Silica, Corning 7980, 3" x 1" x 1 mm, **USA**
FSMS3X1X1, each \$39.90

Coverslips (Cover Glasses)

Coverslips, premium, borosilicate glass, 1 oz/pk, Cs/10 oz

#0 Thickness (0.08-0.12 mm)

Circles, Germany*

25 mm, ~222 pcs/oz, Cs/10 oz. **0-25** \$214.00



#1 Thickness (0.13-0.17 mm)

Squares & Rectangles, Erie (EpreDia)

22 x 22 mm, ~162 pcs/oz, Cs/10 oz. **1-2222** \$70.70

22 x 30 mm, ~119 pcs/oz, Cs/10 oz. **1-2230** \$67.00

22 x 50 mm, ~71 pcs/oz, Cs/10 oz. **1-2250** \$67.00

24 x 50 mm, ~65 pcs/oz, Cs/10 oz. **1-2450** \$67.00

24 x 60 mm, ~54 pcs/oz, Cs/10 oz. **1-2460** \$67.00



#1 Thickness (0.13-0.17 mm)

Circles, Germany*

12 mm, ~663 pcs/oz, Cs/5 oz. **1-12** \$249.00

18 mm, ~296 pcs/oz, Cs/5 oz. **1-18** \$150.00



#1.5 Thickness (0.16-0.19 mm). *Recommended for most users*

Squares & Rectangles, Erie (EpreDia)

18 x 18 mm, ~209 pcs/oz, Cs/10 oz. **1.5-1818** \$71.70

22 x 22 mm, ~140 pcs/oz, Cs/10 oz. **1.5-2222** \$73.35

22 x 50 mm, ~62 pcs/oz, Cs/10 oz. **1.5-2250** \$67.00

24 x 40 mm, ~71 pcs/oz, Cs/10 oz. **1.5-2440** \$67.00

24 x 50 mm, ~56 pcs/oz, Cs/10 oz. **1.5-2450** \$63.80

24 x 60 mm, ~47 pcs/oz, Cs/10 oz. **1.5-2460** \$67.00



#1.5 Thickness (0.16-0.19 mm). *Recommended for most users*

Circles, Germany*

12 mm, ~552 pcs/oz. **1.5-12** cs/5 oz \$198.00 cs/10 oz —

18 mm, ~245 pcs/oz. **1.5-18** \$250.00 —

22 mm, ~164 pcs/oz. **1.5-22** — \$252.00

25 mm, ~127 pcs/oz. **1.5-25** — \$252.00



*German glass is preferred for cell growth/adhesion

Coverslips, specialty

Fused Quartz, **USA**

0.15-0.25 mm Thick, NSG Type OZ®

1" Square, each. **QCS1X1** \$19.00

1" Circle, each. **QCS1X0** \$19.00

0.16 ± 0.01 mm Thick, A1

25 mm Square, each. **R525025** \$32.40

25 mm Circle, each. **R525000** \$23.40

Fused Silica, **USA**

0.15-0.25 mm Thick, Corning 7980

1" Square, each. **FSCS1X1** \$21.50

1" Circle, each. **FSCS1X0** \$21.50

0.16 ± 0.01 mm Thick, S1-UV

25 mm Square, each. **R425025** \$32.00

25 mm Circle, each. **R425000** \$28.80