

## Accessories

### Slide & Coverslip Handling

#### Forceps


##### Spade Tips

- Angled (Kuehne), 11 cm (4¼")  
 Stainless Steel, economy. **33-1041** ..... \$4.15  
 Teflon® Coated, Bel-Art. **379420000** ..... \$35.60  
 Straight, stainless steel, 10.5 cm, economy. **57.926** ..... \$8.45



##### Square Point Tips, straight

- Stainless Steel, 12 cm (4¾"), economy. **1244** ..... \$2.25  
 Teflon® Coated, 11 cm (4¼"), Bel-Art. **379430000** ..... \$39.40



Minimize contamination during slide & coverslip handling. Will not scratch; corrosion and heat resistant.

- Hemostat**, Polypropylene, 12.5 cm. **HEM-300.00** ..... \$1.95

### Immersion Oil & Bottles



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities.  
**USA**

#### Immersion Oil

##### Type A, low viscosity (150 cSt). Short W.D. optics

- 7.5 mL in LDPE dispensing bottle. **18101** ..... \$10.00  
 20 mL in HDPE dispensing bottle. **18102** ..... \$15.35  
 480 mL in amber bottle. **18104** ..... \$110.00

##### Type B, high viscosity (1250 cSt). Long W.D. optics

- 7.5 mL in LDPE dispensing bottle. **18105** ..... \$10.00  
 20 mL in HDPE dispensing bottle. **18106** ..... \$15.35  
 480 mL in amber bottle. **18108** ..... \$110.00

##### Type NVH, very high viscosity (21,000 cSt)

- 480 mL in amber bottle. **18113** ..... \$67.00


##### Type HF, low fluorescence (700 cSt)

- 480 mL in amber bottle. **18118** ..... \$169.00

##### Type FF, fluorescence free (170 cSt)

- 480 mL in amber bottle. **18124** ..... \$163.00

A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



- Bottle, Dispensing**, HDPE, 20 mL, each. **18145** ..... \$8.50

### Slide Storage

#### PLASTIC SLIDE BOXES



25-Slide Squeeze-Release Cover



100-Slide Steel-Hinged Cover

50-Slide Hinged Cover



#### Slide Boxes, hold 75 x 25 mm (3 x 1") slides

- 25-Slide, squeeze-release cover, cork bottom, 9.5 x 8.5 x 3 cm  
 White, ABS plastic. **97-0025** ..... \$2.50  
 50-Slide, hinged cover, cork bottom, 20 x 9.5 x 3.5 cm  
 Blue, ABS plastic. **97-0056** ..... \$6.80  
 Red, ABS plastic. **97-0058** ..... \$6.80  
 100-Slide, steel-hinged cover, cork bottom  
 Blue, ABS plastic. **97-0126** ..... \$8.90  
 Red, ABS plastic. **97-0128** ..... \$8.90  
 White, ABS plastic. **97-0125** ..... \$8.90



Store 4500 75 x 25 mm slides in 6 divided drawers. Each cabinet comes with a sample of 6 compressor follow blocks and 6 drawer labels. Safely stack up to 12 cabinet units on a single base.  
 18¾"(D) x 15¼"(W) x 5"(H)

#### Slide Cabinets, steel

- 6-Drawer, stackable  
 Green. **513500G** ..... \$255.00  
 Tan. **513500T** ..... \$255.00  
 White. **513500W** ..... \$255.00  
 Memorandum Markers, paperboard, Pk/100. **C108** ..... \$12.75  
 Compressor Follow Block, rubber, each. **C109** ..... \$1.90  
 Base, raises cabinets 4¼" off floor  
 Green. **513503G** ..... \$169.00  
 Tan. **513503T** ..... \$169.00  
 White. **513503W** ..... \$169.00



Cardboard Folder



Plastic Folder

#### Slide Folders, 20-slide, hold 75 x 25 mm (3 x 1") slides

- Cardboard, black tape, economy. **HS9902** ..... \$6.35 (12-587-10)  
 Plastic, gray tray w/2 clear doors. **M750-20GY** ..... \$11.60  
 Transparent covers hold visible slides securely in place, but can fold flat against back of tray when opened. Slides are slightly elevated for easy retrieval.