


Slide Holders, perfect for slide drying or staining
 52-Slide (approx.), 12 x 4", HDPE, **USA. L500-C** \$30.95
 91-Slide (approx.), 12 x 7", HDPE, **USA. L500-D** \$41.35




5-Slide mailer can also be used as a staining jar.

Slide Mailer, 5-Slide, plastic, end-opening, Cs/200. **62706** \$87.30



SL-3-264

SL-3-182



Slide Sleeves, 10-gauge vinyl
 4 x 5½" Card Holder Plus 4 Slide Pockets. **SL-3-264** \$0.45
 8½ x 11" Page Holds 18 Slides. **SL-3-182** \$0.85

Aluminum



Plastic

Slide Trays, 20-slide
 Aluminum, Eberbach, **USA. E2870** \$22.00
 Plastic, white. **445810000** \$17.30

Slide & Cassette Pens, Marker & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker.
 Pilot pens are ideal for cytologists and pathologists.



Pens, Pilot®, permanent, extra fine point
 Black, Bx/12. **44102** \$17.45
 ★ Blue, Bx/12. **44103** \$17.45
 Green, Bx/12. **44105** \$17.45
 Red, Bx/12. **44104** \$17.45



The 0.75 mm tip is ideal for use on small surfaces, e.g., histology cassettes, microscope slides and cryogenic vials. The permanent ink is waterproof and resistant to most laboratory chemicals and solvents.

Marker, Diamond PerfectMark™, black, Bx/10. **1000** \$33.30



Scribes, retractable
 Carbide-Tipped (sharpenable). **RSC** \$24.30
 Carbide Tip (replacement). **RCT-1** \$8.45
 Diamond-Tipped. **RSD** \$28.35
 Diamond Tip (replacement). **RDT-1** \$14.35

Slide Warmers



Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.

Slide Warmers
 Small, 26 x 18 cm (10 x 7"). **XH-2002** \$371.00
 Large, 63 x 20 cm (25 x 8"). **XH-2001** \$487.00