


## Staining



Ten slide compartments, separated by 1/2" of space, are additionally isolated by eight barrier dividers when the lid is closed. Individual slides are suspended 1/2" above the chamber floor by four pedestal posts and four corner posts.

### Staining/Moisture Chambers

10-Slide, polystyrene, <b>USA</b>	<u>ea</u>	<u>pk/3</u>
Clear. <b>240-9000-010</b> .....	\$54.95	154.00
Black. <b>240-9020-Z10</b> .....	\$56.10	157.00



20-Slide, Simport

★ Clear Lid (PETG), black base (ABS), each. <b>M920-1</b> .....	\$195.00
★ Black Lid (ABS), black base (ABS), each. <b>M920-2</b> .....	\$195.00



### Staining Dishes


Glass (flint)

20-Slide, holds slides horizontally in metal rack	
Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H)	
Set (dish, rack & lid). <b>SDS-20</b> .....	\$31.85
Dish & Lid. <b>SD-20</b> .....	\$19.40
Rack only, stainless steel. <b>SR-20</b> .....	\$16.00
30-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H)	
Set (dish, rack & lid). <b>SDS-30</b> .....	\$43.50
Dish & Lid. <b>SD-30</b> .....	\$21.70
Rack only, stainless steel. <b>SR-30</b> .....	\$26.65
60-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H)	
Set (dish, rack & lid). <b>SDS-60</b> .....	\$60.80
Dish & Lid. <b>SD-60</b> .....	\$33.20
Rack only, stainless steel. <b>SR-60</b> .....	\$34.36



Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

Plastic (polyoxymethylene [POM]), holds 25 slides	
Set (dish, rack & lid). <b>0990-99</b> .....	\$50.90



### Staining Jars, hold slides (up to 76 x 26 mm) vertically

Coplin, holds 10 slides back-to-back	<u>ea</u>	<u>cs/12</u>
Glass w/Ground Glass Cover, <b>USA. 673</b> .....	\$13.45	149.00
Glass Cover, Fits 673 Jar. <b>673-LID</b> .....	\$4.95	—
Glass w/Screw-On Phenolic Cap, <b>USA. 674</b> .....	\$13.45	149.00
Phenolic Cap, Fits 674 Jar, black. <b>674-LID</b> .....	\$5.05	—
PP w/Screw-On Cap. <b>60062101</b> .....	\$5.05	49.20
Polypropylene jar is unbreakable, microwave proof, stain resistant, and chemically inert. Cap is domed.		
Wide Top, Glass w/Square Lid, <b>USA. 675</b> .....	\$18.00	197.00
Holds 16 slides back-to-back		
Glass Lid, Fits 675 Jar. <b>675-LID</b> .....	\$5.50	—



Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

### Wide Top, Plastic (POM)

12-Slide, EasyDip™, Simport	
Jar only, Blue, Cs/6. <b>M900-12B</b> .....	\$67.65
Jar only, Green, Cs/6. <b>M900-12G</b> .....	\$67.65
Jar only, Pink, Cs/6. <b>M900-12P</b> .....	\$67.65
Jar only, White, Cs/6. <b>M900-12W</b> .....	\$67.65
Jar only, Yellow, Cs/6. <b>M900-12Y</b> .....	\$67.65
Rack/Lid only, gray, Cs/6. <b>M905-12DGY</b> .....	\$74.90



Anodized aluminum staining jar holder

(Jars sold separately)

Staining Jar Holder, 6-Place. <b>M906</b> .....	\$82.90
---	---------

Continued