



Simport



Simport  
Macrosette

### Tissue Cassettes, Attached Lids, slotted pores

Histosette™ I, snap-apart lids, Pk/1500.....\$154.00  
30° front labeling area.

Blue. **M490-6** Pink. **M490-3**  
Green. **M490-4** White. **M490-2**

Histosette™ II, integral lids, Pk/1500.....\$204.00  
45° front labeling area.

Blue. **M492-6** Pink. **M492-3**  
Green. **M492-4** White. **M492-2**

Uniset™, hinged lids, Pk/1500 .....\$165.00  
35° front labeling area. (Alternative to Sakura's Uni-Cassette®)

Blue. **M505-6** Pink. **M505-3**  
Green. **M505-4** White. **M505-2**

Macrosette™, 2X normal depth, white, Pk/750. **M512** .....\$221.00  
Dimensions: 40.1 x 28.5 x 13 mm(H)

## Cassette Storage



Matchbox Style

**Box, Matchbox Style (Pill Box), cardboard, white**

2 1/8 x 1 1/8 x 5/8", Pk/72. **20047** ..... \$35.45

**Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. 8342** ..... \$121.00

Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks  
in 4 drawers, each divided into 12 rows.  
Cabinets can be stacked up to 8 high.  
Dimensions: 49.5 x 38 x 21 cm  
(19 1/2 x 15 x 8 1/4")



**Tissue Cassette Cabinet, Steel, 4-drawer. 485400** ..... \$386.00

Store up to 2100 cassettes in 6 drawers.  
Each drawer includes sponge block  
spacers and an identification label.  
Cabinets can be stacked up to 4 high.  
Dimensions: 44.5 x 24 x 39.5 cm(H)  
(17 1/2 x 9 1/2 x 15 1/2")



**Tissue Cassette Cabinet, ABS. 97-0600** ..... \$204.00



Simport


Store 250 cassettes or 165 embedding  
rings in each stackable drawer.  
High-impact polystyrene resists  
temperatures from -80°C to +80°C.  
Drawer dimensions: 40.5 x 23 x 5.1 cm.

**Tissue Cassette Drawers, Polystyrene**

Pk/2. **M495-2** ..... \$77.55

Pk/6. **M495-6** ..... \$176.00

## Brushes



LB-30

### Brushes

Flat Bristles (paint style)

Camel Hair, 1" wide, 6 1/2" total length. **BR665B (LB-37)** ..... \$6.65

Round Bristles, camel hair

3/4" Diameter, 8 1/2" total length. **BR661 (LB-30)** ..... \$11.15