

## Centrifuges



GCH-24

**Centrifuge, Hematocrit, 24-place. GCH-24**..... \$2,439.00  
 24-place capillary rotor and reading disk lid are included.



C1012



C1008-C

### Centrifuges, Mini

5,500 RPM (2000 x g), 12-place, MyFuge™12. **C1012** ..... \$448.00  
 (1-year warranty)  
 6,000 RPM (2000 x g), 8-place, MyFuge™Mini. **C1008-\_\_** \* \$333.00  
 \*Specify clear (-C), green (-G), or purple (-P) for lid color.  
 (1-year warranty)

## Clamps & Supports

Rectangular  
Steel Base



Triangular  
Cast-Iron Base




### Support (Ring) Stands, finished metal base/plated support rod

**Rectangular Steel Base**  
 6 x 4" base, 5/16 x 18" rod. **3104SSS66** ..... \$15.45  
 8 x 5" base, 3/8 x 20" rod. **3104SSS67** ..... \$19.70  
 9 x 6" base, 1/2 x 24" rod. **3104SSS68** ..... \$25.00

**Rectangular Cast-Iron Base**  
 Provides greater stability than steel base.  
 6 x 4" base, 5/16 x 18" rod. **3104SSC14** ..... \$18.80  
 8 x 5" base, 3/8 x 20" rod. **3104SSC15** ..... \$20.80  
 9 x 6" base, 1/2 x 24" rod. **3104SSC16** ..... \$30.40

**Triangular Cast-Iron Base**  
 5" Legs, 1/2 x 24" rod. **1257G2324** ..... \$25.00  
 6" Legs, 1/2 x 36" rod. **1257G23** ..... \$30.65



600COBR3FR




600THCL01



600CLBURDD



97-4023



1257G49

### Support Clamps

Burette, plastic coated jaws, 360° rotatable. **600COBR3FR** .... \$8.70  
 Burette, Double, deluxe. **600CLBURDD** ..... \$22.95  
 Plastic covered spring jaws on a metal frame securely grip  
 burettes up to 30 mm OD. Fits rods up to 3/4" (19 mm) OD.  
 Hook Collar, nickel plated. **1257G49** ..... \$6.45  
 Fits rods up to 1/2" (13 mm) in diameter.  
 Right-Angle, aluminum. **97-4023** ..... \$8.60  
 Fits rods up to 5/8" in diameter.  
 Thermometer, aluminum, rotates 360°. **600THCL01** ..... \$11.10  
 Fits rods up to 1/2" (13 mm) in diameter.



### Support Rings, plated cast-iron

3" OD (2 1/2" ID). **600SRC103** ..... \$9.55  
 4" OD (3 1/2" ID). **600SRC104** ..... \$12.45  
 5" OD (4 1/2" ID). **600SRC105** ..... \$13.25



Plain



Ceramic  
Center

### Wire Gauze Squares, galvanized iron wire w/turned edges

**Plain**  
 4 x 4". **31040004** ..... \$1.45  
 5 x 5". **31040005** ..... \$1.30  
 6 x 6". **31040006** ..... \$1.55

**Ceramic Center**  
 4 x 4". **31040001** ..... \$1.60  
 5 x 5". **31040002** ..... \$1.70  
 6 x 6". **31040003** ..... \$2.25

Continued