

Drying Racks



188050000

189330023



HS23243A

Drying Racks, Benchtop, Bel-Art

Single-sided. **188050000**..... \$174.00

Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 43 x 17 x 48 cm(H).

Double-sided

2-Tier. **189330022** \$433.00

3-Tier. **189330023** \$589.00

Drying Rack, Benchtop/Wallmount, Heathrow

Single-sided. **HS23243A** \$148.00

Dry up to 90 pieces of labware with this self-standing or wall-mounted rack. Constructed of vinyl-coated wire hairpin loops with a polystyrene tray and rack.

Dimensions: 46.2 x 18.2 x 52.5 cm(H).

Electrophoresis



Electrophoresis System, MyGel™ Mini. **E1101** \$649.00

Eyewear

Absorbs
99.9% UV



Glasses, Safety, integral side shield. **43574002**..... \$6.00



Flexible frame and adjustable elastic strap provide a tight seal.

Goggles, Safety, vented (chemical splash). **97-2001**..... \$4.40

Filter & Blotting Papers



Ahlstrom - Munksjo
Filter Paper

Qualitative, Standard Grades (for gravity filtration)

Slow flow (4 mL/min), 1.5 µm retention, Grade 610

(Whatman 5; S&S 602)

9 cm, Pk/100. **6100-0900** \$12.90

12.5 cm, Pk/100. **6100-1250** \$16.10

15 cm, Pk/100. **6100-1500** \$19.40

20.5 cm, Pk/100. **6100-2050** \$30.65

24 cm, Pk/100. **6100-2400** \$42.70

Medium flow (30 mL/min), 2 µm retention, Grade 642

(Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)

9 cm, Pk/100. **6420-0900** \$10.05

15 cm, Pk/100. **6420-1500** \$18.05

Medium flow (55 mL/min), 2.5 µm retention, Grade 601

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

4.25 cm, Pk/100. **6010-0425** \$7.70

9 cm, Pk/100. **6010-0900** \$10.35

12.5 cm, Pk/100. **6010-1250** \$14.85

15 cm, Pk/100. **6010-1500** \$18.10

18.5 cm, Pk/100. **6010-1850** \$26.25

20.5 cm, Pk/100. **6010-2050** \$28.30

33 cm, Pk/100. **6010-3300** \$76.40

Medium flow (60 mL/min), 6 µm retention, Grade 613

(Whatman Student; S&S 595, 0903)

9 cm, Pk/100. **6130-0900** \$8.80

15 cm, Pk/100. **6130-1500** \$15.35

20.5 cm, Pk/100. **6130-2050** \$23.00

33 cm, Pk/100. **6130-3300** \$51.15

Medium flow (200 mL/min), 10 µm retention, Grade 631

(Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)

20.5 cm, Pk/100. **6310-2050** \$28.50

33 cm, Pk/100. **6310-3300** \$77.20

Fast flow (235 mL/min), 25 µm retention, creped, Grade 615

9 cm, Pk/100. **6150-0900** \$9.10

12.5 cm, Pk/100. **6150-1250** \$13.20

15 cm, Pk/100. **6150-1500** \$15.85

20.5 cm, Pk/100. **6150-2050** \$23.70

33 cm, Pk/100. **6150-3300** \$52.70

Very Fast flow (360 mL/min), 35 µm retention, creped, Grade 617

9 cm, Pk/50. **6170-0900** \$8.30

12.5 cm, Pk/50. **6170-1250** \$11.05

15 cm, Pk/50. **6170-1500** \$12.85

20.5 cm, Pk/50. **6170-2050** \$19.60

33 cm, Pk/50. **6170-3300** \$42.00

Qualitative, Pleated (Fluted)

Slow flow (30 mL/min), 2 µm retention, Grade 542

(Whatman 2V)

12.5 cm, Pk/100. **5420-1250** \$82.45

18.5 cm, Pk/100. **5420-1850** \$93.30