

**Qualitative, Wet-Strengthened Grades** (for vacuum filtration)

- Medium flow** (65 mL/min), 3 µm retention, **Grade 909**  
 (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)  
 32 cm, Pk/100. **9090-3200** ..... \$83.20
- Fast flow** (245 mL/min), 27 µm retention, creped, **Grade 962**  
 (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)  
 11 cm, Pk/100. **9620-1100** ..... \$12.35

**Quantitative, Ashless**

(Designed for use in analytical and gravimetric analysis)

- Slow flow** (4 mL/min), 1.5 µm retention, **Grade 94**  
 (Whatman 42; S&S 589/3)  
 9 cm, Pk/100. **0940-0900** ..... \$34.85  
 15 cm, Pk/100. **0940-1500** ..... \$60.00
- Medium flow** (30 mL/min), 2 µm retention, **Grade 74**  
 (Whatman 40; S&S 589/2, 589/4)  
 9 cm, Pk/100. **0740-0900** ..... \$34.85  
 15 cm, Pk/100. **0740-1500** ..... \$60.00
- Fast flow** (220 mL/min), 15 µm retention, **Grade 55**  
 (Whatman 54, 541; S&S 1505)  
 9 cm, Pk/100. **0550-0900** ..... \$40.00  
 15 cm, Pk/100. **0550-1500** ..... \$63.25

**Blotting/Chromatography/Electrophoresis Papers**


- 0.19 mm thick**, 125 mm/30 min absorption, **Grade 601**  
 (Whatman 1, 1CHR; S&S 593, 597, 2043A)  
 46 x 57 cm, Pk/100. **6018-4657** ..... \$157.00
- 0.34 mm thick**, 140 mm/30 min absorption, **Grade 238**  
 (Whatman 3MM, S&S 598, 231<sup>o</sup>)
- 20 x 20 cm, Pk/100. **2388-2020** ..... #1 \$80.00  
 46 x 57 cm, Pk/100. **2388-4657** ..... \$305.00
- 0.42 mm thick**, 150 mm/30 min absorption, **Grade 237**  
 (Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703; Albet 2040B, FP591; Munktell 3S/H, FN7a)
- NEW** 46 x 57 cm, Pk/100. **2378-4657** ..... \$352.00
- 0.95 mm thick**, 180 mm/30 min absorption, **Grade 243**  
 (Whatman 17CHR; S&S 2668)  
 11 x 14 cm, Pk/100. **2438-1114** ..... \$31.75  
 20 x 20 cm, Pk/100. **2438-2020** ..... \$66.15
- 2.48 mm thick**, 240 mm/30 min absorption, **Grade 320**  
 (S&S 2727)  
 11 x 14 cm, Pk/50. **3208-1114** ..... \$63.45  
 20 x 20 cm, Pk/50. **3208-2020** ..... \$93.05



Smooth tips & polypropylene grips

- Forceps, Filter**, 11.5 cm. **300-025** ..... \$4.40

**Freezer Racks**



**Freezer Racks, stainless steel, USA**

175V-13-2L

Restocking fee of 20% on all returned freezer racks listed below.  
 No returns on custom racks.

**Horizontal** (fits upright freezers)

- Horizontal, Holds 2" Boxes**
- |   | 1-9      | ≥10    |
|---|----------|--------|
| 9-Place (3L x 3H). <b>175HS-9-2</b> .....   | \$59.00  | 53.00  |
| Dimensions: 16.31(L) x 5.5(W) x 6.5" (H)    |          |        |
| 12-Place (3L x 4H). <b>175HS-12-2</b> ..... | \$78.00  | 71.00  |
| Dimensions: 16.31(L) x 5.5(W) x 8.75" (H)   |          |        |
| 15-Place (3L x 5H). <b>175HS-15-2</b> ..... | \$98.00  | 89.00  |
| Dimensions: 16.31(L) x 5.5(W) x 11" (H)     |          |        |
| 16-Place (4L x 4H). <b>175HL-16-2</b> ..... | \$105.00 | 95.00  |
| Dimensions: 22(L) x 5.5(W) x 8.75" (H)      |          |        |
| 20-Place (4L x 5H). <b>175HL-20-2</b> ..... | \$132.00 | 119.00 |
| Dimensions: 22(L) x 5.5(W) x 11" (H)        |          |        |

- Horizontal, Holds 3" Boxes**
- |   |         |       |
|---|---------|-------|
| 9-Place (3L x 3H). <b>175HS-9-3</b> .....   | \$61.00 | 55.00 |
| Dimensions: 16.31(L) x 5.5(W) x 9.5"(H)     |         |       |
| 12-Place (4L x 3H). <b>175HL-12-3</b> ..... | \$80.00 | 73.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H)        |         |       |

**★ Horizontal Drawer** (fits upright freezers)

- Horizontal Drawer, Holds 2" Boxes**
- |  | 1-9      | ≥10    |
|--|----------|--------|
| 9-Place (3L x 3H). <b>175HD-332</b> .....  | \$85.00  | 77.00  |
| Dimensions: 16.5(L) x 5.5(W) x 7.125"(H)   |          |        |
| 12-Place (3L x 4H). <b>175HD-342</b> ..... | \$111.00 | 101.00 |
| Dimensions: 16.5(L) x 5.5(W) x 9.5" (H)    |          |        |
| 15-Place (3L x 5H). <b>175HD-352</b> ..... | \$141.00 | 128.00 |
| Dimensions: 16.5(L) x 5.5(W) x 11.812" (H) |          |        |
| 16-Place (4L x 4H). <b>175HD-442</b> ..... | \$151.00 | 137.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5" (H)      |          |        |
| 20-Place (5L x 4H). <b>175HD-542</b> ..... | \$188.00 | 171.00 |
| Dimensions: 27.25(L) x 5.5(W) x 9.5" (H)   |          |        |

- Horizontal Drawer, Holds 5" Boxes**
- |   | 1-9     | ≥10   |
|---|---------|-------|
| 6-Place (3L x 2H). <b>175RD-325</b> ..... | \$74.00 | 67.00 |
| Dimensions: 18.125(L) x 6(W) x 10.75" (H) |         |       |
| 8-Place (4L x 2H). <b>175RD-425</b> ..... | \$99.00 | 90.00 |
| Dimensions: 24(L) x 6(W) x 10.75" (H)     |         |       |

- Horizontal Drawer, Holds 15 mL tubes**
- |   |          |        |
|---|----------|--------|
| 160-Place (80L x 2H). <b>175HD-32ML-15</b> .....  | \$99.00  | 90.00  |
| Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)            |          |        |
| 220-Place (110L x 2H). <b>175HD-42ML-15</b> ..... | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5" (H)             |          |        |

- Horizontal Drawer, Holds 50 mL tubes**
- |   |          |        |
|---|----------|--------|
| 84-Place (42L x 2H). <b>175HD-42ML-50</b> ..... | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H)            |          |        |

Continued