




Safely grip hot or cold labware with this silicone rubber hand protector partially lined with mini suction cups. Withstands temperatures from 260°C (500°F) to -57°C (-70°F).

**Hand Protector, red. 600HNDPTC** ..... \$20.40

### Homogenizers



Rapidly homogenize, emulsify, suspend and disrupt 0.1 mL to 250 mL samples. Includes a 5 x 50 mm flat-bottom probe & 7 x 50 mm saw tooth probe for use with micro tubes.

**Homogenizer, Handheld. D1000**..... \$1,269.00



Benchmark BeadBlaster™ 24

Homogenize a variety of sample types, from soft tissue to bone, in 35 seconds or less.

**Homogenizer, Benchtop, Microtube. D2400**..... \$8,404.00

- Holds (24) 2.0 mL tubes, or (6) 5 mL tubes w/optional tube holder
- Speed adjusts from 4.00 to 7.00 m/s in 0.05 m/s increments
- Stores 50 programs in memory, each with a maximum of 10 cycles
- 2-year warranty

#### Homogenizer Prefilled Tubes

2.0 mL screw-capped tubes are prefilled with your choice of the following beads:

Glass, acid washed

- 0.1 mm, Pk/50. **D1031-01**.....\$135.00
- 0.5 mm, Pk/50. **D1031-05**.....\$135.00

Zirconium Beads, acid washed/heat treated/nuclease-free

- 1.5 mm, Pk/50. **D1032-15**.....\$189.00
- 3.0 mm, Pk/50. **D1032-30**.....\$189.00

Stainless Steel Beads, acid washed

- 2.8 mm, Pk/50. **D1033-28**.....\$189.00

### Hot Plates, Hot Plate Stirrer & Magnetic Stirrers



**Hot Plates, analog**

- 6" Diameter Cast Iron Burner, 1000W. **310HOTPLT** ..... \$32.40
- 7.5 x 7.5" Ceramic Platform, 80-380°C. **H4000-H** ..... \$312.00

H4000-HS



**Hot Plate Stirrer, 80-380°C, Benchmark. H4000-HS** ..... \$449.00



S1005

H4000-S



**Magnetic Stirrers, analog, Benchmark**

- 3.4 x 3.4" Platform, 100-1500 rpm, MiniMag™. **S1005**..... \$148.00 (1-year warranty)
- 7.1 x 7.1" Ceramic Platform, 60-1500 rpm. **H4000-S**..... \$312.00 (2-year warranty)