

Incubators



Incubators, Mini Digital, MyTemp™, Benchmark (2-year warranty)

Heat Only, 1° above ambient to 60°C. **H2200-H**\$609.00

Heat/Cool, 15° below ambient to 60°C. **H2200-HC**\$917.00

Heat/Cool/CO₂, 10° below ambient to 55°C. **H2300-HC2** \$4,169.00

- 20 L capacity (0.8 cu ft)
- Digital temperature control with LED display
- Accepts bottles and flasks up to 2 L
- Two adjustable/removable shelves
- 9.3"(D) x 10.3"(W) x 12.8"(H) interior;
14.5"(D) x 13.2"(W) x 18.7"(H) exterior



Incubators, Dual-Convection, SureTemp™ (2-year warranty)

40 L. **H2505-40**\$2,139.00

13.7"(D) x 13.7"(W) x 13.7"(H) interior

70 L. **H2505-70**\$2,699.00

13.7"(D) x 17.7"(W) x 17.7"(H) interior

130 L. **H2505-130**\$3,449.00

17.7"(D) x 21.7"(W) x 21.7"(H) interior

- Precisely control temperature up to 75°C
- Record temperatures with SureCheck™ data logging
- Select either gravity or mechanical convection mode



Incubator, Shaking, Incu-Shaker™ Mini. **H1001-M**\$3,199.00

9.5" x 11.5" platform w/ rubber mat

- 5° above ambient to 70°C
- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Domed lid accepts flasks up to 1 L
- 16.2"(D) x 11"(W) x 13"(H) exterior

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

Steel Platform, Magic Clamp™, 9.5 x 11". **H1000-MR**\$326.00

Clamps, magnetic

250 mL Flask. **H1000-MR-250**\$49.85

500 mL Flask. **H1000-MR-500**\$49.85

1000 mL Flask. **H1000-MR-1000**\$76.25

Mortar and Pestle Sets

Mortar and Pestle Sets

Glass, **USA**

59 mL (2 oz) Capacity. **826-2**\$16.35

118 mL (4 oz) Capacity. **826-3**\$17.20

237 mL (8 oz) Capacity. **826-4**\$20.40

Porcelain, outside wall and rim glazed

30 mL Capacity, 63 mm top OD. **LPC-154**\$5.00

70 mL Capacity. **600JMD070**\$13.35

150 mL Capacity. **600JMD150**\$17.60

320 mL Capacity, 135 mm top OD. **LPC-156**\$11.50



Parafilm



Parafilm, thermoplastic sealing film, **USA**

2" x 250' (5 cm x 76 m) Roll. **HS234526A**\$34.00

4" x 125' (10 cm x 38 m) Roll. **HS234526B**\$34.90

Petri Dishes



Petri Dishes, polystyrene, Kord-Valmark

60 x 15 mm, stackable

Cs/500, sterile sleeve/20. **2901**\$118.00

100 x 15 mm

Slippable (rimless lid for automation)

Standard, Cs/500, sterile sleeve/25. **2917**\$126.00

★ Stackable

Standard, Cs/500, sterile sleeve/25. **2900**\$126.00

100 x 25 mm, stackable

Cs/340, sterile sleeve/17. **2906**\$137.00

pH Paper



Measure pH with 1 pH unit accuracy from 1-14.

Package contains 100 paper strips, 2 x 1/4", in a plastic vial with color chart.

pH Paper, pH 1-14, Pk/100. 600PHP100\$4.00