



These mini, air displacement pipettors provide an economical and accurate method for pipetting fixed volumes. Pipettors feature a two-step plunger and are fully autoclavable.

Pipettors, Single Channel

Fixed Volume, pocket-size, Diamond® Jr.

5 µL, red. 3322-5	\$32.95
10 µL, pink. 3322-10	\$32.95
20 µL, green. 3322-20	\$32.95
50 µL, purple. 3322-50	\$32.95
100 µL, orange. 3322-100	\$32.95

Single channel, adjustable volume pipettors are fully autoclavable.



Adjustable Volume, Diamond® Apex™ (3-year warranty)

0.5-10 µL, 0.02 µL increments, red. 3351-10	\$219.00
5-50 µL, 0.1 µL increments, yellow. 3351-50	\$219.00
10-100 µL, 0.1 µL increments, blue. 3351-100	\$219.00
100-1000 µL, 1 µL increments, orange. 3351-1000	\$219.00

Pitchers



Polypropylene



Stainless Steel

Pitchers

Polypropylene, graduated (mL), tall form, tapered	
1000 mL, 50 mL graduations, each. 289910000	\$7.20
2000 mL, 100 mL graduations, each. 289920000	\$8.35
3000 mL, 100 mL graduations, each. 289930000	\$13.80
Stainless Steel, graduated (mL/oz), tapered, Polar Ware	
1000 mL/32 oz, each. T1063	\$73.20

Rockers, Shakers & Mixers



Rockers, Fixed Speed, Benchmark (2-year warranty)

2-D (seesaw), 12 x 12" platform w/flat mat	
BenchBlotter™. BR1000	\$557.00
<ul style="list-style-type: none"> • Fixed speed of 12 rocks/minute • Platform angle is ±8° • Cold room & incubator safe (4-45°C) • Load capacity is 2.3 kg (5 lbs) 	
3-D (nutating)	
Small, 10.5 x 7.5" platform w/dimpled mat	
MiniMixer™. B3D1020	\$459.00
Large, 12 x 12" platform w/dimpled mat	
BioMixer™. B3D1320	\$566.00
<ul style="list-style-type: none"> • Fixed speed of 24 rocks/minute • Platform angle is ±20° • Cold room & incubator safe (4-45°C) • Load capacity is 1.8 kg (4 lbs) 	

Continued