


## Rotating Mixers



**Tube Rotator**, programmable, Benchmark (2-year warranty)  
 Roto-Bot Mini. **R4040** ..... \$766.00  
 Roto-Bot. **R4045** ..... \$1,115.00


- Speed adjusts from 1-80 rpm
- Tube capacity (starter set)
  - R4040: 12 x 1.5 mL, 2 x 15 mL & 2 x 50 mL
  - R4045: 30 x 1.5 mL, 8 x 15 mL & 2 x 50 mL
- Timer is programmable from 1 sec – 99 hrs, 59 min / Cont.

## Vortexers



**Vortex Mixers**, Benchmark (2-year warranty)  
 Vornado™, mini vortexer, fixed 2800 rpm. **BV101** \* ..... \$233.00  
 \*Specify Gray, Blue or Purple cup head  
 BenchMixer™, single cup head, 200-3200 rpm. **BV1000** ..... \$302.00  
 Flat Head Attachment, 3". **BV1000-FLAT** ..... \$46.55

## Rulers, Meter Stick & Tape Measure



**Rulers**, 15 cm/6" long, scaled in both metric and English units  
 Vinyl, flexible  
 Transparent. **4461** ..... \$0.95  
 White. **4460** ..... \$0.55  
 Stainless Steel, 12 mm wide, 0.5 mm thick, **USA. 605-065** ..... \$6.70

**Meter Stick**, lacquered wood, plain ends. **605-050** ..... \$4.90  
 One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

**Tape Measure**, 152.4 cm/60", reinforced fiberglass. **11036** ..... \$3.70  
 One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.

## Sharps Containers



4402269



440321391

**Sharps Containers**, red polypropylene  
 5.4 Quart, horizontal entry lid, Prevent™, each. **4402269** ..... \$9.90  
 Dimensions: 12" x 4.75" x 11"(H)  
 12 Gallon, hinged vertical entry lid, each. **440321391** ..... \$32.20  
 Dimensions: 18.25" x 12.75" x 18.75"(H)

## Specimen Containers

Screw-Cap, 90 mL



Screw-Cap,  
480 & 960 mL



Snap-Cap

## Specimen Containers

Screw-Cap, polypropylene (PP)  
 90 mL, sterile, Cs/400. **840175** ..... \$150.00  
 480 mL (16 oz), Cs/24. **59014510** ..... \$71.95  
 83 mm (3 1/4") opening  
 960 mL (32 oz), Cs/12. **59014511** ..... \$59.05  
 111 mm (4 3/8") opening  
 Snap-Cap, white HDPE bottom/white LDPE lid  
 480 mL (16 oz), Cs/100. **271016W** ..... \$95.75  
 9.9 cm opening x 11.1 cm high  
 1000 mL (34 oz), Cs/100. **271032W** ..... \$144.00  
 11.7 cm opening x 13.7 cm high

## Spot Plate



**Spot Plate**, glazed white porcelain, 21 mm diameter wells  
 12-Well, 115 x 95 x 7 mm, each. **3104187701** ..... \$4.50