

Plastic culture tubes with dual-position caps.



Tubes, Culture,

- 12 x 75 mm (5 mL)
 - Polystyrene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110428** \$98.70
 - Polystyrene, Nonsterile, Cs/1000. **110410** \$64.00
 - Polypropylene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110438** \$98.90
- 13 x 100 mm (10 mL)
 - Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1510** \$106.00
- 17 x 100 mm (16 mL)
 - Polypropylene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110178** \$159.00
- 16 x 125 mm (19 mL)
 - Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1515** \$175.00
- 16 x 150 mm (23 mL)
 - Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1517** \$180.00



Culture Tube Caps, polypropylene
Up to 16 mm Tubes

- Blue, Cs/1000. **118152B** \$110.00
- Green, Cs/1000. **118152G** \$110.00
- Red, Cs/1000. **118152R** \$110.00
- Yellow, Cs/1000. **118152Y** \$110.00



15 mL



Black

RNase/DNase-free

Tubes, Microcentrifuge, Economy, polypropylene

- 0.5 mL, 0.1 mL graduations, natural, Pk/1000. **HS4422** \$39.10
- 1.5 mL, 0.25 mL graduations, Heathrow
 - Natural, max RCF is 16,400 x g
Pk/500. **HS4323** \$19.00
 - Pk/1000. **HS4325** \$35.95
 - Black, max RCF is 18,000 x g, Pk/500. **HS4323K** \$21.15
- 5 mL, natural, max RCF is 25,000 x g, Benchmark
 - Nonsterile, 4 bags/50, Pk/200. **C1005T5** \$68.45
 - Sterile, 10 bags/20, Pk/200. **C1005T5ST** \$75.10



0.6 mL



1.5 mL



2.0 mL

Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable. Free of RNase/DNase.

Tubes, Microcentrifuge, Premium, PP, G-Tube™, Bio Plas, USA

Amber color protects light sensitive specimens.

- 0.6 mL, Pk/500 \$28.20
 - Amber. **4048** Lavender. **4047** Natural. **4040**
 - Mixed (blue, green, orange, pink, & yellow). **4049** \$34.00
- 1.5 mL, 0.5 mL graduations, frosted area on side, Pk/500 \$28.20
 - Amber. **4038** Blue. **4033** Lavender. **4037** Natural. **4030**
 - Mixed (blue, green, orange, pink, & yellow). **4039** \$34.00
- 2.0 mL, 0.5 mL graduations, frosted area on side
 - Natural, Pk/500. **4050** \$30.30



Tubes, PCR, 0.2 mL thin-wall, polypropylene, Simport

- Individual, flat caps, natural, Pk/1000. **T325-2VN** \$51.20



Vials, Cryogenic, round bottom, internal thread

- O-Ring Seal, 4 mL, Cs/1000. **T301-4** \$302.00
- Washer Seal, 4 mL, Cs/1000. **T311-4** \$302.00



Vinyl grip, **USA**

Vial Decapping Pliers

- 13 mm Crimp Seals. **D-13** \$72.00
- 20 mm Crimp Seals. **D-20** \$72.00