

Tube Storage

Cryogenic Storage Boxes

SB2CC-64 Cryogenic
Tube Storage Box



Tube Storage Boxes, Cardboard, vented

| | |
|---|---------|
| 2" Box Only, 5¼ x 5¼ x 2", each. SB2CC | \$1.95 |
| 2" Box w/dividers, white | |
| 100-Cell (10 x 10), Pk/12. SB2CC-100 | \$31.00 |
| ★ 81-Cell (9 x 9), Pk/12. SB2CC-81 | \$30.75 |
| 64-Cell (8 x 8), Pk/12. SB2CC-64 | \$30.75 |
| 2" Mini box w/dividers, 3 x 3 x 2" | |
| 25-Cell (5 x 5), Pk/12. MB2C-25 | \$33.40 |
| Holds 13 mm vials and 4.0-5.0 mL cryovials. | |
| 3" Box Only, 5¼ x 5¼ x 3", each. SB3C | \$1.95 |
| 3" Box w/dividers | |
| 100-Cell (10 x 10), Pk/12. SB3C-100 | \$31.00 |
| 81-Cell (9 x 9), Pk/12. SB3C-81 | \$30.75 |
| 64-Cell (8 x 8), Pk/12. SB3C-64 | \$30.75 |
| Dividers only, Fit 2 & 3" Boxes | |
| 100-Cell (10 x 10), Pk/12. CD-100 | \$12.90 |
| Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials. | |
| 81-Cell (9 x 9), Pk/12. CD-81 | \$12.60 |
| Holds 13 mm vials and 4.0-5.0 mL cryovials. | |
| 64-Cell (8 x 8), Pk/12. CD-64 | \$12.60 |
| Holds 14 mm vials. | |
| 49-Cell (7 x 7), Pk/12. CD-49 | \$12.90 |
| Holds 14 mm vials. | |
| 25-Cell (5 x 5), 25 mm cell, Pk/12. CD-25 | \$12.60 |
| 16-Cell (4 x 4), 31 mm cell, Pk/12. CD-16 | \$12.60 |
| 5" Box Only, 5¾ x 5¾ x 5", each. SB4C | \$3.60 |
| 5" Box w/dividers | |
| 64-Cell (8 x 8), Pk/12. SB4C-64 | \$47.40 |
| 36-Cell (6 x 6), 15 mL tubes, Pk/12. SB4C-36 | \$47.40 |
| 16-Cell (4 x 4), 50 mL tubes, Pk/12. SB4C-16 | \$47.40 |
| Dividers only, Fit 5" Box | |
| 64-Cell (8 x 8), Pk/12. CD4-64 | \$12.60 |
| 36-Cell (6 x 6), Pk/12. CD4-36 | \$12.60 |
| 16-Cell (4 x 4), Pk/12. CD4-16 | \$12.60 |



Holds 1.0-2.0 mL tubes and cryovials.
Features a locking hinged lid, and base
with molded alphanumeric coordinates.
Dimensions: 12.8 x 12.8 x 5.4 cm.

Tube Storage Boxes, Polypropylene


| | |
|---|---------|
| 81-Well (9x9), 2 blue/2 yellow boxes, Pk/4. 60066501 | \$36.50 |
|---|---------|

Tube Racks




Tube Racks, Grip, polypropylene, No-Wire™, Bel-Art

| | |
|--|---------|
| 18-20 mm Tubes, 40-place, each. 187490002 | \$36.95 |
|--|---------|



Tube Racks, Interlocking, polypropylene, Heathrow

| | |
|---|---------|
| 4-Way, mixed colors, Pk/5. HS29022G | \$47.40 |
| Holds (4) 50 mL tubes, (12) 15 mL tubes, (32) 1.5 mL tubes, or (32) 0.5 mL tubes. | |
| 4-Way Rotatable, Rota-Rack™, four connected modules | |
| Large Tubes, 22.4 x 11 x 11 cm, each. HS29040A | \$49.30 |
| Each module holds (2) 50 mL tubes, (8) 15 mL tubes, (6) 20 mm tubes, or (10) 10-13 mm tubes. Mixed colors. | |



Tube Racks, Mega Rack, polypropylene, fully-assembled

| | |
|------------------------------------|---------|
| 10-13 mm Tubes | |
| 216-Place, 33.7 x 26.7 x 9.1 cm(H) | |
| Blue, each. 34524322A | \$51.60 |
| White, each. 34524322E | \$51.60 |
| 432-Place, 67.4 x 26.7 x 9.1 cm(H) | |
| Blue, each. 345120113 | \$90.25 |
| White, each. 345120112 | \$90.25 |



Tube Racks, Microtube, polypropylene, Heathrow

| | |
|--|---------|
| 80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm | |
| Holds 1.5-2.0 mL tubes in conical wells. | |
| Blue, Pk/5. HS29025B | \$28.50 |
| Pink, Pk/5. HS29025D | \$28.50 |
| Mixed, Pk/5. HS29025G | \$28.50 |