

**Mercury Short Arc**

50W Mercury Arc, 100 Hrs, Osram. <b>HBO50W/AC/L1</b> .....	\$115.00
<i>(78564, 78588, 8-B199, 38-16-19, 414140, 2055, 3040-0090, 500-137, SHL-50AC1, BA-001, 5000487)</i>	
50W Mercury Arc, 100 Hrs, Osram. <b>HBO50W/AC/L2</b> .....	\$115.00
<i>(31-31-65, SHL-50AC2, BA-002, 5000857)</i>	
50W Mercury Arc, 200 Hrs, Osram. <b>HBO50W/3</b> .....	\$140.00
<i>(BA-003)</i>	
100W Mercury Arc, 200 Hrs, Ushio. <b>USH-102D</b> .....	\$109.00
<i>(78589, 8-B192U, 500-138, 3800-18-4060/D, 2087-001, 3040-0380, BA-005A, HBO100W/2, M/5000/133, 5000274)</i>	
100W Mercury Arc, 300 Hrs, Ushio. <b>USH-103D</b> .....	\$145.00
<i>(78590, 8-B1920, 380301-9350, USH-103OL, 5001329, HBO103W/2)</i>	
200W Mercury Arc, 1000 hrs, Ushio. <b>USH-200DP</b> .....	\$135.80
<i>(HBO200W/2, 31-01-01, 649, 860014, 78580, 8-B190, 3800-18-4070, 38-16-15, 050-320, MA885, 91-00-50, 505-620, 78590, 5000276)</i>	

**Xenon Short Arc**

75W Xenon Arc, 400 hrs, Ushio. <b>UXL-75XE</b> .....	\$195.00
<i>(XBO75W/2, 8-B198, 910881-2013, 38-16-02, 500-139, 91-01-42, 3800-79-9190, 78682, 620-125, 903-050, BA-013, 5000371)</i>	
150W Xenon Arc, Short Gap. <b>XBO150W/CROFR</b> .....	\$260.00
<i>(8-B195, BA-018, UXL-150S-0, 3801-22-2450, 500-008, 5000326)</i>	
180W Xenon Arc, Osram. <b>XBOR180W45COFR</b> .....	\$450.00
<i>(69183, 407.309)</i>	

**Fluorescent**

4W Fluorescent, Tubular, Eiko. <b>F4T5/D</b> .....	\$7.00
<i>(31-31-03, MA309/05, 3000108)</i>	
5W Fluorescent, G23, Sylvania. <b>CF5DS/841</b> .....	\$11.70
<i>(21278, 800-138, MA68, 48920-81, 20303)</i>	
9W Fluorescent, G23, Sylvania. <b>CF9DS/841</b> .....	\$4.95
<i>(20305)</i>	
22W Fluorescent, 8" OD Circle, Sylvania. <b>FC8T9/D/RS</b> .....	\$19.00
<i>(29980)</i>	
32W Fluorescent, 12" OD Circle, Sylvania. <b>FC12T9/D/RS</b> .....	\$15.90
<i>(20030)</i>	

**Germicidal**

4W Germicidal (UV-C), G5, Osram. <b>G4T5/OF</b> .....	\$10.35
<i>(23397)</i>	
15W Germicidal (UV-C), G13, Sylvania. <b>G15T8</b> .....	\$17.00
<i>(21612, 3000007)</i>	
25W Germicidal (UV-C), G13, Ushio. <b>G25T8</b> .....	\$19.80
<i>(3000008)</i>	
39W Germicidal, single-pin, Ushio. <b>G36T5L</b> .....	\$35.80
<i>(23443, 00206, 3000312)</i>	
65W Germicidal, single-pin, 155.4 cm, Higuchi. <b>G64T5L</b> .....	\$83.00

**Heat (Infrared)**

250W, R40, shatter resistant, GE. <b>250BR40/HI/TC</b> .....	\$14.00
--	---------

\* Recommend Bulb Removal/Installation Tool (shown below)



Remove and install prefocused/flanged bulbs safely with this simple tool. Push rubber-lined tube directly onto top of bulb and turn clockwise to install, or counterclockwise to remove.

**Bulb Removal/Installation Tool. 3283**..... \$16.00

**Batteries**

**Batteries**, button cell, 11.6 mm diameter

Silver Oxide, 1.55V, 5.4 mm high, Energizer. <b>357</b> .....	\$4.00
<i>(303, A76, AG13, EPX76, LR1154, LR44)</i>	
Alkaline, 1.5V, 3.05 mm high. <b>389</b> .....	\$2.25
<i>(LR1130, LR54, SR1130, AG10, 189)</i>	

**Fuses**



**Fuses**

$\frac{1}{4} \times \frac{1}{4}$ "	
Fast-Acting, AGC	
0.1A/250V. <b>312-100</b> .....	\$2.00
0.2A/250V. <b>312-200</b> .....	\$2.00
1A/250V. <b>312-001</b> .....	\$2.00
2A/250V. <b>312-002</b> .....	\$2.00
3A/250V. <b>312-003</b> .....	\$2.00
4A/250V. <b>312-004</b> .....	\$2.00
Time-Delay (Slow-Blow), MDL	
0.7A/250V. <b>313-700</b> .....	\$2.00
1A/250V. <b>313-001</b> .....	\$2.00
2.5A/250V.....	\$2.00
5 x 20 mm	
Fast-Acting, GMA	
0.5A/250V. <b>235-500</b> .....	\$2.00
1A/250V. <b>235-001</b> .....	\$2.00
3A/250V. <b>235-003</b> .....	\$2.00
5A/125V. <b>235-005</b> .....	\$2.00
Time-Delay (Slow-Blow), Glass tube	
0.5A/250V. <b>239-500</b> .....	\$2.00
1A/250V. <b>239-001</b> .....	\$2.00
2A/250V. <b>239-002</b> .....	\$2.00
2.5A/250V.....	\$2.00
3A/250V. <b>239-003</b> .....	\$2.00
5A/250V. <b>239-005</b> .....	\$2.00
Time-Delay (Slow-Blow), Ceramic	
1A/250V. <b>215-001</b> .....	\$2.00