


Microscope Filters, glass



| | |
|--|---------|
| 29.8 mm Diameter | |
| Blue, clear. MA289/05 | \$50.00 |
| Blue Frosted. MA289 | \$50.00 |
| Diffusion (frosted). MA288 | \$50.00 |
| 32 mm Diameter | |
| Blue, clear. 923-160 | \$15.00 |
| Blue Frosted. 3236 | \$7.50 |
| Diffusion (frosted). 3235-W | \$7.50 |
| Green, clear. CM-0321-G | \$7.00 |
| 40 mm Diameter | |
| Blue Frosted. MA563/05 | \$70.00 |
| Diffusion (frosted). MA566 | \$70.00 |
| 45 mm Diameter | |
| Blue, clear. 14671 | \$10.40 |
| Diffusion (frosted). 14673 | \$35.00 |
| Green Interference, 546/60, Nikon. MBN11200 | \$52.00 |
| Neutral Density, Nikon | |
| 50% Transmittance, ND2A. MBN21802 | \$42.00 |
| 12.5% Transmittance, ND8A. MBN21808 | \$42.00 |



Renew and improve stage surface. Cleans easily.
Light surface provides contrast for previewing specimens.
Additionally, PTFE coating aids moving slides on the stage.

| | |
|--|--------|
| Microscope Stage Laminate, USA. 16251 | \$8.00 |
| 6.25 x 7.75" sheet provided with a #10 scalpel. | |
| White Teflon® (PTFE) on self-adhesive vinyl is 0.010" thick. | |

Object Marker Accessories

| | |
|--|---------|
| Refill Ink, black, 10 mL, Shachihata. 22112 | \$7.00 |
| Adapters | |
| RMS/45 to CFI60, 15 mm long. 14111 | \$43.00 |
| RMS ID to 25 mm OD Thread, brass. 14115 | \$52.00 |

Microscopes with Nikon CFI60 infinity-corrected (∞) objectives require a 14111 adapter.
Microscopes with Leica 25 mm thread/45 mm parfocal length objectives require a 14115 adapter.
Microscopes with 160 mm tube length objectives do not require an adapter.



Objectives, 160 mm tube length, DIN, Meiji, Motic, National
Super High Contrast Achromats, National

| | |
|--|----------|
| 4X, N.A. 0.10. 704-160ASC | \$34.00 |
| 10X, N.A. 0.25. 710-160ASC | \$43.00 |
| 40X, N.A. 0.65. 740-160ASC | \$79.00 |
| 100X Oil, N.A. 1.25. 799-160ASC | \$117.00 |

Semi-Plan Achromats, National (unless noted)

| | |
|---|----------|
| 2X, N.A. 0.05, Accu-Scope. 3193SP-00 | \$70.00 |
| 4X, N.A. 0.10. 704-160SP | \$39.00 |
| 10X, N.A. 0.25. 710-160SP | \$40.00 |
| 40X, N.A. 0.65. 740-160SP | \$51.00 |
| 50X Oil, N.A. 0.95, Meiji. MA694 | \$330.00 |
| 60X, Strain Free, N.A. 0.85. MA645 | \$360.00 |
| 100X Oil, N.A. 1.25. 799-160SP | \$82.00 |

Plan Achromats, National (unless noted)

| | |
|--|----------|
| 4X, N.A. 0.10. 704-160P | \$40.00 |
| 10X, N.A. 0.25. 710-160P | \$65.00 |
| 20X, N.A. 0.40. 720-160P | \$147.00 |
| 40X, N.A. 0.65. 740-160P | \$123.00 |
| 50X Oil, N.A. 0.95, Accu-Scope. 3190PL-00 | \$315.00 |
| 60X, N.A. 0.85. 760-160P | \$232.00 |
| 100X Oil, N.A. 1.25. 799-160P | \$197.00 |

Protect spare microscope objectives.



| | |
|--|--------|
| Objective Case , DIN, polystyrene. 15601 | \$0.50 |
|--|--------|



Extender tubes do not contain lenses.

Objective Extender Tubes, RMS thread, brass, **USA**

| | |
|--------------------------------|---------|
| 10 mm Long. 14105 | \$44.20 |
| 18 mm Long. 14106 | \$42.75 |

Objective Shims, fit RMS thread (20.32 mm Dia.), **USA**

Parfocalize objectives with these stainless steel rings.

| | |
|---|---------|
| 25 µm (0.001") Thick, Pk/10. 10032061 | \$16.15 |
| 50 µm (0.002") Thick, Pk/10. 10032062 | \$16.15 |
| 75 µm (0.003") Thick, Pk/10. 10032063 | \$16.15 |
| 100 µm (0.004") Thick, Pk/10. 10032064 | \$16.15 |
| 125 µm (0.005") Thick, Pk/10. 10032065 | \$16.15 |
| Assortment Pack, pk/10. 100320630 | \$19.80 |

Includes two each of the following thicknesses: 25, 50, 75, 100, and 125 µm.

Highly efficient, linear polarizing film has excellent color reproducibility. Transmission is 38% for a single polarizer, while extinction is better than 99.99% for crossed polars. Film is 0.76 mm thick.



| | |
|--|---------|
| Polarizing Film , 50 x 50 mm Pair w/Plastic Case, USA. 15751 | \$15.00 |
|--|---------|