



Reticles & Stage Micrometers



Reticles, eyepiece, 1.5 mm thick, Klarmann Rulings, USA

*Specify reticle diameter using 2-digit suffix,
for example, use **KR-207-21** for 21 mm diameter reticle

★ 5 mm/100 Parts. KR-209-*	\$105.00
★ 10 mm/100 Parts. KR-207-*	\$105.00
Crossline 90°. KR-301-*	\$105.00
Crossline 90° with 10 mm/100 Parts Scale. KR-257-*	\$105.00
Indexed Square, 100 sqs, each 1 x 1 mm. KR-406B-*	\$105.00
1-10 on x-axis; A-J on y-axis	
Miller Disc	
★ Corner Square. KR-802-*	\$105.00
Central Square. KR-802A-*	\$105.00
Walton Beckett, G-22. KR-827-1-*	\$259.00
Whipple Disc. KR-413-*	\$105.00



Stage Micrometers, Transmitted Light

★ 1 mm/100 Parts Scale (0.01 mm), Meiji. MA285	\$130.00
1 mm-0.01 mm/div on 1 mm thick slide, protected with microcoverslip. Preferred for compound microscopes.	
10 mm/100 Parts Scale (0.1 mm), Meiji. MA187	\$130.00
10 mm-0.1 mm/div on 1 mm thick slide, protected with 1 mm thick coverslip. Designed for stereo microscopes.	
25 mm Scale, USA. KR-812	\$230.00
25mm-1 mm/div, 5 mm-0.1 mm/div, 2 mm-0.02 mm/div on 3 mm thick slide. Preferred at magnifications below 400X. Ideal for large field stereo microscopes.	
0.2 mm Square/400 Squares, Meiji. MA287	\$130.00
20 vertical lines, 1 mm long and spaced 0.01 mm apart, intersect 20 similar horizontal lines to create a central 0.2 mm square consisting of 400 squares, each 0.01 mm. 1 mm thick slide is protected with microcoverslip. Useful for x-y calibration.	
0.04"/40 Parts Scale (0.001"), Meiji. MA286	\$130.00
0.04"-0.001"/div on 1 mm thick slide, protected with microcoverslip.	

Stage Micrometer, Reflected Light

1 mm/100 Parts Scale (0.01 mm), Meiji. MA292	\$250.00
1 mm-0.01 mm/div on 3 mm thick metal slide	

Digital Cameras



National

Moticam A5 Digital Camera	\$495.00
<ul style="list-style-type: none"> 5.0 MP (2592 x 1944) resolution with 1/2.8" sCMOS Sensor (scientific Complementary Metal Oxide Semiconductor) CS ring adapter 4-dot calibration slide MotiConnect imaging software 	
<i>Requires Windows 10 operating system or higher</i>	



National

Moticam A8 Digital Camera	\$529.00
<ul style="list-style-type: none"> 8.0 MP (3840 x 2160) resolution with 1/3" sCMOS Sensor (scientific Complementary Metal Oxide Semiconductor) CS ring adapter 4-dot calibration slide MotiConnect imaging software 	
<i>Requires Windows 10 operating system or higher</i>	

Accessories



Adapter, CS-Mount to C-Mount, 5 mm spacer. 17110	\$15.00
<i>Do not use this adapter ring on an existing C-mount camera, as it will negatively impact the quality and parfocality of the imaging system.</i>	



Couplers, Diagnostic Instruments, USA

1X C-Mount for 38 mm ISO Photoport. D10NLC	\$135.00
<i>Fits newer Nikon and Leica microscopes</i>	

Couplers, National

1X C-Mount for 23.2 mm ID Eyetube. 930-023	\$29.00
<i>Lens free - this coupler is not parfocal</i>	