

**Locking**



**Forceps, stainless steel**  
 Self-Closing (cross-locking), squeeze handles to open tips  
 Fine Point, smooth tips  
 Economy, 13 cm (5 $\frac{1}{8}$ " ). **300-066** ..... \$5.45  
 Peer-Vigor, 11 cm (4 $\frac{1}{4}$ " ), Swiss. **57.324** ..... \$42.40  
 Medium Point, serrated tips, 17 cm (6 $\frac{5}{8}$ " ). **TWZ-948.00** ..... \$3.95




Slide-Locking, grid-serrated tips  
 Fine Point  
 Economy, 14 cm (5 $\frac{1}{2}$ " ). **57.0357A** ..... \$9.00  
 Dumont, 16 cm (6 $\frac{1}{4}$ " ). **TWZ-101.72** ..... \$45.75  
 Medium Point  
 Economy, 16 cm (6 $\frac{1}{4}$ " ). **203** ..... \$9.00  
 Dumont, 16 cm (6 $\frac{1}{4}$ " ). **TWZ-101.76** ..... \$43.60  
 Snap-Closing, angled, fine points, 15.5 cm (6" ). **3610** ..... \$5.80  
 serrated tips & handles, guide pin

**Tissue**



**Forceps, stainless steel**  
 Adson-Brown, 9 x 9 side grasping teeth, 12 cm. **17-574** ..... \$16.05  
 Adson Tissue, 1 x 2 teeth, 12 cm (4 $\frac{3}{4}$ " ). **17-560** ..... \$8.55  
 Iris Tissue, 1 x 2 teeth  
 Straight, economy, 10 cm (4" ). **26-295** ..... \$10.75  
 Angled, economy, 11.5 cm (4 $\frac{1}{2}$ " ). **19-350ANG** ..... \$5.00  
 Half Curve, standard, 10 cm (4" ). **45-822** ..... \$17.60





Russian Tissue, circular serrated tips, 15.5 cm (6" ). **2275** ..... \$5.25  
 Tissue, 1 x 2 teeth, 14 cm (5 $\frac{1}{2}$ " ). **17-145** ..... \$8.55

**Tongs**



**Tongs, Crucible, serrated jaws, stainless steel**  
 24.5 cm (9 $\frac{3}{4}$ " ). **320-101** ..... \$12.50  
 51 cm (20" ). **320-120** ..... \$30.15


**Hemostats & Clamps**



**Hemostats, serrated jaws**  
 Hartman Mosquito, stainless steel  
 Straight  
 Economy, 9 cm (3 $\frac{1}{2}$ " ). **12-180S** ..... \$6.05  
 Standard, 9 cm (3 $\frac{1}{2}$ " ). **20-352** ..... \$21.30  
 Curved  
 Economy, 9 cm (3 $\frac{1}{2}$ " ). **306-012** ..... \$9.10  
 Standard, 10 cm (4" ). **20-360** ..... \$18.75  
 Halstead Mosquito, 12.5 cm (5" ), stainless steel  
 Straight  
 Economy. **14-065** ..... \$10.00  
 Standard. **20-501** ..... \$22.45  
 Curved  
 Economy. **14-070** ..... \$10.00  
 Standard. **20-519** ..... \$22.85  
 Kelly, stainless steel  
 Straight  
 Economy, 14.5 cm (5 $\frac{3}{4}$ " ). **306-016** ..... \$8.20  
 Standard, 14.5 cm (5 $\frac{3}{4}$ " ). **20-550** ..... \$23.60  
 Curved  
 Economy, 14 cm (5 $\frac{1}{2}$ " ). **306-017** ..... \$9.90  
 Standard, 14.5 cm (5 $\frac{3}{4}$ " ). **20-568** ..... \$23.95  
 Rochester-Pean, 20.5 cm (8" ), stainless steel  
 Straight  
 Economy. **1920** ..... \$6.50  
 Standard. **21-977** ..... \$33.95  
 Curved  
 Economy. **1930** ..... \$6.50  
 Standard. **21-979** ..... \$33.95  
 Polypropylene, 12.5 cm (5" ). **HEM-300.00** ..... \$3.25



Specimen Forceps, 25 cm (9 $\frac{3}{4}$ " ), serrated jaws. **300-500** ..... \$18.70



**Sponge Forceps, Foerster, serrated jaws, stainless steel**  
 Straight  
 17.5 cm (7" ), economy. **14-660** ..... \$17.50  
 24 cm (9 $\frac{1}{2}$ " ), economy. **14-670** ..... \$21.10  
 Curved, 17.5 cm (7" ), economy. **14-665** ..... \$17.50

## Accessories

### Slide & Coverslip Handling

#### Forceps


##### Spade Tips

- Angled (Kuehne), 11 cm (4¼")  
 Stainless Steel, economy. **33-1041** ..... \$4.15  
 Teflon® Coated, Bel-Art. **379420000** ..... \$35.60  
 Straight, stainless steel, 10.5 cm, economy. **57.926** ..... \$8.45



##### Square Point Tips, straight

- Stainless Steel, 12 cm (4¾"), economy. **1244** ..... \$2.25  
 Teflon® Coated, 11 cm (4¼"), Bel-Art. **379430000** ..... \$39.40



Minimize contamination during slide & coverslip handling. Will not scratch; corrosion and heat resistant.

- Hemostat**, Polypropylene, 12.5 cm. **HEM-300.00** ..... \$1.95

### Immersion Oil & Bottles



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities.  
**USA**

#### Immersion Oil

##### Type A, low viscosity (150 cSt). Short W.D. optics

- 7.5 mL in LDPE dispensing bottle. **18101** ..... \$10.00  
 20 mL in HDPE dispensing bottle. **18102** ..... \$15.35  
 480 mL in amber bottle. **18104** ..... \$110.00

##### Type B, high viscosity (1250 cSt). Long W.D. optics

- 7.5 mL in LDPE dispensing bottle. **18105** ..... \$10.00  
 20 mL in HDPE dispensing bottle. **18106** ..... \$15.35  
 480 mL in amber bottle. **18108** ..... \$110.00

##### Type NVH, very high viscosity (21,000 cSt)

- 480 mL in amber bottle. **18113** ..... \$67.00


##### Type HF, low fluorescence (700 cSt)

- 480 mL in amber bottle. **18118** ..... \$169.00

##### Type FF, fluorescence free (170 cSt)

- 480 mL in amber bottle. **18124** ..... \$163.00


A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



- Bottle, Dispensing**, HDPE, 20 mL, each. **18145** ..... \$8.50

### Slide Storage

#### PLASTIC SLIDE BOXES



#### Slide Boxes, hold 75 x 25 mm (3 x 1") slides

- 25-Slide, squeeze-release cover, cork bottom, 9.5 x 8.5 x 3 cm  
 White, ABS plastic. **97-0025** ..... \$2.50  
 50-Slide, hinged cover, cork bottom, 20 x 9.5 x 3.5 cm  
 Blue, ABS plastic. **97-0056** ..... \$6.80  
 Red, ABS plastic. **97-0058** ..... \$6.80  
 100-Slide, steel-hinged cover, cork bottom  
 Blue, ABS plastic. **97-0126** ..... \$8.90  
 Red, ABS plastic. **97-0128** ..... \$8.90  
 White, ABS plastic. **97-0125** ..... \$8.90



Store 4500 75 x 25 mm slides in 6 divided drawers. Each cabinet comes with a sample of 6 compressor follow blocks and 6 drawer labels. Safely stack up to 12 cabinet units on a single base.  
 18¾"(D) x 15¼"(W) x 5"(H)

#### Slide Cabinets, steel

- 6-Drawer, stackable  
 Green. **513500G** ..... \$255.00  
 Tan. **513500T** ..... \$255.00  
 White. **513500W** ..... \$255.00  
 Memorandum Markers, paperboard, Pk/100. **C108** ..... \$12.75  
 Compressor Follow Block, rubber, each. **C109** ..... \$1.90  
 Base, raises cabinets ¾" off floor  
 Green. **513503G** ..... \$169.00  
 Tan. **513503T** ..... \$169.00  
 White. **513503W** ..... \$169.00



Cardboard Folder



Plastic Folder

#### Slide Folders, 20-slide, hold 75 x 25 mm (3 x 1") slides

- Cardboard, black tape, economy. **HS9902** ..... \$6.35 (12-587-10)  
 Plastic, gray tray w/2 clear doors. **M750-20GY** ..... \$11.60  
 Transparent covers hold visible slides securely in place, but can fold flat against back of tray when opened. Slides are slightly elevated for easy retrieval.