



Bone Shears, Constable, 23 cm (9") long, 50 mm (2") cutting jaws
 Economy. **33-1001** \$45.00
 German. **965-620** \$160.00



Postmortem Hammer, hook handle, 25 cm (10") long
 Economy. **33-1032** \$47.50
 Standard. **74-043** \$80.00



16 mm ($\frac{5}{8}$ ") cutting blade;
 9 cm ($3\frac{1}{2}$ ") long

Skull Breaker, Virchow, economy. **33-1011** \$42.40

Scalpels



Scalpels

Disposable, economy, sterile, gray plastic handle	<u>bx/10</u>	<u>cs/50 bx</u>
★ #10 carbon steel blade. 9410	\$11.65	516.00
#11 carbon steel blade. 9411	\$11.65	516.00
#22 carbon steel blade. 9422	\$11.65	516.00



Miniature, plastic thumb handle, nonsterile, England
 #11P stainless steel blade, each. **2201111** \$1.35




Smooth Handle



Fluted Handle

Reusable, one-piece stainless steel, $1\frac{1}{2}$ " blade
 Smooth Handle, 15.5 cm (6"). **14-350** \$3.15
 Fluted Handle, narrow blade, 15.5 cm (6"). **14-360** \$3.15



Scalpel Blades, individually wrapped

Carbon Steel, Economy, Bx/100	\$23.20
#10 blade. 9310	#21 blade. 9321
#11 blade. 9311	#22 blade. 9322
#15 blade. 9315	#23 blade. 9323
#20 blade. 9320	#24 blade. 9324
Carbon Steel, Germany, nonsterile	
#60 blade, 25 strips of 4, Bx/100. 2200260	\$59.00
#70 blade, 25 strips of 4, Bx/100. 2200270	\$70.00
Carbon Steel, Swann-Morton, nonsterile	
#60 blade, Bx/50. 22002SM60	\$39.60
Stainless Steel, Economy, Bx/100	\$21.95
#10 blade. 36-010S	#21 blade. 052-610
#11 blade. 052-602	#22 blade. 052-612
#20 blade. 36-020S	
Stainless Steel, Lance, Bx/100	\$47.50
#11 blade. 92111	#15 blade. 92115
Stainless steel, Swann-Morton, sterile	
#22 blade, Bx/100. 22001SM22	\$52.00
(RS-9861-19)	

Continued



Safely remove and store up to 100 used blades in this disposable box.

Scalpel Blade Remover. 22048..... \$4.35



Scalpel Handles, reusable

Stainless Steel, solid

- #3 Style, accepts blades #9-16
 - Economy. **02-490** \$3.45
 - Standard. **10-305**..... \$9.25↓
- #4 Style, accepts blades #20-26
 - Economy. **02-492** \$3.45
 - Standard. **10-310**..... \$9.25↓
- #7 Style, accepts blades #9-16
 - Economy. **02-494** \$5.35
 - Standard. **10-315**..... \$18.50

Stainless Steel, Hollow

- #8 Style, accepts blades #60 & #70. **HAN8** \$11.00
- PVC Coated Stainless Steel. *Not autoclavable*
- #3 Style, accepts blades #10-17. **KNF-225.00** \$6.00
- Plastic, contoured to fit palm of hand
- #5 Style, accepts blades #10-17. **KNF-221.00** \$5.25
- #6 Style, accepts blades #20-26. **1605** \$5.00
- #8 Style, accepts blades #60 & #70
- Bard-Parker™. **371080**..... \$44.95

Scissors




Barraquer


Scissors, Spring, ambidextrous, stainless steel

Barraquer (butterfly), 6 cm long, economy. **3737** \$34.50
10 mm angled-on-flat blades

36205/5-855



M3720



Castroviejo, straight, sharp tips, economy

- 7 cm (2¾"), 5 mm blades. **36205** \$45.00
Effective cutting edge is 3 mm; tip-to-tip is 3 mm.
- 8.5 cm (3¼"), 5 mm blades. **5-855** \$45.00
Effective cutting edge is 4 mm; tip-to-tip is 3 mm.
- 11.5 cm (4½"), 17 mm blades. **M3720** \$26.00
Effective cutting edge is 12 mm.



Uniband™

Lifetime guarantee against rust and defects.

Featherlite™



Uniband™, sharp tips, Laschal Surgical, USA

- Straight, 12 mm Blades. **LA-4** \$32.00
- Straight, 20 mm Blades. **LA-1** \$32.00
- Straight, 20 mm Serrated Blades. **LA-2** \$32.00
- Curved, 12 mm Blades. **LA-4C** \$32.00
- Curved, 20 mm Blades. **LA-1C** \$32.00

Featherlite™, Castroviejo, sharp tips, Laschal Surgical, USA

- Greater tactile sensitivity than Uniband™ scissors.
- Straight, 12 mm Blades. **MPF-4** \$80.80
- Curved, 12 mm blades. **MPF-4C** \$80.80

Scissors, Bonn, 15 mm blades, 9 cm (3½"), Germany

- Straight, Roboz. **RS-5840** \$107.00
- Curved, Roboz. **RS-5841** \$107.00

Straight

Curved



Scissors, Delicate, stainless steel

- Straight, sharp tips
 - 9 cm (3½"), 20 mm blades, Roboz. **RS-5880** \$102.00
 - 10.5 cm (4"), 22 mm blades, Standard. **13-003** \$21.30
- Curved, sharp tips
 - 9 cm (3½"), 20 mm blades, Roboz. **RS-5881** \$106.00
 - 10 cm (4"), 22 mm blades, Standard. **13-059** \$22.45

Continued