



Safely remove and store up to 100 used blades in this disposable box.

Scalpel Blade Remover. 22048..... \$4.35



Scalpel Handles, reusable

Stainless Steel, solid

- #3 Style, accepts blades #9-16
 - Economy. **02-490** \$3.45
 - Standard. **10-305**..... \$9.25↓
- #4 Style, accepts blades #20-26
 - Economy. **02-492** \$3.45
 - Standard. **10-310**..... \$9.25↓
- #7 Style, accepts blades #9-16
 - Economy. **02-494** \$5.35
 - Standard. **10-315**..... \$18.50

Stainless Steel, Hollow

- #8 Style, accepts blades #60 & #70. **HAN8** \$11.00
- PVC Coated Stainless Steel. *Not autoclavable*
- #3 Style, accepts blades #10-17. **KNF-225.00** \$6.00
- Plastic, contoured to fit palm of hand
- #5 Style, accepts blades #10-17. **KNF-221.00** \$5.25
- #6 Style, accepts blades #20-26. **1605** \$5.00
- #8 Style, accepts blades #60 & #70
- Bard-Parker™. **371080**..... \$44.95

Scissors




Barraquer


Scissors, Spring, ambidextrous, stainless steel

Barraquer (butterfly), 6 cm long, economy. **3737** \$34.50
10 mm angled-on-flat blades

36205/5-855



M3720



Castroviejo, straight, sharp tips, economy

- 7 cm (2¾"), 5 mm blades. **36205** \$45.00
Effective cutting edge is 3 mm; tip-to-tip is 3 mm.
- 8.5 cm (3¼"), 5 mm blades. **5-855** \$45.00
Effective cutting edge is 4 mm; tip-to-tip is 3 mm.
- 11.5 cm (4½"), 17 mm blades. **M3720** \$26.00
Effective cutting edge is 12 mm.



Uniband™

Lifetime guarantee against rust and defects.

Featherlite™



Uniband™, sharp tips, Laschal Surgical, USA

- Straight, 12 mm Blades. **LA-4** \$32.00
- Straight, 20 mm Blades. **LA-1** \$32.00
- Straight, 20 mm Serrated Blades. **LA-2** \$32.00
- Curved, 12 mm Blades. **LA-4C** \$32.00
- Curved, 20 mm Blades. **LA-1C** \$32.00

Featherlite™, Castroviejo, sharp tips, Laschal Surgical, USA

- Greater tactile sensitivity than Uniband™ scissors.
- Straight, 12 mm Blades. **MPF-4** \$80.80
- Curved, 12 mm blades. **MPF-4C** \$80.80

Scissors, Bonn, 15 mm blades, 9 cm (3½"), Germany

- Straight, Roboz. **RS-5840** \$107.00
- Curved, Roboz. **RS-5841** \$107.00

Straight


Curved



Scissors, Delicate, stainless steel

- Straight, sharp tips
 - 9 cm (3½"), 20 mm blades, Roboz. **RS-5880** \$102.00
 - 10.5 cm (4"), 22 mm blades, Standard. **13-003** \$21.30
- Curved, sharp tips
 - 9 cm (3½"), 20 mm blades, Roboz. **RS-5881** \$106.00
 - 10 cm (4"), 22 mm blades, Standard. **13-059** \$22.45

Continued



Scissors, Dissecting, 11.5 cm (4½"), stainless steel

Straight, sharp tips

25 mm Blades, Economy. 310-040	\$4.95
25 mm Blades, Standard. 13-002	\$19.00
24 mm Blades, Roboz. RS-5912	\$91.10
35 mm Blades, Economy. 310-045	\$4.80

Curved, sharp tips

25 mm Blades, Economy. 310-035	\$4.80
25 mm Blades, Standard. 13-010	\$19.40
24 mm Blades, Roboz. RS-5913	\$94.05




Angled (angled-on-side), stainless steel

11.5 cm (4½"), Economy. 310-066	\$6.50
12 cm (4¾"), 30 mm blades, Roboz. RS-6708	\$99.95



Scissors, Deaver, 14 cm (5½"), stainless steel

Straight, S/S (sharp/sharp) tips, Standard. 10-632	\$21.80
Straight, S/B (sharp/blunt) tips, Standard. 10-630	\$21.30




Scissors, Mayo, semi-blunt tips, stainless steel

Straight

14 cm (5½"), Economy. 310-086	\$9.90
14.5 cm (5¾"), Standard. 15-575	\$22.55
14.5 cm (5¾"), Roboz. RS-6870	\$87.10
17.5 cm (6¾"), Economy. 310-088	\$12.15
17.5 cm (6¾"), Standard. 15-579	\$26.95
17.5 cm (6¾"), Standard-TC. 16-854	\$58.50
23 cm (9"), Standard. 15-585	\$42.35

Curved

14 cm (5½"), Economy. 310-087	\$10.40
14.5 cm (5¾"), Standard. 15-576	\$22.95
14 cm (5½"), Roboz. RS-6871	\$87.10



Scissors, Metzenbaum, stainless steel

Straight

15 cm (5¾"), Economy. 12-225	\$14.00
---	---------

Curved

14.5 cm (5¾"), Economy. 12-230	\$14.00
14.5 cm (5¾"), Standard. 15-552	\$23.05
14 cm (5½"), German. 105-569	\$46.20

Blunt/Blunt Tips 

Sharp/Blunt Tips 

Sharp/Sharp Tips 

Scissors, Operating, stainless steel

Straight, S/S (sharp/sharp) tips

11.5 cm (4½"), Standard. 10-456	\$17.85
14 cm (5½"), Economy. 12-055	\$8.55
14 cm (5½"), LR (large rings), Economy. 11-107SSLR	\$6.05
14 cm (5½"), Standard. 10-552	\$21.30
14.5 cm (5¾"), German. 101-550	\$36.50
15.5 cm (6"), Economy. 12-060	\$10.75
16.5 cm (6½"), Economy. 12-065	\$10.75
16.5 cm (6½"), Standard. 10-602	\$24.75
16.5 cm (6½"), German. 101-600	\$40.20

Straight, S/B (sharp/blunt) tips

11.5 cm (4½"), Economy. 12-015	\$8.55
12 cm (4¾"), Standard. 10-453	\$17.85
12 cm (4¾"), German. 100-453	\$33.60
12 cm (4¾"), 30 mm blades, Roboz. RS-6700	\$99.95
14 cm (5½"), Economy. 12-025	\$8.55
14 cm (5½"), Standard. 10-550	\$21.30
14 cm (5½"), German. 100-552	\$36.50
16.5 cm (6½"), Economy. 310-111	\$11.00
16.5 cm (6½"), Standard. 10-600	\$24.75
16.5 cm (6½"), 60 mm blades, Roboz. RS-6824	\$76.25

Straight, B/B (blunt/blunt) tips

11.5 cm (4½"), Economy. 310-100	\$7.40
14 cm (5½"), Economy. 12-085	\$8.55
14 cm (5½"), Standard. 10-558	\$21.30
14 cm (5½"), German. 102-558	\$36.50
16.5 cm (6½"), Economy. 12-095	\$10.75
16.5 cm (6½"), Standard. 10-608	\$24.75
16.5 cm (6½"), 60 mm blades, Roboz. RS-6828	\$76.25

Curved, S/S (sharp/sharp) tips


14 cm (5½"), Economy. 12-155	\$8.55
14.5 cm (5¾"), Standard. 10-566	\$13.25
16.5 cm (6½"), Economy. 12-165	\$10.75
16.5 cm (6½"), Standard. 10-616	\$25.25

Scissors, Operating (cont.)


- Curved, S/B (sharp/blunt) tips
- 11.5 cm (4½"), Economy. **310-101** \$7.65
 - 14 cm (5½"), Economy. **12-125** \$8.55
 - 14.5 cm (5¾"), Standard. **10-560** \$15.40 ↓
 - 14 cm (5½"), German. **100-560** \$37.90
 - 16.5 cm (6½"), Economy. **12-135** \$10.75
 - 16.5 cm (6½"), Standard. **10-610** \$25.25
 - 16.5 cm (6½"), 60 mm blades, Roboz. **RS-6849** \$81.20
- Curved, B/B (blunt/blunt) tips
- 11.5 cm (4½"), Economy. **310-102** \$7.85
 - 14 cm (5½"), Economy. **12-185** \$8.55
 - 14 cm (5½"), Standard. **10-568** \$21.80
 - 16.5 cm (6½"), Economy. **12-195** \$10.75
 - 16.5 cm (6½"), Standard. **10-618** \$25.25



- Scissors, Stevens Tenotomy, 11.5 cm (4½"), stainless steel**
 Straight, blunt tips, Standard. **13-520** \$24.75



- Scissors, Strabismus, stainless steel**
- Straight
- 10 cm (4"), Economy. **26-105** \$10.60
 - 11.5 cm (4½"), Standard. **13-524** \$23.65
- Curved
- 11 cm (4¼"), Economy. **26-110** \$10.60
 - 11.5 cm (4½"), Standard. **13-526** \$24.15



- Scissors, Bandage, stainless steel, economy**
- Knowles, straight, sharp probe tips, 14 cm (5½"). **1070** \$6.25
 - Lister, angled, blunt probe tips, 14 cm (5½"). **1050** \$5.00
- Can be used as enterotomy scissors.



Cut through everything from wire to bone with these durable scissors. Bottom blade is serrated and has a footplate for safety. Blades are stainless steel & handles are plastic. Autoclavable to 120° C.

- Utility Shears, 18 cm (7¼"), economy. 310-341** \$5.70

Scrapers

Remove paraffin, ice, etc., safely from your laboratory equipment with these plastic scrapers.



- Scrapers, plastic**
- 5 cm (2") Wide Blade. **311-245** \$3.70
 - 10 cm (4") Wide Blade. **311-250** \$4.70
 - 15 cm (6") Wide Blade. **311-255** \$5.70


Spatulas & Scoops



- Spatulas, Hardwood Handle, stainless steel blade**
- 10 cm (4") Blade. **14-365B** \$5.40
 - 15 cm (6") Blade. **14-365D** \$7.05
 - 20 cm (8") Blade. **20-497** \$6.05



- Spatulas, Heavy-Duty, stainless steel**
 Round End/Square End, both flat
- 30.5 cm (12"). **320-425** \$16.10
 - 46 cm (18"). **320-430** \$21.50



- Spatulas, Micro, stainless steel**
- Square End/Spoon End
- 19 cm (7½"). **320-411** \$4.50
 - 23 cm (9"). **320-414** \$4.65
- Round End/Tapered End, 18 cm (7") long. **320-370** \$5.00
 Round end: 8 x 51 mm; Tapered end: 6 x 51 mm
- V-Shaped End/Round End, 16 cm (6¼") long. **HS15906** \$7.45



- V-Shaped (Reagent Digger), 16.5 cm (6½"). **600LSCSS6** \$2.60
 Handle, wooden, 13 cm long, each. **600LSHND1** \$3.00