

Dissecting Kits

Inquire about customizing a kit for your surgical procedure.
Kit components are discounted 10%.

Basic kit. 38501 \$5.30

Contents: Dissecting Case, 6-loop (1257370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

Mouse kit. 38551 \$17.85

Contents: Dissecting Case, 6-loop (1257370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

Dissecting Pans, stainless steel (without wax)

20⁷/₈ x 12¹³/₁₆ x 2¹/₂". **AM20029** \$21.90
25³/₄ x 17³/₄ x 2¹/₄". **AM68357** \$64.90



Insect Pins, 38 mm (1½") long, nylon head, Czech Republic

Black Enameled, Pk/100 \$7.00
#000, 0.25 mm dia. **355425** #1, 0.40 mm dia. **355440**
#00, 0.30 mm dia. **355430** #2, 0.45 mm dia. **355445**
#0, 0.35 mm dia. **355435** #3, 0.50 mm dia. **355450**
Stainless Steel, Pk/100 \$8.75
#000, 0.25 mm dia. **355525** #1, 0.40 mm dia. **355540**
#00, 0.30 mm dia. **355530** #2, 0.45 mm dia. **355545**
#0, 0.35 mm dia. **355535** #3, 0.50 mm dia. **355550**



Instrument Cup, silicone, 3½" opening, 2½" deep. 150-1005. \$57.00



Instrument Storage Boxes, 10 x 4 x 2", 18-8 stainless steel

Small, 7 x 3 x 1½". **35-406** \$26.60
Large, 10 x 4 x 2". **35-408** \$35.00

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible
Transparent. **4461** \$0.95
White. **4460** \$0.55
Stainless Steel, 12 mm wide, **USA. 605-065** \$6.25



Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.



Sterilizer, Hot Bead, Germinator™ 500, USA. 5-1450 \$426.00
Glass Beads (replacement), 1.5 mm, Bg/12 oz. **5-1455** \$30.00



Safely incinerate microorganisms in 5-7 seconds at 815°C (1500°F) within this ceramic, stainless steel-encased chamber. Compact and flameless operation allows this unit to be used in anaerobic and aerobic chambers. Ideal for sterilizing loops and needles, glass tube/pipette mouths, and various glass and metal instruments. 1-year warranty.

Sterilizer, Infrared Heat, BactiZapper™. B1000 \$476.00

6100A



Sterilizing Trays, plastic, autoclavable, PST, USA

B/L/M (base, lid & mat) included.
6 x 2½ x ¾" (15.2 x 6.4 x 1.9 cm). **2600A** \$23.00
(500252)
★ 6½ x 4 x ¾" (16.2 x 10.2 x 1.9 cm). **4600A** \$33.00
10 x 6 x ¾" (25.4 x 15.2 x 1.9 cm). **6100A** \$67.00
(500254)



Instrument Trays

Polypropylene, deep w/flat cover, autoclavable
8½ x 6 x 2¾" (22 x 15 x 7 cm), Pk/2. **60081736** \$25.15

Continued



15F



35-309

Stainless Steel

Shallow (Mayo), **USA**

15¹/₈ x 10⁵/₈ x 3³/₄" (38.4 x 27 x 1.9 cm). **15F** \$51.65

Deep w/ Dome Cover

8³/₄ x 3¹/₄ x 1¹/₂" (22.2 x 8.4 x 3.8 cm). **35-309** \$39.90

Deep – Tray Only, **USA**

8⁷/₈ x 5 x 2" (22.7 x 12.9 x 5.1 cm). **35-314NP** \$41.40



This natural Arkansas stone outlasts and produces a better finish than a man-made stone. Mounted in a wooden box. Dimensions: 4 x 2 x 1/2"

Sharpening Stone, fine grit (hard). ABR-800.10 \$31.95

Includes a plastic cleaning basket.
Requires only 42 kHz of power.
Dimensions: 6 x 3 1/2 x 2"(D).



Ultrasonic Cleaner, 600 mL, mechanical timer. CLN-520.20 \$45.95

Cleaning Solution, 1 Quart. **CLN-850.10** \$9.95

Biodegradable concentrate mixes 1:30.