

Accessories

Slide & Coverslip Handling

Forceps


Spade Tips

- Angled (Kuehne), 11 cm (4¼")
 Stainless Steel, economy. **33-1041** \$4.15
 Teflon® Coated, Bel-Art. **379420000** \$35.60
 Straight, stainless steel, 10.5 cm, economy. **57.926** \$8.45



Square Point Tips, straight

- Stainless Steel, 12 cm (4¾"), economy. **1244** \$2.25
 Teflon® Coated, 11 cm (4¼"), Bel-Art. **379430000** \$39.40



Minimize contamination during slide & coverslip handling. Will not scratch; corrosion and heat resistant.

- Hemostat, Polypropylene, 12.5 cm. **HEM-300.00** \$1.95

Immersion Oil & Bottles



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities.
USA

Immersion Oil

Type A, low viscosity (150 cSt). Short W.D. optics

- 7.5 mL in LDPE dispensing bottle. **18101** \$10.00
 20 mL in HDPE dispensing bottle. **18102** \$15.35
 480 mL in amber bottle. **18104** \$110.00

Type B, high viscosity (1250 cSt). Long W.D. optics

- 7.5 mL in LDPE dispensing bottle. **18105** \$10.00
 20 mL in HDPE dispensing bottle. **18106** \$15.35
 480 mL in amber bottle. **18108** \$110.00

Type NVH, very high viscosity (21,000 cSt)

- 480 mL in amber bottle. **18113** \$67.00


Type HF, low fluorescence (700 cSt)

- 480 mL in amber bottle. **18118** \$169.00

Type FF, fluorescence free (170 cSt)

- 480 mL in amber bottle. **18124** \$163.00

A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



- Bottle, Dispensing, HDPE, 20 mL, each. **18145** \$8.50

Slide Storage

PLASTIC SLIDE BOXES



25-Slide Squeeze-Release Cover



100-Slide Steel-Hinged Cover

50-Slide Hinged Cover



Slide Boxes, hold 75 x 25 mm (3 x 1") slides

- 25-Slide, squeeze-release cover, cork bottom, 9.5 x 8.5 x 3 cm
 White, ABS plastic. **97-0025** \$2.50
 50-Slide, hinged cover, cork bottom, 20 x 9.5 x 3.5 cm
 Blue, ABS plastic. **97-0056** \$6.80
 Red, ABS plastic. **97-0058** \$6.80
 100-Slide, steel-hinged cover, cork bottom
 Blue, ABS plastic. **97-0126** \$8.90
 Red, ABS plastic. **97-0128** \$8.90
 White, ABS plastic. **97-0125** \$8.90



Store 4500 75 x 25 mm slides in 6 divided drawers. Each cabinet comes with a sample of 6 compressor follow blocks and 6 drawer labels. Safely stack up to 12 cabinet units on a single base.
 18¾"(D) x 15¼"(W) x 5"(H)

Slide Cabinets, steel

- 6-Drawer, stackable
 Green. **513500G** \$255.00
 Tan. **513500T** \$255.00
 White. **513500W** \$255.00
 Memorandum Markers, paperboard, Pk/100. **C108** \$12.75
 Compressor Follow Block, rubber, each. **C109** \$1.90
 Base, raises cabinets 4¼" off floor
 Green. **513503G** \$169.00
 Tan. **513503T** \$169.00
 White. **513503W** \$169.00



Cardboard Folder



Plastic Folder

Slide Folders, 20-slide, hold 75 x 25 mm (3 x 1") slides

- Cardboard, black tape, economy. **HS9902** \$6.35 (12-587-10)
 Plastic, gray tray w/2 clear doors. **M750-20GY** \$11.60
 Transparent covers hold visible slides securely in place, but can fold flat against back of tray when opened. Slides are slightly elevated for easy retrieval.

Beakers

Polypropylene Beakers



Beakers

Glass (borosilicate), low form, double scale

Economy

10 mL, no graduations, Pk/12 600BG10	\$35.50
50 mL, Pk/12. 600BG50	\$41.10
100 mL, Pk/12. 600BG100	\$43.15
150 mL, Pk/12. 600BG150	\$46.15
250 mL, Pk/12. 600BG250	\$49.85
400 mL, Pk/12. 600BG400	\$62.75
600 mL, Pk/6. 600BG600	\$45.80
800 mL, Pk/6. 600BG800	\$77.75
1000 mL, Pk/6. 600BG1000	\$74.00
2000 mL, each. 600BG2000	\$23.20
4000 mL, each. 600BG4000	\$69.75

Heavy-Duty

50 mL, Pk/12. 600BG100350	\$44.55
100 mL, Pk/12. 600BG1003100	\$75.50
150 mL, Pk/12. 600BG1003150	\$88.50
250 mL, Pk/12. 600BG1003250	\$102.00
400 mL, Pk/12. 600BG1003400	\$116.00
600 mL, Pk/6. 600BG1003600	\$88.75
1000 mL, each. 600BG100310000	\$20.65

Polypropylene, premium, low form, single scale

100 mL, 25 mL graduations, each. 262110000	\$4.05
250 mL, 50 mL graduations, each. 262130000	\$6.75
600 mL, 100 mL graduations, each. 262150000	\$10.50
1000 mL, 200 mL graduations, each. 262160000	\$10.60
2000 mL, 500 mL graduations, each. 262170000	\$32.15
4000 mL, 1000 mL graduations, each. 262190000	\$53.60

Bench Protector



Protect your bench surface using this 40 mil absorbent liner backed with 0.5 mil polyethylene film.

Bench Liner, Roll, 50' x 20", white, Labmat™. **246750000**..... \$95.45

Bottles




Deliver controlled drops from standard 0.030" tip opening with these LDPE bottles. Sealer cap snaps on easily.

Bottles, Dispensing, squeeze bottles

LDPE (low-density polyethylene), Yorker cap

30 mL (1 oz), Bg/12. 116370001	\$33.15
185 mL (6 oz), Bg/12. 116370006	\$63.85

A long nozzle, textured sides, and slim profile make this chemically inert, HDPE container ideal for dispensing immersion oil.



HDPE (high-density polyethylene)

20 mL, each. 18145	\$8.50
---------------------------------	--------

Amber glass protects light-sensitive liquids. Screw-cap dropper features a black rubber bulb and glass pipet.



Boston Round

Bottles, Dropping, Amber Glass (flint)

30 mL (1 oz), Bx/12. 3104053	\$19.30
Bottle Dropper Assembly (replacement) 20/400, 7 x 76 mm Pipet, Pk/12. 3104053D	\$11.10
60 mL (2 oz), Bx/12. 3104054	\$27.00
Bottle Dropper Assembly (replacement) 20/400, 7 x 89 mm Pipet, Pk/12. 3104054D	\$11.15

French Square



Bottles, Dropping, Clear Glass (flint)

15 mL (1/2 oz), Bx/12. 3104062	\$33.80
Bottle Dropper Assembly (replacement) 20/400, 7 x 60 mm Pipet, Pk/12. 3104074	\$11.55
30 mL (1 oz), Bx/12. 3104063	\$39.20
Bottle Dropper Assembly (replacement) 24/400, 7 x 75 mm Pipet, Pk/12. 3104075	\$14.75

Continued