


## Staining



Ten slide compartments, separated by 1/2" of space, are additionally isolated by eight barrier dividers when the lid is closed. Individual slides are suspended 1/2" above the chamber floor by four pedestal posts and four corner posts.

### Staining/Moisture Chambers

10-Slide, polystyrene, <b>USA</b>	<u>ea</u>	<u>pk/3</u>
Clear. <b>240-9000-010</b> .....	\$54.95	154.00
Black. <b>240-9020-Z10</b> .....	\$56.10	157.00



20-Slide, Simport

★ Clear Lid (PETG), black base (ABS), each. <b>M920-1</b> .....	\$195.00
★ Black Lid (ABS), black base (ABS), each. <b>M920-2</b> .....	\$195.00



### Staining Dishes


Glass (flint)

20-Slide, holds slides horizontally in metal rack	
Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H)	
Set (dish, rack & lid). <b>SDS-20</b> .....	\$31.85
Dish & Lid. <b>SD-20</b> .....	\$19.40
Rack only, stainless steel. <b>SR-20</b> .....	\$16.00
30-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H)	
Set (dish, rack & lid). <b>SDS-30</b> .....	\$43.50
Dish & Lid. <b>SD-30</b> .....	\$21.70
Rack only, stainless steel. <b>SR-30</b> .....	\$26.65
60-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H)	
Set (dish, rack & lid). <b>SDS-60</b> .....	\$60.80
Dish & Lid. <b>SD-60</b> .....	\$33.20
Rack only, stainless steel. <b>SR-60</b> .....	\$34.36



Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

Plastic (polyoxymethylene [POM]), holds 25 slides	
Set (dish, rack & lid). <b>0990-99</b> .....	\$50.90



### Staining Jars, hold slides (up to 76 x 26 mm) vertically

Coplin, holds 10 slides back-to-back	<u>ea</u>	<u>cs/12</u>
Glass w/Ground Glass Cover, <b>USA. 673</b> .....	\$13.45	149.00
Glass Cover, Fits 673 Jar. <b>673-LID</b> .....	\$4.95	—
Glass w/Screw-On Phenolic Cap, <b>USA. 674</b> .....	\$13.45	149.00
Phenolic Cap, Fits 674 Jar, black. <b>674-LID</b> ...	\$5.05	—
PP w/Screw-On Cap. <b>60062101</b> .....	\$5.05	49.20
Polypropylene jar is unbreakable, microwave proof, stain resistant, and chemically inert. Cap is domed.		
Wide Top, Glass w/Square Lid, <b>USA. 675</b> .....	\$18.00	197.00
Holds 16 slides back-to-back		
Glass Lid, Fits 675 Jar. <b>675-LID</b> .....	\$5.50	—



Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

### Wide Top, Plastic (POM)

12-Slide, EasyDip™, Simport	
Jar only, Blue, Cs/6. <b>M900-12B</b> .....	\$67.65
Jar only, Green, Cs/6. <b>M900-12G</b> .....	\$67.65
Jar only, Pink, Cs/6. <b>M900-12P</b> .....	\$67.65
Jar only, White, Cs/6. <b>M900-12W</b> .....	\$67.65
Jar only, Yellow, Cs/6. <b>M900-12Y</b> .....	\$67.65
Rack/Lid only, gray, Cs/6. <b>M905-12DGY</b> .....	\$74.90



Anodized aluminum staining jar holder

(Jars sold separately)

Staining Jar Holder, 6-Place. <b>M906</b> .....	\$82.90
---	---------

Continued




- 24-Slide, 4 x 2 x 4", PP (microwaveable), light green  
 Jar & Lid. **4501030** ..... \$11.65  
 Rack only, gray. **4501042** ..... \$21.45  
 Staining Jar Station, 12-Place. **4501011** ..... \$221.00  
 Includes: metal frame, (12) 24-slide dishes  
 w/ lids (4501030) & (1) rack (4501042).  
 Dimensions: 29(L) x 5(W) x 5"(H)


**Tissue (Embedding) & Biopsy Cassettes**



- Base Molds**, Plastic, disposable, Simport  
 15 x 15 x 5 mm, Cs/1000. **M475-2** ..... \$156.00  
 24 x 24 x 5 mm, Cs/1000. **M475-3** ..... \$156.00  
 30 x 24 x 5 mm, Cs/1000. **M475-4** ..... \$156.00  
 37 x 24 x 5 mm, Cs/1000. **M475-5** ..... \$156.00  
 37 x 24 x 10 mm, Cs/500. **M475-10** ..... \$80.30  
 Designed for Simport Macrosette™.




Solvent resistant,  
0.2 x 0.2 mm mesh.




- Biopsy Bags**  
 Nylon, 30 x 45 mm, 5 bags/200, Pk/1000. **4502011** ..... \$235.00  
 Paper (Tea Bags), 2¼ x 2¼", Pk/500. **56100** ..... \$46.50

Made of 60 PPI cellular  
 polyester urethane foam.  
 Use rectangles with cassettes,  
 circles with capsules.



- Biopsy Pads**  
 Rectangular  
 Black Foam, 30.2 x 25.4 mm, Bg/1000. **339014** ..... \$42.00  
 Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. **M476-1** ..... \$222.00  
 Circular, Simport  
 Blue Foam, 25.4 x 2.7 mm, Cs/10,000. **M476-2** ..... \$273.00

**Biopsy Cassettes**



- Biopsy Cassettes, Attached Lids**, acetyl polymer, Simport  
 Histosette™ I, 1 mm square pores, Pk/1500 ..... \$154.00  
 30° front labeling area.  
 Blue. **M491-6**                      Pink. **M491-3**  
 Green. **M491-4**                      White. **M491-2**
- Uniset™, hinged lids, Pk/1500 ..... \$165.00  
 35° front labeling area. (Alternative to Sakura's Uni-Cassette®)  
 Blue. **M506-6**                      Pink. **M506-3**  
 Green. **M506-4**                      White. **M506-2**
- Micromesh™, 0.38 mm pores, Pk/1000 ..... \$251.00  
 45° front labeling area.  
 Blue. **M507-6**                      Pink. **M507-3**                      Yellow. **M507-5**  
 Green. **M507-4**                      White. **M507-2**

**Tissue Cassettes**



- Tissue Cassettes, Detached Lids**, slotted pores  
 Histosette™ II, Simport, Pk/1000 ..... \$107.00  
 45° front labeling area.  
 Blue. **M485-6**                      Pink. **M485-3**  
 Green. **M485-4**                      White. **M485-2**
- Uni-Cassettes™, **USA**, Pk/500 ..... \$77.00  
 Aqua. **6059**                      Lime. **6062**  
 Blue. **6054**                      Peach. **6064**  
 Burnt Orange. **6056**                      Pink. **6052**  
 Chartreuse. **6061**                      Red. **6060**  
 Gold. **6058**                      Rose. **6063**  
 Green. **6053**                      Tan. **6055**  
 Lavender. **6065**                      White. **6050**  
 Lilac. **6057**                      Yellow. **6051**

Continued