



Simport



Simport
Macrosette

Tissue Cassettes, Attached Lids, slotted pores

Histosette™ I, snap-apart lids, Pk/1500..... \$154.00
30° front labeling area.

Blue. **M490-6** Pink. **M490-3**
Green. **M490-4** White. **M490-2**

Histosette™ II, integral lids, Pk/1500..... \$204.00
45° front labeling area.

Blue. **M492-6** Pink. **M492-3**
Green. **M492-4** White. **M492-2**

Unisette™, hinged lids, Pk/1500 \$165.00
35° front labeling area. (Alternative to Sakura's Uni-Cassette®)

Blue. **M505-6** Pink. **M505-3**
Green. **M505-4** White. **M505-2**

Macrosette™, 2X normal depth, white, Pk/750. **M512** \$221.00
Dimensions: 40.1 x 28.5 x 13 mm(H)

Cassette Storage



Matchbox Style

Box, Matchbox Style (Pill Box), cardboard, white

2 1/8 x 1 1/8 x 5/8", Pk/72. **20047** \$35.45

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. 8342 \$121.00

Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks
in 4 drawers, each divided into 12 rows.
Cabinets can be stacked up to 8 high.
Dimensions: 49.5 x 38 x 21 cm
(19 1/2 x 15 x 8 1/4")



Tissue Cassette Cabinet, Steel, 4-drawer. 485400 \$386.00

Store up to 2100 cassettes in 6 drawers.
Each drawer includes sponge block
spacers and an identification label.
Cabinets can be stacked up to 4 high.
Dimensions: 44.5 x 24 x 39.5 cm(H)
(17 1/2 x 9 1/2 x 15 1/2")



Tissue Cassette Cabinet, ABS. 97-0600 \$204.00



Simport


Store 250 cassettes or 165 embedding
rings in each stackable drawer.
High-impact polystyrene resists
temperatures from -80°C to +80°C.
Drawer dimensions: 40.5 x 23 x 5.1 cm.

Tissue Cassette Drawers, Polystyrene

Pk/2. **M495-2** \$77.55

Pk/6. **M495-6** \$176.00

Brushes



LB-37



LB-30

Brushes

Flat Bristles (paint style)

Camel Hair, 1" wide, 6 1/2" total length. **BR665B (LB-37)** \$6.65

Round Bristles, camel hair

3/4" Diameter, 8 1/2" total length. **BR661 (LB-30)** \$11.15

Borosilicate glass bottle;
blue PP cap (GL45)



Bottles, Media, graduated, autoclavable

100 mL, Cs/10. 600BM0100	\$96.85
250 mL, Cs/10. 600BM0250	\$110.00
500 mL, Cs/10. 600BM0500	\$125.00
1000 mL, Cs/10. 600BM1000	\$178.00
Replacement Items	
Bottle Caps, GL45, blue, Pk/10. B3000-CAP	\$47.60
Bottle Drip-Free Rings, GL45, Pk/10. B3000-RIN	\$25.65



Globally Harmonized System (GHS) wash bottles are vented to prevent solvent drips. Polyethylene bottles (LDPE) with polypropylene (PP) caps.

Bottles, Wash, GHS Labeled, 500 mL (empty)

Acetone, Pk/6. 600370051	\$41.30
Bleach, Pk/6. 600370059	\$41.30
Deionized Water, Pk/6. 600370064	\$41.30
Ethanol, Pk/6. 600370052	\$41.30
Methanol, Pk/6. 600370053	\$41.30
Water, Pk/6. 600370061	\$41.30

Bottles, Wash, Safety Labeled, 500 mL (empty)

4-Color diamond, LDPE bottle/PP cap	
Soap, Pk/4. 117160014	\$41.30
Universal (label your own bottle), Pk/4. 117160009	\$41.30



Bottles, Wash, Integral Spout, LDPE bottle/PP cap

250 mL (8 oz), Cs/10. 601633	\$76.50
500 mL (16 oz), Cs/10. 601634	\$87.45



Bottles, Narrow-Mouth, HDPE bottle/polypropylene cap

125 mL, Pk/12. 106200015	\$25.45
250 mL, Pk/12. 106200016	\$36.70
500 mL, 28 mm cap, Pk/12. 106200007	\$42.95
1000 mL, Pk/6. 106200008	\$33.60



Bottles, Wide-Mouth

HDPE bottle/polypropylene cap	
125 mL, 38 mm cap, Pk/12. 60033407	\$25.45
250 mL, Pk/12. 106250006	\$36.70
500 mL, 53 mm cap, Pk/12. 60033409	\$42.75
1000 mL, Pk/6. 60033410	\$35.95
Polypropylene bottle and cap, autoclavable	
125 mL, 38 mm cap, Pk/12. 60033307	\$25.15
250 mL, 43 mm cap, Pk/12. 60033308	\$36.70
500 mL, 53 mm cap, Pk/12. 106320007	\$42.45
1000 mL, 63 mm cap, Pk/6. 60033310	\$37.40



Bottle Top Dispensers

0.5-5 mL. 600BTDR2	\$312.00
5-60 mL. 600BTDR5	\$372.00

Boxes



Box, Matchbox Style (Pill Box), cardboard, white

2 1/8 x 1 1/8 x 5/8", Pk/72. 20047	\$35.45
---	---------

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. **8342**

Each cardboard box can store approximately 540 paraffin blocks.



Boxes, Glass Disposal, cardboard, polyethylene bag included

Bench Model, 8 x 8 x 10", Pk/6. 246530002	\$122.00
Floor Model, 12 x 12 x 27", Pk/6. 246530001	\$163.00
Glass Disposal Box Cover, reusable, PP. 132040001	\$53.65
<i>Replacement item for the floor model. Nonrefundable.</i>	
<i>Polypropylene cover is autoclavable.</i>	

Continued