



Simport



Simport
Macrosette

Tissue Cassettes, Attached Lids, slotted pores

Histosette™ I, snap-apart lids, Pk/1500..... \$154.00
30° front labeling area.

Blue. **M490-6** Pink. **M490-3**
Green. **M490-4** White. **M490-2**

Histosette™ II, integral lids, Pk/1500..... \$204.00
45° front labeling area.

Blue. **M492-6** Pink. **M492-3**
Green. **M492-4** White. **M492-2**

Unisette™, hinged lids, Pk/1500 \$165.00
35° front labeling area. (Alternative to Sakura's Uni-Cassette®)

Blue. **M505-6** Pink. **M505-3**
Green. **M505-4** White. **M505-2**

Macrosette™, 2X normal depth, white, Pk/750. **M512** \$221.00
Dimensions: 40.1 x 28.5 x 13 mm(H)

Cassette Storage



Matchbox Style

Box, Matchbox Style (Pill Box), cardboard, white

2 1/8 x 1 1/8 x 5/8", Pk/72. **20047** \$35.45

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. 8342 \$121.00

Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks
in 4 drawers, each divided into 12 rows.
Cabinets can be stacked up to 8 high.
Dimensions: 49.5 x 38 x 21 cm
(19 1/2 x 15 x 8 1/4")



Tissue Cassette Cabinet, Steel, 4-drawer. 485400 \$386.00

Store up to 2100 cassettes in 6 drawers.
Each drawer includes sponge block
spacers and an identification label.
Cabinets can be stacked up to 4 high.
Dimensions: 44.5 x 24 x 39.5 cm(H)
(17 1/2 x 9 1/2 x 15 1/2")



Tissue Cassette Cabinet, ABS. 97-0600 \$204.00



Simport


Store 250 cassettes or 165 embedding
rings in each stackable drawer.
High-impact polystyrene resists
temperatures from -80°C to +80°C.
Drawer dimensions: 40.5 x 23 x 5.1 cm.

Tissue Cassette Drawers, Polystyrene

Pk/2. **M495-2** \$77.55

Pk/6. **M495-6** \$176.00

Brushes



LB-37



LB-30

Brushes

Flat Bristles (paint style)

Camel Hair, 1" wide, 6 1/2" total length. **BR665B (LB-37)** \$6.65

Round Bristles, camel hair

3/4" Diameter, 8 1/2" total length. **BR661 (LB-30)** \$11.15

Hinged
Polystyrene



Boxes, Plastic, (see following chart)

Plastic Boxes				
Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
H-1	1 x 1 x ¼"	-	Reticles, Pointers	\$0.40
H-1/11B	1 x 1 x ¼"	1	Reticles, Pointers	\$0.50
H-54	2 9/16 x 2 9/16 x ½"	-	Filters	\$1.20
H-54/2916G	2 9/16 x 2 9/16 x ½"	1	Filters	\$1.35
H-28/334B2	3 ¾ x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40
H-10	2 7/8 x 1 3/16 x 3/8"	-	Small Parts	\$0.60
H-21/512B	5 ½ x 1 ¾ x ½"	1	Forceps & Scissors	\$1.80

Brushes



Brushes

- Bench, horse hair, 10" long. **BRS-969.10** \$9.95
Bristle head: 4¾"(L) x 1½"(W); bristle length: 1"
- Flat Bristles (paint style)
Camel Hair, 1" wide, 6½" total length. **BR665B (LB-37)** \$6.65
- Horsehair, heavy-duty bristles, metal handle. **33761880** \$0.50 ↓
- Round Bristles, camel hair, 8½" long. **BR661 (LB-30)** \$11.80
- Swivel-Head, black nylon bristles. **BR565 (LB-14)** \$5.00
Bristle head: 4¼"(L) x 1½"(W); total length: 12"

Burner



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing. Height is 15.6 cm (6 1/8")



H-5550 gridded top, 3.8 cm (1½") diameter, produces a uniform blue cone of light. Height is 21.6 cm (8½")

Burners, Bunsen, natural gas, Humboldt Mfg, USA

- Standard, flame stabilizer. **H-5870** \$32.65
- High Temperature, Meeker-type. **H-5550** \$98.80
- Accessory
Bunsen Burner Natural Gas Connector
3' Long. **H-2087** \$39.15



- Burner Lighter**, replaceable single flint, round file. **SOL-290.00** ... \$4.40
- Flints (replacement), Pk/5. **SOL-290.01** \$2.95

Carboys



Carboys with Spigot, polypropylene, round

- 10 L (2.5 gal). **60034051** \$84.45
- 20 L (5 gal). **60034052** \$117.00