

Drying Racks



188050000

189330023



HS23243A

Drying Racks, Benchtop, Bel-Art

Single-sided. **188050000**.....\$174.00

Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 43 x 17 x 48 cm(H).

Double-sided

2-Tier. **189330022**\$433.00

3-Tier. **189330023**\$589.00

Drying Rack, Benchtop/Wallmount, Heathrow

Single-sided. **HS23243A**\$148.00

Dry up to 90 pieces of labware with this self-standing or wall-mounted rack. Constructed of vinyl-coated wire hairpin loops with a polystyrene tray and rack. Dimensions: 46.2 x 18.2 x 52.5 cm(H).

Electrophoresis



Electrophoresis System, MyGel™ Mini. **E1101**\$649.00

Eyewear

Absorbs
99.9% UV



Glasses, Safety, integral side shield. **43574002**.....\$6.00



Flexible frame and adjustable elastic strap provide a tight seal.

Goggles, Safety, vented (chemical splash). **97-2001**.....\$4.40

Filter & Blotting Papers



Ahlstrom - Munksjo
Filter Paper

Qualitative, Standard Grades (for gravity filtration)

Slow flow (4 mL/min), 1.5 µm retention, Grade 610

(Whatman 5; S&S 602)

9 cm, Pk/100. **6100-0900**\$12.90

12.5 cm, Pk/100. **6100-1250**\$16.10

15 cm, Pk/100. **6100-1500**\$19.40

20.5 cm, Pk/100. **6100-2050**\$30.65

24 cm, Pk/100. **6100-2400**\$42.70

Medium flow (30 mL/min), 2 µm retention, Grade 642

(Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)

9 cm, Pk/100. **6420-0900**\$10.05

15 cm, Pk/100. **6420-1500**\$18.05

Medium flow (55 mL/min), 2.5 µm retention, Grade 601

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

4.25 cm, Pk/100. **6010-0425**\$7.70

9 cm, Pk/100. **6010-0900**\$10.35

12.5 cm, Pk/100. **6010-1250**\$14.85

15 cm, Pk/100. **6010-1500**\$18.10

18.5 cm, Pk/100. **6010-1850**\$26.25

20.5 cm, Pk/100. **6010-2050**\$28.30

33 cm, Pk/100. **6010-3300**\$76.40

Medium flow (60 mL/min), 6 µm retention, Grade 613

(Whatman Student; S&S 595, 0903)

9 cm, Pk/100. **6130-0900**\$8.80

15 cm, Pk/100. **6130-1500**\$15.35

20.5 cm, Pk/100. **6130-2050**\$23.00

33 cm, Pk/100. **6130-3300**\$51.15

Medium flow (200 mL/min), 10 µm retention, Grade 631

(Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)

20.5 cm, Pk/100. **6310-2050**\$28.50

33 cm, Pk/100. **6310-3300**\$77.20

Fast flow (235 mL/min), 25 µm retention, creped, Grade 615

9 cm, Pk/100. **6150-0900**\$9.10

12.5 cm, Pk/100. **6150-1250**\$13.20

15 cm, Pk/100. **6150-1500**\$15.85

20.5 cm, Pk/100. **6150-2050**\$23.70

33 cm, Pk/100. **6150-3300**\$52.70

Very Fast flow (360 mL/min), 35 µm retention, creped, Grade 617

9 cm, Pk/50. **6170-0900**\$8.30

12.5 cm, Pk/50. **6170-1250**\$11.05

15 cm, Pk/50. **6170-1500**\$12.85

20.5 cm, Pk/50. **6170-2050**\$19.60

33 cm, Pk/50. **6170-3300**\$42.00

Qualitative, Pleated (Fluted)

Slow flow (30 mL/min), 2 µm retention, Grade 542

(Whatman 2V)

12.5 cm, Pk/100. **5420-1250**\$82.45

18.5 cm, Pk/100. **5420-1850**\$93.30

Qualitative, Wet-Strengthened Grades (for vacuum filtration)

- Medium flow** (65 mL/min), 3 µm retention, **Grade 909**
 (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)
 32 cm, Pk/100. **9090-3200** \$83.20
- Fast flow** (245 mL/min), 27 µm retention, creped, **Grade 962**
 (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)
 11 cm, Pk/100. **9620-1100** \$12.35

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

- Slow flow** (4 mL/min), 1.5 µm retention, **Grade 94**
 (Whatman 42; S&S 589/3)
 9 cm, Pk/100. **0940-0900** \$34.85
 15 cm, Pk/100. **0940-1500** \$60.00
- Medium flow** (30 mL/min), 2 µm retention, **Grade 74**
 (Whatman 40; S&S 589/2, 589/4)
 9 cm, Pk/100. **0740-0900** \$34.85
 15 cm, Pk/100. **0740-1500** \$60.00
- Fast flow** (220 mL/min), 15 µm retention, **Grade 55**
 (Whatman 54, 541; S&S 1505)
 9 cm, Pk/100. **0550-0900** \$40.00
 15 cm, Pk/100. **0550-1500** \$63.25

Blotting/Chromatography/Electrophoresis Papers


- 0.19 mm thick**, 125 mm/30 min absorption, **Grade 601**
 (Whatman 1, 1CHR; S&S 593, 597, 2043A)
 46 x 57 cm, Pk/100. **6018-4657** \$157.00
- 0.34 mm thick**, 140 mm/30 min absorption, **Grade 238**
 (Whatman 3MM, S&S 598, 231^o)
 20 x 20 cm, Pk/100. **2388-2020** \$80.00
 46 x 57 cm, Pk/100. **2388-4657** \$305.00
- 0.42 mm thick**, 150 mm/30 min absorption, **Grade 237**
 (Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703; Albet 2040B, FP591; Munktell 3S/H, FN7a)
 46 x 57 cm, Pk/100. **2378-4657** \$352.00
- NEW** **0.95 mm thick**, 180 mm/30 min absorption, **Grade 243**
 (Whatman 17CHR; S&S 2668)
 11 x 14 cm, Pk/100. **2438-1114** \$31.75
 20 x 20 cm, Pk/100. **2438-2020** \$66.15
- 2.48 mm thick**, 240 mm/30 min absorption, **Grade 320**
 (S&S 2727)
 11 x 14 cm, Pk/50. **3208-1114** \$63.45
 20 x 20 cm, Pk/50. **3208-2020** \$93.05



Smooth tips & polypropylene grips

- Forceps, Filter, 11.5 cm. 300-025** \$4.40

Freezer Racks



Freezer Racks, stainless steel, USA

175V-13-2L

Restocking fee of 20% on all returned freezer racks listed below.
 No returns on custom racks.

Horizontal (fits upright freezers)

- Horizontal, Holds 2" Boxes**
- | | 1-9 | ≥10 |
|---|----------|--------|
| 9-Place (3L x 3H). 175HS-9-2 | \$59.00 | 53.00 |
| Dimensions: 16.31(L) x 5.5(W) x 6.5" (H) | | |
| 12-Place (3L x 4H). 175HS-12-2 | \$78.00 | 71.00 |
| Dimensions: 16.31(L) x 5.5(W) x 8.75" (H) | | |
| 15-Place (3L x 5H). 175HS-15-2 | \$98.00 | 89.00 |
| Dimensions: 16.31(L) x 5.5(W) x 11" (H) | | |
| 16-Place (4L x 4H). 175HL-16-2 | \$105.00 | 95.00 |
| Dimensions: 22(L) x 5.5(W) x 8.75" (H) | | |
| 20-Place (4L x 5H). 175HL-20-2 | \$132.00 | 119.00 |
| Dimensions: 22(L) x 5.5(W) x 11" (H) | | |

- Horizontal, Holds 3" Boxes**
- | | | |
|---|---------|-------|
| 9-Place (3L x 3H). 175HS-9-3 | \$61.00 | 55.00 |
| Dimensions: 16.31(L) x 5.5(W) x 9.5"(H) | | |
| 12-Place (4L x 3H). 175HL-12-3 | \$80.00 | 73.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H) | | |

Horizontal Drawer (fits upright freezers)

- Horizontal Drawer, Holds 2" Boxes**
- | | 1-9 | ≥10 |
|--|----------|--------|
| 9-Place (3L x 3H). 175HD-332 | \$85.00 | 77.00 |
| Dimensions: 16.5(L) x 5.5(W) x 7.125"(H) | | |
| 12-Place (3L x 4H). 175HD-342 | \$111.00 | 101.00 |
| Dimensions: 16.5(L) x 5.5(W) x 9.5" (H) | | |
| 15-Place (3L x 5H). 175HD-352 | \$141.00 | 128.00 |
| Dimensions: 16.5(L) x 5.5(W) x 11.812" (H) | | |
| 16-Place (4L x 4H). 175HD-442 | \$151.00 | 137.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5" (H) | | |
| 20-Place (5L x 4H). 175HD-542 | \$188.00 | 171.00 |
| Dimensions: 27.25(L) x 5.5(W) x 9.5" (H) | | |

- Horizontal Drawer, Holds 5" Boxes**
- | | 1-9 | ≥10 |
|---|---------|-------|
| 6-Place (3L x 2H). 175RD-325 | \$74.00 | 67.00 |
| Dimensions: 18.125(L) x 6(W) x 10.75" (H) | | |
| 8-Place (4L x 2H). 175RD-425 | \$99.00 | 90.00 |
| Dimensions: 24(L) x 6(W) x 10.75" (H) | | |

- Horizontal Drawer, Holds 15 mL tubes**
- | | | |
|---|----------|--------|
| 160-Place (80L x 2H). 175HD-32ML-15 | \$99.00 | 90.00 |
| Dimensions: 16.5(L) x 5.5(W) x 9.5"(H) | | |
| 220-Place (110L x 2H). 175HD-42ML-15 | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5" (H) | | |

- Horizontal Drawer, Holds 50 mL tubes**
- | | | |
|---|----------|--------|
| 84-Place (42L x 2H). 175HD-42ML-50 | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H) | | |

Continued