

**Qualitative, Wet-Strengthened Grades** (for vacuum filtration)

- Medium flow** (65 mL/min), 3 µm retention, **Grade 909**  
 (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)  
 32 cm, Pk/100. **9090-3200** ..... \$83.20
- Fast flow** (245 mL/min), 27 µm retention, creped, **Grade 962**  
 (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)  
 11 cm, Pk/100. **9620-1100** ..... \$12.35

**Quantitative, Ashless**

(Designed for use in analytical and gravimetric analysis)

- Slow flow** (4 mL/min), 1.5 µm retention, **Grade 94**  
 (Whatman 42; S&S 589/3)  
 9 cm, Pk/100. **0940-0900** ..... \$34.85  
 15 cm, Pk/100. **0940-1500** ..... \$60.00
- Medium flow** (30 mL/min), 2 µm retention, **Grade 74**  
 (Whatman 40; S&S 589/2, 589/4)  
 9 cm, Pk/100. **0740-0900** ..... \$34.85  
 15 cm, Pk/100. **0740-1500** ..... \$60.00
- Fast flow** (220 mL/min), 15 µm retention, **Grade 55**  
 (Whatman 54, 541; S&S 1505)  
 9 cm, Pk/100. **0550-0900** ..... \$40.00  
 15 cm, Pk/100. **0550-1500** ..... \$63.25

**Blotting/Chromatography/Electrophoresis Papers**

- 0.19 mm thick**, 125 mm/30 min absorption, **Grade 601**  
 (Whatman 1, 1CHR; S&S 593, 597, 2043A)  
 46 x 57 cm, Pk/100. **6018-4657** ..... \$157.00
- 0.34 mm thick**, 140 mm/30 min absorption, **Grade 238**  
 (Whatman 3MM, S&S 598, 231<sup>o</sup>)  
 20 x 20 cm, Pk/100. **2388-2020** ..... \$80.00  
 46 x 57 cm, Pk/100. **2388-4657** ..... \$305.00
- 0.42 mm thick**, 150 mm/30 min absorption, **Grade 237**  
 (Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703; Albet 2040B, FP591; Munktell 3S/H, FN7a)  
 46 x 57 cm, Pk/100. **2378-4657** ..... \$352.00
- NEW** **0.95 mm thick**, 180 mm/30 min absorption, **Grade 243**  
 (Whatman 17CHR; S&S 2668)  
 11 x 14 cm, Pk/100. **2438-1114** ..... \$31.75  
 20 x 20 cm, Pk/100. **2438-2020** ..... \$66.15
- 2.48 mm thick**, 240 mm/30 min absorption, **Grade 320**  
 (S&S 2727)  
 11 x 14 cm, Pk/50. **3208-1114** ..... \$63.45  
 20 x 20 cm, Pk/50. **3208-2020** ..... \$93.05



Smooth tips & polypropylene grips

- Forceps, Filter, 11.5 cm. 300-025** ..... \$4.40

**Freezer Racks**



175V-13-2L

**Freezer Racks, stainless steel, USA**

Restocking fee of 20% on all returned freezer racks listed below.  
 No returns on custom racks.

**Horizontal** (fits upright freezers)

- Horizontal, Holds 2" Boxes**
- |   | 1-9      | ≥10    |
|---|----------|--------|
| 9-Place (3L x 3H). <b>175HS-9-2</b> .....   | \$59.00  | 53.00  |
| Dimensions: 16.31(L) x 5.5(W) x 6.5" (H)    |          |        |
| 12-Place (3L x 4H). <b>175HS-12-2</b> ..... | \$78.00  | 71.00  |
| Dimensions: 16.31(L) x 5.5(W) x 8.75" (H)   |          |        |
| 15-Place (3L x 5H). <b>175HS-15-2</b> ..... | \$98.00  | 89.00  |
| Dimensions: 16.31(L) x 5.5(W) x 11" (H)     |          |        |
| 16-Place (4L x 4H). <b>175HL-16-2</b> ..... | \$105.00 | 95.00  |
| Dimensions: 22(L) x 5.5(W) x 8.75" (H)      |          |        |
| 20-Place (4L x 5H). <b>175HL-20-2</b> ..... | \$132.00 | 119.00 |
| Dimensions: 22(L) x 5.5(W) x 11" (H)        |          |        |

- Horizontal, Holds 3" Boxes**
- |   |         |       |
|---|---------|-------|
| 9-Place (3L x 3H). <b>175HS-9-3</b> .....   | \$61.00 | 55.00 |
| Dimensions: 16.31(L) x 5.5(W) x 9.5"(H)     |         |       |
| 12-Place (4L x 3H). <b>175HL-12-3</b> ..... | \$80.00 | 73.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H)        |         |       |

**Horizontal Drawer** (fits upright freezers)

- Horizontal Drawer, Holds 2" Boxes**
- |  | 1-9      | ≥10    |
|--|----------|--------|
| 9-Place (3L x 3H). <b>175HD-332</b> .....  | \$85.00  | 77.00  |
| Dimensions: 16.5(L) x 5.5(W) x 7.125"(H)   |          |        |
| 12-Place (3L x 4H). <b>175HD-342</b> ..... | \$111.00 | 101.00 |
| Dimensions: 16.5(L) x 5.5(W) x 9.5" (H)    |          |        |
| 15-Place (3L x 5H). <b>175HD-352</b> ..... | \$141.00 | 128.00 |
| Dimensions: 16.5(L) x 5.5(W) x 11.812" (H) |          |        |
| 16-Place (4L x 4H). <b>175HD-442</b> ..... | \$151.00 | 137.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5" (H)      |          |        |
| 20-Place (5L x 4H). <b>175HD-542</b> ..... | \$188.00 | 171.00 |
| Dimensions: 27.25(L) x 5.5(W) x 9.5" (H)   |          |        |

- Horizontal Drawer, Holds 5" Boxes**
- |   | 1-9     | ≥10   |
|---|---------|-------|
| 6-Place (3L x 2H). <b>175RD-325</b> ..... | \$74.00 | 67.00 |
| Dimensions: 18.125(L) x 6(W) x 10.75" (H) |         |       |
| 8-Place (4L x 2H). <b>175RD-425</b> ..... | \$99.00 | 90.00 |
| Dimensions: 24(L) x 6(W) x 10.75" (H)     |         |       |

- Horizontal Drawer, Holds 15 mL tubes**
- |   |          |        |
|---|----------|--------|
| 160-Place (80L x 2H). <b>175HD-32ML-15</b> .....  | \$99.00  | 90.00  |
| Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)            |          |        |
| 220-Place (110L x 2H). <b>175HD-42ML-15</b> ..... | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5" (H)             |          |        |

- Horizontal Drawer, Holds 50 mL tubes**
- |   |          |        |
|---|----------|--------|
| 84-Place (42L x 2H). <b>175HD-42ML-50</b> ..... | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H)            |          |        |

Continued

**Horizontal Drawer-Tray** (fits upright freezers)


<b>Drawer-Tray, Holds 2" Boxes</b> , aluminum trays	1-9	≥10
12-Place (3L x 4H). <b>175HDT-342</b> .....	\$96.00	86.00
Dimensions: 16.63(L) x 5.5(W) x 8.5" (H)		
18-Place (3L x 6H). <b>175HDT-362</b> .....	\$143.00	129.00
Dimensions: 16.63(L) x 5.5(W) x 12.75" (H)		
24-Place (4L x 6H). <b>175HDT-462</b> .....	\$189.00	172.00
Dimensions: 22(L) x 5.5(W) x 12.75" (H)		

**Vertical** (fits chest freezers)

<b>Vertical, Holds 2" Boxes</b>		
13-Place. <b>175V-13-2L</b> .....	\$88.00	79.00
Dimensions: 5.5(L) x 5.625(W) x 28.562" (H)		

**Tube Storage Boxes** (cryogenic) – see in GENERAL LAB SUPPLIES & EQUIPMENT section

**Funnels**



**Funnels**

Disposable, polypropylene, short stems, Simport		
★ 55 mm Top ID, takes 11 cm filter, Cs/100. <b>F490-3</b> .....	\$40.70	
★ 65 mm Top ID, takes 12.5 cm filter, Cs/100. <b>F490-4</b> .....	\$40.70	
Glass (borosilicate), long stem (150 mm)		
50 mm Top ID, Pk/6. <b>600GF614050</b> .....	\$49.95	
65 mm Top ID, Pk/6. <b>600GF614065</b> .....	\$60.00	
75 mm Top ID, Pk/6. <b>600GF614075</b> .....	\$62.00	
Heavy-Duty, Molded-In Handle, polypropylene		
250 mm Top ID, 2.5 L capacity, each. <b>60057232</b> .....	\$26.55	
Micro, polypropylene, Bel-Art		
24 mm Top ID, 2.6 mL capacity, Pk/12. <b>1146850024</b> .....	\$42.30	
28 mm Top ID, 5.2 mL capacity, Pk/12. <b>146850028</b> .....	\$34.65	
43 mm Top ID, 15.7 mL capacity, Pk/12. <b>146850043</b> .....	\$72.95	
Powder, polypropylene		
Parallel Stem		
78 mm Top ID, Pk/12. <b>60057222</b> .....	\$16.85	
100 mm Top ID, Pk/6. <b>146600100</b> .....	\$16.35	
145 mm Top ID, Pk/6. <b>60057224</b> .....	\$21.50	

**Funnels** (cont.)

Standard Stem, Premium, polypropylene, Bel-Art	
Exclusive filter-step internal design provides speedier filtration; external ribs prevent air-lock.	
55 mm Top ID, takes 9 cm filter, Pk/12. <b>146760055</b> .....	\$47.75
65 mm Top ID, takes 11 cm filter, Pk/12. <b>146760065</b> .....	\$49.60
75 mm Top ID, takes 12.5 cm filter, Pk/6. <b>146760075</b> .....	\$25.65
100 mm Top ID, takes 18.5 cm filter, Pk/4. <b>146760100</b> .....	\$32.60
Standard Stem, Economy, polypropylene	
50 mm Top ID, each. <b>310FUNPP50</b> .....	\$0.80
62 mm Top ID, each. <b>310FUNPP62</b> .....	\$0.95
75 mm Top ID, each. <b>310FUNPP75</b> .....	\$1.15
100 mm Top ID, each. <b>310FUNPP100</b> .....	\$1.65

**Gauze Pads**

These 4 x 4" nonsterile gauze sponges are 12-ply.



<b>Gauze Pads</b> , 10 pk/200, Cs/2000. <b>44122000</b> .....	\$88.45
---	---------

**Germicidal Wipes**

Sani-Cloth Plus



<b>Germicidal Wipes</b> , Tub/160 <b>440Q89072</b> .....	\$10.00
--	---------

**Gloves & Hand Protection**



Nitrile Gloves

**Gloves**

Purple Nitrile, powder-free, 5.5 mil, textured fingers, Halyard	
Nitrile gloves are hypoallergenic and provide greater puncture and chemical resistance than latex or vinyl gloves.	
Small, bx/100, Cs/10 Bx. <b>55081</b> .....	\$139.00
Medium, bx/100, Cs/10 Bx. <b>55082</b> .....	\$139.00
Large, bx/100, Cs/10 Bx. <b>55083</b> .....	\$139.00
X-Large, bx/90, Cs/10 Bx. <b>55084</b> .....	\$139.00
Purple Nitrile-XTRA, 12" Long, powder-free, 5.9 mil, Halyard	
Offers additional protection against chemotherapy drugs and common chemicals.	
Small, bx/50, Cs/10 Bx. <b>50601</b> .....	\$213.00
Medium, bx/50, Cs/10 Bx. <b>50602</b> .....	\$213.00
Large, bx/50, Cs/10 Bx. <b>50603</b> .....	\$213.00
X-Large, bx/50, Cs/10 Bx. <b>50604</b> .....	\$213.00

Continued