


Accessories & Parts

Counting Chambers

Bright-Line®
Hemocytometer



- Hemocytometer**, Bright-Line®, molded borosilicate glass, **USA**
Reduced Thickness/phase counting chamber (3.5 mm thick)
Chamber with 2 Coverslips, Hausser. **1475** \$305.00
- Howard Mold Chamber**, Hausser, **USA**
Requires Howard Disc - see RETICLES & STAGE MICROMETERS
Chamber with 2 Coverslips (5080 & 5090). **3820** \$383.00
Coverslip only, 28 x 33 x 0.5 mm (thin), each. **5080** \$46.80

England Field Finder, 1 mm thick, **USA. KR-834** \$253.00

Eyepiece Pointer, wire. **14301** \$2.00



Eyetebe Sleeve, 23.2 mm increased to 30 mm. **14320** \$35.00

Lens Tissue

- Fine, 10 x 15 cm sheets, Bk/50 Sheets. **LENSF-46-50** \$3.85
Type 1, Class 1, lintless
- Premium, 7.5 x 15 cm sheets, Pk/50 Sheets. **77605** \$4.00



- Microscope Cleaning Kit. 77600** \$30.00
- Replacement Items
- Blower Bulb. **77602** \$10.50
 - Brush, 1½", polyester. **77603** \$3.00
 - Cleaning Stick. **77604** \$2.00
 - Lens Tissue, Premium, 7.5 x 15 cm, 50 sheets. **77605** \$4.00
 - Solvent Jar, (Empty). **77607** \$2.50
 - Utility Box. **77608** \$9.00



Protect your microscope from dust, light, and moisture with a durable rip-stop nylon, washable cover.

Microscope Dust Covers, rip-stop nylon, burgundy, **USA**

- AO100** \$23.00
- 4BIN** \$10.00
Fits Binocular Head Only - 6½" x 8"(W)
- NBIN** \$10.00
Fits Binocular Head Only - 11" x 8"(W)
- NSML** \$25.00
Fits Nikon Alphaphot (1&2), SC; Stereos; National 420,162,168; Swift M15B-P
- NDFC** \$44.00
Fits Nikon Diaphot w/EPI-FL & Camera; Jenco BC-2A1, BC-364
- D200** \$44.00
Fits Nikon Diaphot 200/300; Olympus Vanox
- D2FC** \$44.00
Fits Nikon Diaphot 200/300 w/EPI-FL & Camera, Ti
- N400** \$30.00
Fits Nikon E400/E600/50i/55i, TS100; Jenco CP-2A1; Meiji IM7100; Olympus BX 40/41, CK2 w/mech. stage
- N4CA** \$30.00
Fits Nikon E400/E600 w/Camera; Meiji IM7200; Olympus BX 40/41 w/Camera, Zeiss Axioskop 2 w/Camera
- N4FL** \$30.00
Fits Nikon E400/E600 with EPI-FL
- NLAB** \$23.00
Fits Nikon Labophot (1&2), Optiphot (1&2), TMS/TMS-F, E200; Olympus CK2; Leica DME, DM500/DM750
- NLCA** \$30.00
Fits Nikon Labophot (1&2)/Optiphot (1&2) w/Camera; Meiji MT4300/5300
- NLFC** \$30.00
Fits Nikon Labophot (1&)/Optiphot (1&2) w/EPI-FL & Camera
- NLTH- *** \$33.00
Fits Nikon Labophot (1&2)/Optiphot (1&2) w/TH-L
- NLTC- *** \$38.00
Fits Nikon Labophot (1&2)/Optiphot (1&2) w/TH-L & Camera
**Both TH-L covers include a cover for the student head.*
Specify opening as Right, Left, or Both sides using suffix.
- NLMV** \$28.00
Fits Nikon Labophot (1&2)/Optiphot (1&2) w/Multi-View
- NLMC** \$30.00
Fits Nikon Labophot (1&2)/Optiphot (1&2) w/M-V & Camera

SLIDE & CASSETTE PROCESSING

Microscope Slides

Slides, general, 75 x 25 x 1 mm

Plain

Globe, ground edges, 90° corners
Cs/20 Bx, bx/72. **1301** \$152.00

Frosted (etched)

Globe, beveled edges, clipped corners
Cs/20 Bx, bx/72. **1334** \$204.00

White End

Globe, beveled edges, clipped corners
Cs/20 Bx, bx/72. **1334W** \$238.00

Superfrost®, "SPECIMEN". **4951-001T** \$327.00
(12-550-143, 48311-600)

Superfrost® Plus



Concavity, 2-Well



Slides, specialty

Charged (silanated), 75 x 25 x 1 mm

Recommended for adhering tissue sections
and cell preparations.

Diamond™, white end, ground edges, 90° corners
Cs/20 Bx, bx/72. **1358W** \$475.00

Superfrost® Plus, white end, Erie (EpreDia)
Cs/20 Bx, bx/72. **4951PLUS-001** \$682.00 ↓
Preferred slide for Ventana instrumentation.
(12-550-15, 48311-703)

Concavity

75 x 25 x 1.3 mm, 16 mm dia., 0.8 mm deep, economy
1-Well, Bx/12. **600CS3X11** \$7.75

2-Well, Bx/12. **600CS3X12** \$10.75

76 x 26 x 3.2 mm, 18 mm dia., 1.75 mm deep, Erie (EpreDia)
1-Well, Bx/12. **1526** \$97.35

Fused Quartz, NSG Type OZ®, 3" x 1" x 1 mm, USA

QMS3X1X1, each \$20.00

Fused Silica, Corning 7980, 3" x 1" x 1 mm, USA

FSMS3X1X1, each \$39.90

Coverslips (Cover Glasses)

Coverslips, premium, borosilicate glass, 1 oz/pk, Cs/10 oz

#0 Thickness (0.08-0.12 mm)

Circles, Germany*

25 mm, ~222 pcs/oz, Cs/10 oz. **0-25** \$214.00



#1 Thickness (0.13-0.17 mm)

Squares & Rectangles, Erie (EpreDia)


22 x 22 mm, ~162 pcs/oz, Cs/10 oz. **1-2222** \$70.70

22 x 30 mm, ~119 pcs/oz, Cs/10 oz. **1-2230** \$67.00

22 x 50 mm, ~71 pcs/oz, Cs/10 oz. **1-2250** \$67.00

24 x 50 mm, ~65 pcs/oz, Cs/10 oz. **1-2450** \$67.00

24 x 60 mm, ~54 pcs/oz, Cs/10 oz. **1-2460** \$67.00



#1 Thickness (0.13-0.17 mm)

Circles, Germany*

12 mm, ~663 pcs/oz, Cs/5 oz. **1-12** \$249.00

18 mm, ~296 pcs/oz, Cs/5 oz. **1-18** \$150.00



#1.5 Thickness (0.16-0.19 mm). Recommended for most users

Squares & Rectangles, Erie (EpreDia)

18 x 18 mm, ~209 pcs/oz, Cs/10 oz. **1.5-1818** \$71.70

22 x 22 mm, ~140 pcs/oz, Cs/10 oz. **1.5-2222** \$73.35

22 x 50 mm, ~62 pcs/oz, Cs/10 oz. **1.5-2250** \$67.00

24 x 40 mm, ~71 pcs/oz, Cs/10 oz. **1.5-2440** \$67.00

24 x 50 mm, ~56 pcs/oz, Cs/10 oz. **1.5-2450** \$63.80

24 x 60 mm, ~47 pcs/oz, Cs/10 oz. **1.5-2460** \$67.00



#1.5 Thickness (0.16-0.19 mm). Recommended for most users

Circles, Germany*

12 mm, ~552 pcs/oz. **1.5-12** cs/5 oz \$198.00 cs/10 oz —

18 mm, ~245 pcs/oz. **1.5-18** \$250.00 —

22 mm, ~164 pcs/oz. **1.5-22** — \$252.00

25 mm, ~127 pcs/oz. **1.5-25** — \$252.00



*German glass is preferred for cell growth/adhesion

Coverslips, specialty

Fused Quartz, USA

0.15-0.25 mm Thick, NSG Type OZ®

1" Square, each. **QCS1X1** \$19.00

1" Circle, each. **QCS1X0** \$19.00

0.16 ± 0.01 mm Thick, A1

25 mm Square, each. **R525025** \$32.40

25 mm Circle, each. **R525000** \$23.40

Fused Silica, USA

0.15-0.25 mm Thick, Corning 7980

1" Square, each. **FSCS1X1** \$21.50

1" Circle, each. **FSCS1X0** \$21.50

0.16 ± 0.01 mm Thick, S1-UV

25 mm Square, each. **R425025** \$32.00

25 mm Circle, each. **R425000** \$28.80