


LAB, SURGICAL & DENTAL INSTRUMENTS

Legend

- Economy:** stainless steel; forged and finished in Pakistan
Standard: surgical-grade (OR) German steel; forged in Germany and finished in Pakistan and/or Germany
German: surgical-grade (OR) German steel; made entirely (forged and finished) in Germany
Roboz: Roboz Surgical. (RS-)
-TC: Tungsten Carbide inlays; identified by gold handles

Forceps (Tweezers) & Tongs



General

Forceps

Fine Point, serrated handles, stainless steel

Straight, serrated tips & handles

- 10 cm (4"), Economy, guide pin. **26-275** \$10.75
- 10 cm (4"), Standard, guide pin. **40-720** \$16.50
- 10 cm (4"), German, guide pin. **420-802** \$32.80
- 10 cm (4"), Roboz, 0.8 mm tip, guide pin. **RS-5130** \$69.30
- ★ 11.5 cm (4½"), Economy, guide pin. **300-130** \$2.80
- 11.5 cm (4½"), Economy, No Pin. **20-373** \$3.35
- 20.5 cm (8"), Economy, guide pin. **300-141** \$8.10
- 25.4 cm (10"), Economy, guide pin. **300-143** \$9.20

Straight, Smooth Tips, serrated handles

- 10 cm (4"), Roboz, 0.8 mm tip, guide pin. **RS-5140** \$69.30
- 20.5 cm (8"), Economy, guide pin. **300-111** \$9.20

Angled, serrated tips

- 10 cm (4"), Roboz, 0.8 mm tip, guide pin. **RS-5139** \$73.25
- 15.5 cm (6"), Economy, guide pin. **3600** \$4.00

Half Curve (45°), serrated tips & handles

- 10 cm (4"), Economy, guide pin. **26-280** \$10.75
- 10 cm (4"), Standard, guide pin. **40-722** \$16.50
- 10 cm (4"), Roboz, 0.8 mm tip, guide pin. **RS-5135** \$73.25
- ★ 11.5 cm (4½"), Student, No Pin. **4190** \$3.30
- ★ 11.5 cm (4½"), Economy, guide pin. **300-131** \$2.90

Full Curve (90°), serrated tips & handles

- 10 cm (4"), Economy, guide pin. **24-665S** \$6.10
- 10 cm (4"), Roboz, 0.8 mm tip, guide pin. **RS-5137** \$73.25

Medium-Fine Point, serrated handles, stainless steel

- Curved, smooth tips, no pin, 15.2 cm (6"). **300-011** \$7.50

Medium Point (Thumb Dressing), serrated handles, stainless steel

Straight, serrated tips, no guide pin

- 9 cm (3½"), Economy. **17-315** \$2.25
- 11.5 cm (4½"), Economy. **300-038** \$2.85
- 11.5 cm (4½"), Standard. **40-700** \$15.80
- 11.5 cm (4½"), Roboz, 1.5 mm tip. **RS-8100** \$33.00
- 15.5 cm (6"), Economy. **300-041** \$3.40
- 15.5 cm (6"), Standard. **40-710** \$12.60
- 15.5 cm (6"), Roboz. **RS-8106** \$34.00
- 16 cm (6¼"), Economy, black matte. **TWZ-703.74** \$4.50
- 20.5 cm (8"), Economy. **2240** \$4.05
- 20.5 cm (8"), Standard. **540-055** \$21.00
- 25.4 cm (10"), Economy. **2250** \$4.80
- 25.4 cm (10"), Standard. **540-060** \$31.50
- 30.5 cm (12"), Economy. **2260** \$5.85
- 51 cm (20"), Economy. **19-380** \$17.90

Exclusive! Extra wide tips ideal for removing specimens from storage containers.

Straight, Smooth Tips, guide pin

- 15.5 cm (6"), economy. **17-320SM** \$3.60

Curved, serrated tips, guide pin

- 11.5 cm (4½"), economy. **17-655** \$3.30



Smooth handles do not have guide pin.

Gooseneck, blunt serrated tips, 15 cm (6"). **300-001** \$5.30



Student Grade, stainless steel, 11.5 cm (4½"). **18-915** \$3.00


straight, serrated tips & handles; does not have guide pin

Corrosion & heat resistant ABS plastic is nonconductive & noncontaminating. Will not scratch optics.



Plastic, straight, serrated tips & handles, 13 cm. **97-3901** \$0.45

Micro



#5 Style w/Hook

Forceps, Micro, smooth tips* & handles, no guide pin

Economy, stainless steel

- #2A Style, straight, 12 cm. **57.0752A** \$7.45
- #3 Style, straight, 12 cm. **19-361** \$5.20
- ★ #5 Style, straight, 11.5 cm. **19-364** \$5.20
- #5 Style w/Hook, straight, 11 cm. **19-364HK** \$22.95
- #7 Style, curved, 11.5 cm. **19-366** \$5.20
- #7B Style, curved, *serrated, nonmag, 11.5 cm. **57.838**... \$14.35



Style #3 (very sharp)



Style #5 (needle sharp)



Style #7 (sharp)

Forceps, Micro (cont)

Dumont Economy, Inox, magnetic, Swiss

- #3 Style, straight, 12 cm (4 3/4"). **TWZ-301.16** \$33.00
- #5 Style, straight, 11 cm (4 1/4"). **TWZ-301.22** \$40.40
- #7 Style, curved, 11.5 cm (4 1/2"). **TWZ-301.28** \$37.55
- SS Style, nonmag, straight, 13.5 cm (5 1/4"). **TWZ-302.56** \$31.00

Dumont, Inox, high-precision, magnetic, Swiss

- #3 Style, straight, 12 cm (4 3/4")
 - Standard Tips (0.17W x 0.10 mm). **RS-5042** \$44.55
- #5 Style, straight, 11 cm (4 1/4")
 - Standard Tips (0.10W x 0.06 mm). **RS-5045** \$38.60
 - Biologie Tips (0.05W x 0.01 mm). **RS-5015** \$44.55
- #55 Style, light shanks, straight, 11 cm (4 1/4")
 - Biologie Tips (0.05W x 0.01 mm). **RS-5010** \$51.50
- #7 Style, curved, 11.5 cm (4 1/2")
 - Standard Tips (0.17W x 0.10 mm). **RS-5047** \$48.50



Dumont, Dumoxel, standard tips (0.10W x 0.06 mm), Swiss
#5/45 Style, tips angled 45°, 11 cm (4 1/4"). **RS-5058** \$46.55

Ideal-tek, stainless steel, nonmagnetic, Swiss

- SS Style, 14 cm (5 1/2"). **TWZ-312.56** \$17.95

Specialty



Forceps, stainless steel

- Adson Dressing, serrated tips & handles, 12 cm. **17-555** \$8.55



- Balance, curved plastic tips, 10 cm (4"). **300-006** \$5.25



- Bishop Harmon, fine serrated tips, 9 cm (3 1/2"). **4211** \$18.00



Cutting, economy

- Style 15A, 11.5 cm (4 1/2"). **TWZ-801.60** \$9.95



- Filter, smooth tips & polypropylene grips, 11.5 cm. **300-025** \$4.80



- Gerald Dressing, 17.5 cm (7"). **19-367** \$6.35



Graefe Eye Dressing, serrated tips & handles, 7 cm (2 1/2")

- Straight, standard. **4215** \$18.00
- Curved, standard. **4216** \$18.00



- Insect Pinning, rounded & serrated tips, 11 cm. **17-682** \$4.15



Paddle Tip, grid-serrated pads, economy

- 22 x 22 mm Pad, standard, 16 cm (6 1/4"). **33-1045** \$14.50
- 30 x 35 mm(W) Pad, 15.2 cm (6"). **33-1045B** \$21.00
- 30 x 40 mm(W) Pad, 15.2 cm (6"). **33-1045A** \$21.00



Wide Handle (Paddle Handle), economy

serrated and perforated grips

Straight, Serrated tips

- 14 cm (5 1/2"), 1.5 mm tips. **15-600** \$25.00

Straight, Smooth tips

- 14 cm (5 1/2"), 0.6 mm tips. **15-625** \$25.00
- 14 cm (5 1/2"), 1.5 mm tips. **15-585** \$25.00
- 17.5 cm (7"), 1.5 mm tips. **15-550** \$25.00
- 17.5 cm (7"), 3 mm tips. **15-500** \$25.00

Angled, Serrated tips

- 14 cm (5 1/2"), 1.5 mm tips. **15-595** \$25.00

Angled, Smooth tips

- 14 cm (5 1/2"), 1.5 mm tips. **15-590** \$25.00
- 17.5 cm (7"), 1.5 mm tips. **15-580** \$25.00

Spade Tips



Angled (Kuehne), 11 cm (4 1/4")

- Stainless Steel, economy. **33-1041** \$4.15
- Teflon® Coated, Bel-Art. **379420000** \$35.60
- Straight, stainless steel, 10.5 cm, economy. **57.926** \$8.45



Square Point Tips, straight

- Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379430000** \$39.40

Locking



Forceps, stainless steel

Self-Closing (cross-locking), squeeze handles to open tips

Fine Point, smooth tips

Economy, 13 cm (5 $\frac{1}{8}$ "). **300-066** \$5.45

Peer-Vigor, 11 cm (4 $\frac{1}{4}$ "), Swiss. **57.324** \$42.40

Medium Point, serrated tips, 17 cm (6 $\frac{5}{8}$ "). **TWZ-948.00** \$3.95



Slide-Locking, grid-serrated tips

Fine Point

Economy, 14 cm (5 $\frac{1}{2}$ "). **57.0357A** \$9.00

Dumont, 16 cm (6 $\frac{1}{4}$ "). **TWZ-101.72** \$45.75

Medium Point

Economy, 16 cm (6 $\frac{1}{4}$ "). **203** \$9.00

Dumont, 16 cm (6 $\frac{1}{4}$ "). **TWZ-101.76** \$43.60

Snap-Closing, angled, fine points, 15.5 cm (6"). **3610** \$5.80

serrated tips & handles, guide pin

Tissue



Forceps, stainless steel

Adson-Brown, 9 x 9 side grasping teeth, 12 cm. **17-574** \$16.05


Adson Tissue, 1 x 2 teeth, 12 cm (4 $\frac{3}{4}$ "). **17-560** \$8.55

Iris Tissue, 1 x 2 teeth

Straight, economy, 10 cm (4"). **26-295** \$10.75

Angled, economy, 11.5 cm (4 $\frac{1}{2}$ "). **19-350ANG** \$5.00

Half Curve, standard, 10 cm (4"). **45-822** \$17.60



Russian Tissue, circular serrated tips, 15.5 cm (6"). **2275** \$5.25

Tissue, 1 x 2 teeth, 14 cm (5 $\frac{1}{2}$ "). **17-145** \$8.55

Tongs





Tongs, Crucible, serrated jaws, stainless steel

24.5 cm (9 $\frac{3}{4}$ "). **320-101** \$12.50

51 cm (20"). **320-120** \$30.15

Hemostats & Clamps



Hemostats, serrated jaws

Hartman Mosquito, stainless steel

Straight

Economy, 9 cm (3 $\frac{1}{2}$ "). **12-180S** \$6.05

Standard, 9 cm (3 $\frac{1}{2}$ "). **20-352** \$21.30

Curved

Economy, 9 cm (3 $\frac{1}{2}$ "). **306-012** \$9.10

Standard, 10 cm (4"). **20-360** \$18.75↓

Halstead Mosquito, 12.5 cm (5"), stainless steel

Straight

Economy. **14-065** \$10.00

Standard. **20-501** \$22.45

Curved

Economy. **14-070** \$10.00

Standard. **20-519** \$22.85

Kelly, stainless steel

Straight

Economy, 14.5 cm (5 $\frac{3}{4}$ "). **306-016** \$8.20

Standard, 14.5 cm (5 $\frac{3}{4}$ "). **20-550** \$23.60

Curved

Economy, 14 cm (5 $\frac{1}{2}$ "). **306-017** \$9.90

Standard, 14.5 cm (5 $\frac{3}{4}$ "). **20-568** \$23.95

Rochester-Pean, 20.5 cm (8"), stainless steel

Straight

Economy. **1920** \$6.50

Standard. **21-977** \$33.95

Curved


Economy. **1930** \$6.50

Standard. **21-979** \$33.95

Polypropylene, 12.5 cm (5"). **HEM-300.00** \$3.25



Specimen Forceps, 25 cm (9 $\frac{3}{4}$ "), serrated jaws. **300-500** \$18.70



Sponge Forceps, Foerster, serrated jaws, stainless steel

Straight

17.5 cm (7"), economy. **14-660** \$17.50

24 cm (9 $\frac{1}{2}$ "), economy. **14-670** \$21.10

Curved, 17.5 cm (7"), economy. **14-665** \$17.50

Clamps



300-120



300-121

Dieffenbach Serrefine, stainless steel, economy

Straight, serrated jaws and handles

60 mm (2³/₈"). **300-120** \$7.40

Curved, serrated jaws and handles

48 mm (1⁷/₈"). **12-248** \$9.15

60 mm (2³/₈"). **300-121** \$7.40



Riveted handle
opens for cleaning.

Johns Hopkins Bulldog, Roboz

Straight, 38 mm (1¹/₂"). **RS-7440** \$88.10



RS-5452

Schwartz Vessel Clips, 1.7 x 8 mm jaws, 2.5 cm long, Germany

Straight Jaws, Roboz. **RS-5452** \$36.00


Slightly Angled Jaws, Roboz. **RS-5457** \$36.00

8.6 cm shaft; 5 mm serrated jaw



Micro Alligator Forceps, stainless steel. **2769** \$50.00

Needle Holders & Wound Closure




Crile-Wood



Olsen-Hegar



Webster



Needle Holders, stainless steel

Baumgartner, cross-serrated jaws, 14 cm (5¹/₂"). **16-285S** \$11.05

Castroviejo, cross-serrated jaws w/lock, 14 cm. **15-302** \$50.00

Crile-Wood, cross-serrated jaws, 15 cm (6")

Economy, central groove. **16-293S** \$13.45

Standard. **32-485** \$33.75

German. **320-366** \$60.80

Derf, 12 cm (4³/₄"), economy. **16-238** \$9.70

serrated jaws with central groove

Olsen-Hegar, cross-serrated jaws with scissors, 14 cm (5¹/₂")

Economy, central groove. **16-296** \$14.30

Economy, No Groove. **16-296A** \$14.45

Standard, central groove. **32-500** \$35.25

Webster, smooth jaws, 12 cm (4³/₄")


Economy. **15-035** \$17.45

Standard. **32-101** \$30.75


Tungsten Carbide. **0817** \$45.00

German. **320-101** \$60.50

Michel Wound Closure System



Michel Clip Applier/Remover, double-ended. **3700** \$13.75
11.5 cm (4¹/₂"), stainless steel



Michel Wound Clips, stainless steel

11 mm, Bx/1000. **380-519** \$74.00

14 mm, Bx/1000. **380-543** \$89.00

16 mm, Bx/1000. **380-568** \$97.00

18 mm, Bx/1000. **380-584** \$98.00

22 mm, Bx/1000. **380-626** \$116.00

Reflex Wound Closure System

CellPoint Scientific, USA

Reflex Clip Applier

For 7 mm Clips. **204-1000** \$165.00

★ For 9 mm Clips. **202-1000** \$165.00


Reflex Wound Clips, nonsterile, stainless steel

7 mm, 5 packs of 20 clips, Bx/100. **203-1000** \$56.00

★ 9 mm, 5 packs of 20 clips, Bx/100. **201-1000** \$46.00

Reflex Clip Remover, fits 7 mm & 9 mm clips. **205-1000** \$99.00

Inoculating Loops/Needle & Cell Spreaders



Inoculating Loops

Disposable

1 µL, Sterile Pk/20, PP, Cs/1000. **HS81121A** \$69.75

1 µL, Sterile Pk/1, PS, Cs/250. **L200-11** \$106.00

Nichrome, 6 mm twisted loop of 26-gauge wire. **600INL703** \$2.00

21 cm aluminum handle

Inoculating Needle. 600INN803 \$2.50

75 mm, 26-gauge nichrome wire; 20 cm aluminum handle

16 mm bar width;
20 cm handle



Cell Spreader, Stainless Steel. 377360006 \$26.10

Microtome Blades

Microtome Blades - see Microtome Blades in Tissue (Embedding) & Biopsy Cassettes section


Accessories

Slide & Coverslip Handling

Forceps


Spade Tips

- Angled (Kuehne), 11 cm (4 1/4")
 Stainless Steel, economy. **33-1041** \$4.15
 Teflon® Coated, Bel-Art. **379420000** \$35.60
 Straight, stainless steel, 10.5 cm, economy. **57.926** \$8.45



Square Point Tips, straight

- Stainless Steel, 12 cm (4 3/4"), economy. **1244** \$2.25
 Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379430000** \$39.40



Minimize contamination during slide & coverslip handling. Will not scratch; corrosion and heat resistant.

Hemostat, Polypropylene, 12.5 cm. **HEM-300.00** \$1.95

Immersion Oil & Bottles



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities.
USA

Immersion Oil

Type A, low viscosity (150 cSt). Short W.D. optics

- 7.5 mL in LDPE dispensing bottle. **18101** \$10.00
 20 mL in HDPE dispensing bottle. **18102** \$15.35
 480 mL in amber bottle. **18104** \$110.00

Type B, high viscosity (1250 cSt). Long W.D. optics

- 7.5 mL in LDPE dispensing bottle. **18105** \$10.00
 20 mL in HDPE dispensing bottle. **18106** \$15.35
 480 mL in amber bottle. **18108** \$110.00

Type NVH, very high viscosity (21,000 cSt)

- 480 mL in amber bottle. **18113** \$67.00


Type HF, low fluorescence (700 cSt)

- 480 mL in amber bottle. **18118** \$169.00

Type FF, fluorescence free (170 cSt)

- 480 mL in amber bottle. **18124** \$163.00

A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



Bottle, Dispensing, HDPE, 20 mL, each. **18145** \$8.50

Slide Storage

PLASTIC SLIDE BOXES




25-Slide
Squeeze-Release Cover



100-Slide
Steel-Hinged Cover

50-Slide
Hinged Cover



Slide Boxes, hold 75 x 25 mm (3 x 1") slides

- 25-Slide, squeeze-release cover, cork bottom, 9.5 x 8.5 x 3 cm
 White, ABS plastic. **97-0025** \$2.50
 50-Slide, hinged cover, cork bottom, 20 x 9.5 x 3.5 cm
 Blue, ABS plastic. **97-0056** \$6.80
 Red, ABS plastic. **97-0058** \$6.80
 100-Slide, steel-hinged cover, cork bottom
 Blue, ABS plastic. **97-0126** \$8.90
 Red, ABS plastic. **97-0128** \$8.90
 White, ABS plastic. **97-0125** \$8.90



Store 4500 75 x 25 mm slides in 6 divided drawers. Each cabinet comes with a sample of 6 compressor follow blocks and 6 drawer labels. Safely stack up to 12 cabinet units on a single base.
 18 3/4"(D) x 15 1/4"(W) x 5"(H)

Slide Cabinets, steel

- 6-Drawer, stackable
 Green. **513500G** \$255.00
 Tan. **513500T** \$255.00
 White. **513500W** \$255.00
 Memorandum Markers, paperboard, Pk/100. **C108** \$12.75
 Compressor Follow Block, rubber, each. **C109** \$1.90
 Base, raises cabinets 4 3/4" off floor
 Green. **513503G** \$169.00
 Tan. **513503T** \$169.00
 White. **513503W** \$169.00



Cardboard
Folder



Plastic
Folder

Slide Folders, 20-slide, hold 75 x 25 mm (3 x 1") slides

- Cardboard, black tape, economy. **HS9902** \$6.35
 (12-587-10)
 Plastic, gray tray w/2 clear doors. **M750-20GY** \$11.60
 Transparent covers hold visible slides securely in place, but can fold flat against back of tray when opened. Slides are slightly elevated for easy retrieval.

Qualitative, Wet-Strengthened Grades (for vacuum filtration)

Medium flow (65 mL/min), 3 µm retention, **Grade 909**

(Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)

32 cm, Pk/100. **9090-3200**\$83.20

Fast flow (245 mL/min), 27 µm retention, creped, **Grade 962**

(Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)

11 cm, Pk/100. **9620-1100**\$12.35

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

Slow flow (4 mL/min), 1.5 µm retention, **Grade 94**

(Whatman 42; S&S 589/3)

9 cm, Pk/100. **0940-0900**\$34.85

15 cm, Pk/100. **0940-1500**\$60.00

Medium flow (30 mL/min), 2 µm retention, **Grade 74**

(Whatman 40; S&S 589/2, 589/4)

9 cm, Pk/100. **0740-0900**\$34.85

15 cm, Pk/100. **0740-1500**\$60.00

Fast flow (220 mL/min), 15 µm retention, **Grade 55**

(Whatman 54, 541; S&S 1505)

9 cm, Pk/100. **0550-0900**\$40.00

15 cm, Pk/100. **0550-1500**\$63.25

Blotting/Chromatography/Electrophoresis Papers

0.19 mm thick, 125 mm/30 min absorption, **Grade 601**

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

46 x 57 cm, Pk/100. **6018-4657**\$157.00

0.34 mm thick, 140 mm/30 min absorption, **Grade 238**

(Whatman 3MM, S&S 598, 231^{cr})

20 x 20 cm, Pk/100. **2388-2020**\$80.00

46 x 57 cm, Pk/100. **2388-4657**\$305.00

0.42 mm thick, 150 mm/30 min absorption, **Grade 237**

(Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703;

Albet 2040B, FP591; Munktel 3S/H, FN7a)

46 x 57 cm, Pk/100. **2378-4657**\$352.00

0.95 mm thick, 180 mm/30 min absorption, **Grade 243**

(Whatman 17CHR; S&S 2668)

11 x 14 cm, Pk/100. **2438-1114**\$31.75

20 x 20 cm, Pk/100. **2438-2020**\$66.15

2.48 mm thick, 240 mm/30 min absorption, **Grade 320**

(S&S 2727)

11 x 14 cm, Pk/50. **3208-1114**\$63.45

20 x 20 cm, Pk/50. **3208-2020**\$93.05



Smooth tips &
polypropylene grips

Forceps, Filter, 11.5 cm. **300-025**\$4.40

Freezer Racks



Freezer Racks, stainless steel, **USA**

175V-13-2L

Restocking fee of 20% on all returned freezer racks listed below.

No returns on custom racks.

Horizontal (fits upright freezers)

Horizontal, Holds 2" Boxes

	1-9	≥10
9-Place (3L x 3H). 175HS-9-2	\$59.00	53.00
Dimensions: 16.31(L) x 5.5(W) x 6.5" (H)		
12-Place (3L x 4H). 175HS-12-2	\$78.00	71.00
Dimensions: 16.31(L) x 5.5(W) x 8.75" (H)		
15-Place (3L x 5H). 175HS-15-2	\$98.00	89.00
Dimensions: 16.31(L) x 5.5(W) x 11" (H)		
16-Place (4L x 4H). 175HL-16-2	\$105.00	95.00
Dimensions: 22(L) x 5.5(W) x 8.75" (H)		
20-Place (4L x 5H). 175HL-20-2	\$132.00	119.00
Dimensions: 22(L) x 5.5(W) x 11" (H)		

Horizontal, Holds 3" Boxes

9-Place (3L x 3H). 175HS-9-3	\$61.00	55.00
Dimensions: 16.31(L) x 5.5(W) x 9.5"(H)		
12-Place (4L x 3H). 175HL-12-3	\$80.00	73.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		



Horizontal Drawer (fits upright freezers)

Horizontal Drawer, Holds 2" Boxes

	1-9	≥10
9-Place (3L x 3H). 175HD-332	\$85.00	77.00
Dimensions: 16.5(L) x 5.5(W) x 7.125"(H)		
12-Place (3L x 4H). 175HD-342	\$111.00	101.00
Dimensions: 16.5(L) x 5.5(W) x 9.5" (H)		
15-Place (3L x 5H). 175HD-352	\$141.00	128.00
Dimensions: 16.5(L) x 5.5(W) x 11.812" (H)		
16-Place (4L x 4H). 175HD-442	\$151.00	137.00
Dimensions: 22(L) x 5.5(W) x 9.5" (H)		
20-Place (5L x 4H). 175HD-542	\$188.00	171.00
Dimensions: 27.25(L) x 5.5(W) x 9.5" (H)		

Horizontal Drawer, Holds 5" Boxes

	1-9	≥10
6-Place (3L x 2H). 175RD-325	\$74.00	67.00
Dimensions: 18.125(L) x 6(W) x 10.75" (H)		
8-Place (4L x 2H). 175RD-425	\$99.00	90.00
Dimensions: 24(L) x 6(W) x 10.75" (H)		

Horizontal Drawer, Holds 15 mL tubes

160-Place (80L x 2H). 175HD-32ML-15	\$99.00	90.00
Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)		
220-Place (110L x 2H). 175HD-42ML-15 ...	\$129.00	118.00
Dimensions: 22(L) x 5.5(W) x 9.5" (H)		

Horizontal Drawer, Holds 50 mL tubes

84-Place (42L x 2H). 175HD-42ML-50	\$129.00	118.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		

Continued