


LAB, SURGICAL & DENTAL INSTRUMENTS

Legend

- Economy:** stainless steel; forged and finished in Pakistan
Standard: surgical-grade (OR) German steel; forged in Germany and finished in Pakistan and/or Germany
German: surgical-grade (OR) German steel; made entirely (forged and finished) in Germany
Roboz: Roboz Surgical. (RS-)
-TC: Tungsten Carbide inlays; identified by gold handles

Forceps (Tweezers) & Tongs



General

Forceps

Fine Point, serrated handles, stainless steel

Straight, serrated tips & handles

- 10 cm (4"), Economy, guide pin. **26-275** \$10.75
- 10 cm (4"), Standard, guide pin. **40-720** \$16.50
- 10 cm (4"), German, guide pin. **420-802** \$32.80
- 10 cm (4"), Roboz, 0.8 mm tip, guide pin. **RS-5130** \$69.30
- ★ 11.5 cm (4½"), Economy, guide pin. **300-130** \$2.80
- 11.5 cm (4½"), Economy, No Pin. **20-373** \$3.35
- 20.5 cm (8"), Economy, guide pin. **300-141** \$8.10
- 25.4 cm (10"), Economy, guide pin. **300-143** \$9.20

Straight, Smooth Tips, serrated handles

- 10 cm (4"), Roboz, 0.8 mm tip, guide pin. **RS-5140** \$69.30
- 20.5 cm (8"), Economy, guide pin. **300-111** \$9.20

Angled, serrated tips

- 10 cm (4"), Roboz, 0.8 mm tip, guide pin. **RS-5139** \$73.25
- 15.5 cm (6"), Economy, guide pin. **3600** \$4.00

Half Curve (45°), serrated tips & handles

- 10 cm (4"), Economy, guide pin. **26-280** \$10.75
- 10 cm (4"), Standard, guide pin. **40-722** \$16.50
- 10 cm (4"), Roboz, 0.8 mm tip, guide pin. **RS-5135** \$73.25
- ★ 11.5 cm (4½"), Student, No Pin. **4190** \$3.30
- ★ 11.5 cm (4½"), Economy, guide pin. **300-131** \$2.90

Full Curve (90°), serrated tips & handles

- 10 cm (4"), Economy, guide pin. **24-665S** \$6.10
- 10 cm (4"), Roboz, 0.8 mm tip, guide pin. **RS-5137** \$73.25

Medium-Fine Point, serrated handles, stainless steel

- Curved, smooth tips, no pin, 15.2 cm (6"). **300-011** \$7.50

Medium Point (Thumb Dressing), serrated handles, stainless steel
 Straight, serrated tips, no guide pin

- 9 cm (3½"), Economy. **17-315** \$2.25
- 11.5 cm (4½"), Economy. **300-038** \$2.85
- 11.5 cm (4½"), Standard. **40-700** \$15.80
- 11.5 cm (4½"), Roboz, 1.5 mm tip. **RS-8100** \$33.00
- 15.5 cm (6"), Economy. **300-041** \$3.40
- 15.5 cm (6"), Standard. **40-710** \$12.60
- 15.5 cm (6"), Roboz. **RS-8106** \$34.00
- 16 cm (6¼"), Economy, black matte. **TWZ-703.74** \$4.50
- 20.5 cm (8"), Economy. **2240** \$4.05
- 20.5 cm (8"), Standard. **540-055** \$21.00
- 25.4 cm (10"), Economy. **2250** \$4.80
- 25.4 cm (10"), Standard. **540-060** \$31.50
- 30.5 cm (12"), Economy. **2260** \$5.85
- 51 cm (20"), Economy. **19-380** \$17.90

Exclusive! Extra wide tips ideal for removing specimens from storage containers.

Straight, Smooth Tips, guide pin

- 15.5 cm (6"), economy. **17-320SM** \$3.60

Curved, serrated tips, guide pin

- 11.5 cm (4½"), economy. **17-655** \$3.30



Smooth handles do not have guide pin.

Gooseneck, blunt serrated tips, 15 cm (6"). **300-001** \$5.30



Student Grade, stainless steel, 11.5 cm (4½"). **18-915** \$3.00


straight, serrated tips & handles; does not have guide pin

Corrosion & heat resistant ABS plastic is nonconductive & noncontaminating. Will not scratch optics.



Plastic, straight, serrated tips & handles, 13 cm. **97-3901** \$0.45

Micro



#5 Style w/Hook

Forceps, Micro, smooth tips* & handles, no guide pin

Economy, stainless steel

- #2A Style, straight, 12 cm. **57.0752A** \$7.45
- #3 Style, straight, 12 cm. **19-361** \$5.20
- ★ #5 Style, straight, 11.5 cm. **19-364** \$5.20
- #5 Style w/Hook, straight, 11 cm. **19-364HK** \$22.95
- #7 Style, curved, 11.5 cm. **19-366** \$5.20
- #7B Style, curved, *serrated, nonmag, 11.5 cm. **57.838**... \$14.35



Style #3 (very sharp)



Style #5 (needle sharp)



Style #7 (sharp)

Forceps, Micro (cont)

Dumont Economy, Inox, magnetic, Swiss

- #3 Style, straight, 12 cm (4 3/4"). **TWZ-301.16** \$33.00
- #5 Style, straight, 11 cm (4 1/4"). **TWZ-301.22** \$40.40
- #7 Style, curved, 11.5 cm (4 1/2"). **TWZ-301.28** \$37.55
- SS Style, nonmag, straight, 13.5 cm (5 1/4"). **TWZ-302.56** \$31.00

Dumont, Inox, high-precision, magnetic, Swiss

- #3 Style, straight, 12 cm (4 3/4")
 - Standard Tips (0.17W x 0.10 mm). **RS-5042** \$44.55
- #5 Style, straight, 11 cm (4 1/4")
 - Standard Tips (0.10W x 0.06 mm). **RS-5045** \$38.60
 - Biologie Tips (0.05W x 0.01 mm). **RS-5015** \$44.55
- #55 Style, light shanks, straight, 11 cm (4 1/4")
 - Biologie Tips (0.05W x 0.01 mm). **RS-5010** \$51.50
- #7 Style, curved, 11.5 cm (4 1/2")
 - Standard Tips (0.17W x 0.10 mm). **RS-5047** \$48.50




Dumont, Dumoxel, standard tips (0.10W x 0.06 mm), Swiss
#5/45 Style, tips angled 45°, 11 cm (4 1/4"). **RS-5058** \$46.55

Ideal-tek, stainless steel, nonmagnetic, Swiss

- SS Style, 14 cm (5 1/2"). **TWZ-312.56** \$17.95

Specialty



Forceps, stainless steel

- Adson Dressing, serrated tips & handles, 12 cm. **17-555** \$8.55



- Balance, curved plastic tips, 10 cm (4"). **300-006** \$5.25



- Bishop Harmon, fine serrated tips, 9 cm (3 1/2"). **4211** \$18.00



Cutting, economy

- Style 15A, 11.5 cm (4 1/2"). **TWZ-801.60** \$9.95



- Filter, smooth tips & polypropylene grips, 11.5 cm. **300-025** \$4.80



- Gerald Dressing, 17.5 cm (7"). **19-367** \$6.35



Graefe Eye Dressing, serrated tips & handles, 7 cm (2 1/2")

- Straight, standard. **4215** \$18.00
- Curved, standard. **4216** \$18.00



- Insect Pinning, rounded & serrated tips, 11 cm. **17-682** \$4.15



Paddle Tip, grid-serrated pads, economy

- 22 x 22 mm Pad, standard, 16 cm (6 1/4"). **33-1045** \$14.50
- 30 x 35 mm(W) Pad, 15.2 cm (6"). **33-1045B** \$21.00
- 30 x 40 mm(W) Pad, 15.2 cm (6"). **33-1045A** \$21.00



Wide Handle (Paddle Handle), economy

serrated and perforated grips

Straight, Serrated tips

- 14 cm (5 1/2"), 1.5 mm tips. **15-600** \$25.00

Straight, Smooth tips

- 14 cm (5 1/2"), 0.6 mm tips. **15-625** \$25.00
- 14 cm (5 1/2"), 1.5 mm tips. **15-585** \$25.00
- 17.5 cm (7"), 1.5 mm tips. **15-550** \$25.00
- 17.5 cm (7"), 3 mm tips. **15-500** \$25.00

Angled, Serrated tips

- 14 cm (5 1/2"), 1.5 mm tips. **15-595** \$25.00

Angled, Smooth tips

- 14 cm (5 1/2"), 1.5 mm tips. **15-590** \$25.00
- 17.5 cm (7"), 1.5 mm tips. **15-580** \$25.00

Spade Tips



Angled (Kuehne), 11 cm (4 1/4")

- Stainless Steel, economy. **33-1041** \$4.15
- Teflon® Coated, Bel-Art. **379420000** \$35.60
- Straight, stainless steel, 10.5 cm, economy. **57.926** \$8.45



Square Point Tips, straight

- Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379430000** \$39.40

Locking



Forceps, stainless steel

Self-Closing (cross-locking), squeeze handles to open tips

Fine Point, smooth tips

Economy, 13 cm (5 $\frac{1}{8}$ "). **300-066** \$5.45

Peer-Vigor, 11 cm (4 $\frac{1}{4}$ "), Swiss. **57.324** \$42.40

Medium Point, serrated tips, 17 cm (6 $\frac{5}{8}$ "). **TWZ-948.00** \$3.95



Slide-Locking, grid-serrated tips

Fine Point

Economy, 14 cm (5 $\frac{1}{2}$ "). **57.0357A** \$9.00

Dumont, 16 cm (6 $\frac{1}{4}$ "). **TWZ-101.72** \$45.75

Medium Point

Economy, 16 cm (6 $\frac{1}{4}$ "). **203** \$9.00

Dumont, 16 cm (6 $\frac{1}{4}$ "). **TWZ-101.76** \$43.60

Snap-Closing, angled, fine points, 15.5 cm (6"). **3610** \$5.80

serrated tips & handles, guide pin

Tissue



Forceps, stainless steel

Adson-Brown, 9 x 9 side grasping teeth, 12 cm. **17-574** \$16.05


Adson Tissue, 1 x 2 teeth, 12 cm (4 $\frac{3}{4}$ "). **17-560** \$8.55

Iris Tissue, 1 x 2 teeth

Straight, economy, 10 cm (4"). **26-295** \$10.75

Angled, economy, 11.5 cm (4 $\frac{1}{2}$ "). **19-350ANG** \$5.00

Half Curve, standard, 10 cm (4"). **45-822** \$17.60



Russian Tissue, circular serrated tips, 15.5 cm (6"). **2275** \$5.25

Tissue, 1 x 2 teeth, 14 cm (5 $\frac{1}{2}$ "). **17-145** \$8.55

Tongs





Tongs, Crucible, serrated jaws, stainless steel

24.5 cm (9 $\frac{3}{4}$ "). **320-101** \$12.50

51 cm (20"). **320-120** \$30.15

Hemostats & Clamps



Hemostats, serrated jaws

Hartman Mosquito, stainless steel

Straight

Economy, 9 cm (3 $\frac{1}{2}$ "). **12-180S** \$6.05

Standard, 9 cm (3 $\frac{1}{2}$ "). **20-352** \$21.30

Curved

Economy, 9 cm (3 $\frac{1}{2}$ "). **306-012** \$9.10

Standard, 10 cm (4"). **20-360** \$18.75↓

Halstead Mosquito, 12.5 cm (5"), stainless steel

Straight

Economy. **14-065** \$10.00

Standard. **20-501** \$22.45

Curved

Economy. **14-070** \$10.00

Standard. **20-519** \$22.85

Kelly, stainless steel

Straight

Economy, 14.5 cm (5 $\frac{3}{4}$ "). **306-016** \$8.20

Standard, 14.5 cm (5 $\frac{3}{4}$ "). **20-550** \$23.60

Curved

Economy, 14 cm (5 $\frac{1}{2}$ "). **306-017** \$9.90

Standard, 14.5 cm (5 $\frac{3}{4}$ "). **20-568** \$23.95

Rochester-Pean, 20.5 cm (8"), stainless steel

Straight

Economy. **1920** \$6.50

Standard. **21-977** \$33.95

Curved


Economy. **1930** \$6.50

Standard. **21-979** \$33.95

Polypropylene, 12.5 cm (5"). **HEM-300.00** \$3.25



Specimen Forceps, 25 cm (9 $\frac{3}{4}$ "), serrated jaws. **300-500** \$18.70



Sponge Forceps, Foerster, serrated jaws, stainless steel

Straight

17.5 cm (7"), economy. **14-660** \$17.50

24 cm (9 $\frac{1}{2}$ "), economy. **14-670** \$21.10

Curved, 17.5 cm (7"), economy. **14-665** \$17.50

Clamps



300-120



300-121

Dieffenbach Serrefine, stainless steel, economy

Straight, serrated jaws and handles

60 mm (2³/₈"). **300-120** \$7.40

Curved, serrated jaws and handles

48 mm (1⁷/₈"). **12-248** \$9.15

60 mm (2³/₈"). **300-121** \$7.40



Riveted handle
opens for cleaning.

Johns Hopkins Bulldog, Roboz

Straight, 38 mm (1¹/₂"). **RS-7440** \$88.10



RS-5452

Schwartz Vessel Clips, 1.7 x 8 mm jaws, 2.5 cm long, Germany

Straight Jaws, Roboz. **RS-5452** \$36.00


Slightly Angled Jaws, Roboz. **RS-5457** \$36.00

8.6 cm shaft; 5 mm serrated jaw



Micro Alligator Forceps, stainless steel. **2769** \$50.00

Needle Holders & Wound Closure




Crile-Wood



Olsen-Hegar



Webster



Needle Holders, stainless steel

Baumgartner, cross-serrated jaws, 14 cm (5¹/₂"). **16-285S** \$11.05

Castroviejo, cross-serrated jaws w/lock, 14 cm. **15-302** \$50.00

Crile-Wood, cross-serrated jaws, 15 cm (6")

Economy, central groove. **16-293S** \$13.45

Standard. **32-485** \$33.75

German. **320-366** \$60.80

Derf, 12 cm (4³/₄"), economy. **16-238** \$9.70

serrated jaws with central groove

Olsen-Hegar, cross-serrated jaws with scissors, 14 cm (5¹/₂")

Economy, central groove. **16-296** \$14.30

Economy, No Groove. **16-296A** \$14.45

Standard, central groove. **32-500** \$35.25

Webster, smooth jaws, 12 cm (4³/₄")

Economy. **15-035** \$17.45

Standard. **32-101** \$30.75


Tungsten Carbide. **0817** \$45.00

German. **320-101** \$60.50

Michel Wound Closure System



Michel Clip Applier/Remover, double-ended. **3700** \$13.75
11.5 cm (4¹/₂"), stainless steel



Michel Wound Clips, stainless steel

11 mm, Bx/1000. **380-519** \$74.00

14 mm, Bx/1000. **380-543** \$89.00

16 mm, Bx/1000. **380-568** \$97.00

18 mm, Bx/1000. **380-584** \$98.00

22 mm, Bx/1000. **380-626** \$116.00

Reflex Wound Closure System

CellPoint Scientific, USA

Reflex Clip Applier

For 7 mm Clips. **204-1000** \$165.00

★ For 9 mm Clips. **202-1000** \$165.00


Reflex Wound Clips, nonsterile, stainless steel

7 mm, 5 packs of 20 clips, Bx/100. **203-1000** \$56.00

★ 9 mm, 5 packs of 20 clips, Bx/100. **201-1000** \$46.00

Reflex Clip Remover, fits 7 mm & 9 mm clips. **205-1000** \$99.00

Inoculating Loops/Needle & Cell Spreaders



Inoculating Loops

Disposable

1 µL, Sterile Pk/20, PP, Cs/1000. **HS81121A** \$69.75

1 µL, Sterile Pk/1, PS, Cs/250. **L200-11** \$106.00

Nichrome, 6 mm twisted loop of 26-gauge wire. **600INL703** \$2.00
21 cm aluminum handle

Inoculating Needle. 600INN803 \$2.50

75 mm, 26-gauge nichrome wire; 20 cm aluminum handle

16 mm bar width;
20 cm handle



Cell Spreader, Stainless Steel. 377360006 \$26.10


Microtome Blades

Microtome Blades - see Microtome Blades in Tissue (Embedding) &
Biopsy Cassettes section

Punches & Tags



Ear Punch, thumb grip, 2 mm hole, plated steel, **USA. 1538**..... \$10.35



Ear Tag Applicator, USA. 1005S1 \$65.80

Ear Tags, nickel-copper alloy

*Specify 3-digit number sequence, e.g., 101-200

**Minimum order of \$60.00 required for tags (any combination)

Stamped

Pk/100. **1005-1-** - * \$41.30
 Pk/500. **1005-1-** - * \$133.00
 Pk/1000. **1005-1-** - * \$188.00

Laser Etched

Pk/100. **1005-1L1-** - * \$60.05
 Pk/500. **1005-1L1-** - * \$226.00
 Pk/1000. **1005-1L1-** - * \$274.00


Razor Blades



Single edge blades are 2.25" long.
 Manufactured by Medical
 Specialties; formerly Pilling Weck
 Surgical Inc. (WECPREP™).


Razor Blades, prep, **USA**, Bx/50. **450205MS** \$45.00

Retractors & Specula




Retractors

Alm, self-retaining, blunt 4 x 4 prongs
 8 cm (3"), economy. **44-1519** \$30.00
 10 cm (4"), economy. **44-1520** \$30.00



Barraquer, wire loop, 4 cm (1½")
 4 mm Loop, economy. **33-1222** \$20.00
 8 mm Loop, economy. **33-1221** \$20.00



Weitlaner, self-retaining, 11.5 cm (4½")
 Sharp 2 x 3 Prongs, economy. **44-1505** \$60.00
 Blunt 2 x 3 Prongs, economy. **44-1501** \$60.00


Vienna Nasal
 Speculum



Specula

Graves Vaginal, medium, economy. **3310** \$18.00
 Vienna Nasal, child, 13 cm (5⅛"). **2722** \$10.50
 lightweight, stainless steel

Rongeurs & Bone Instruments



Blumenthal Rongeur

Rongeurs, stainless steel

Blumenthal, 15 cm (5¾")
 Economy. **410-571** \$48.40
 Standard. **650-571** \$134.00
 Blumenthal, Slightly Angled Jaws, 15 cm (5¾"). **9571** \$37.50
 Mini-Blumenthal, 11.5 cm (4½"). **410-560** \$58.80
 Friedman, 14 cm (5½")
 Economy. **410-574** \$48.40
 Standard. **650-574** \$134.00

Liston Bone Cutter




Bone Cutters, stainless steel, economy

Liston
 14 cm (5½"), 20 mm blades. **320-086** \$25.55
 21.5 cm (8½"), 30 mm blades. **320-087** \$31.25
 Scissor Style, 10 cm (4"), 18 mm blades. **320-025** \$15.15
 straight, strong tips




Bone Nipper, 10 cm (4") long, 6 mm cutting edge. **B286** \$12.00
 angled offset jaws


Continued



Bone Shears, Constable, 23 cm (9") long, 50 mm (2") cutting jaws
 Economy. **33-1001**\$45.00
 German. **965-620**\$160.00



Postmortem Hammer, hook handle, 25 cm (10") long
 Economy. **33-1032**\$47.50
 Standard. **74-043**\$80.00



16 mm ($\frac{5}{8}$ ") cutting blade;
 9 cm ($3\frac{1}{2}$ ") long

Skull Breaker, Virchow, economy. **33-1011**\$42.40

Scalpels



Scalpels

	<u>bx/10</u>	<u>cs/50 bx</u>
Disposable, economy, sterile, gray plastic handle		
★ #10 carbon steel blade. 9410	\$11.65	516.00
#11 carbon steel blade. 9411	\$11.65	516.00
#22 carbon steel blade. 9422	\$11.65	516.00



Miniature, plastic thumb handle, nonsterile, England
 #11P stainless steel blade, each. **2201111**\$1.35




Smooth Handle



Fluted Handle

Reusable, one-piece stainless steel, $1\frac{1}{2}$ " blade
 Smooth Handle, 15.5 cm (6"). **14-350**\$3.15
 Fluted Handle, narrow blade, 15.5 cm (6"). **14-360**\$3.15



Scalpel Blades, individually wrapped

Carbon Steel, Economy, Bx/100	\$23.20
#10 blade. 9310	#21 blade. 9321
#11 blade. 9311	#22 blade. 9322
#15 blade. 9315	#23 blade. 9323
#20 blade. 9320	#24 blade. 9324
Carbon Steel, Germany, nonsterile	
#60 blade, 25 strips of 4, Bx/100. 2200260	\$59.00
#70 blade, 25 strips of 4, Bx/100. 2200270	\$70.00
Carbon Steel, Swann-Morton, nonsterile	
#60 blade, Bx/50. 22002SM60	\$39.60
Stainless Steel, Economy, Bx/100	\$21.95
#10 blade. 36-010S	#21 blade. 052-610
#11 blade. 052-602	#22 blade. 052-612
#20 blade. 36-020S	
Stainless Steel, Lance, Bx/100	\$47.50
#11 blade. 92111	#15 blade. 92115
Stainless steel, Swann-Morton, sterile	
#22 blade, Bx/100. 22001SM22	\$52.00
(RS-9861-19)	

Continued



Safely remove and store up to 100 used blades in this disposable box.

Scalpel Blade Remover. 22048..... \$4.35




Stainless Steel, Solid, #3 Style



Stainless Steel, Solid, #4 Style



Stainless Steel, Solid, #7 Style



Stainless Steel, Hollow



PVC Coated Handle, #3 Style



Plastic Handle, #5/6/8 Styles

Scalpel Handles, reusable

Stainless Steel, solid

#3 Style, accepts blades #9-16

Economy. **02-490** \$3.45

Standard. **10-305**..... \$9.25↓

#4 Style, accepts blades #20-26

Economy. **02-492** \$3.45

Standard. **10-310**..... \$9.25↓

#7 Style, accepts blades #9-16

Economy. **02-494** \$5.35

Standard. **10-315**..... \$18.50

Stainless Steel, Hollow

#8 Style, accepts blades #60 & #70. **HAN8** \$11.00

PVC Coated Stainless Steel. *Not autoclavable*

#3 Style, accepts blades #10-17. **KNF-225.00** \$6.00

Plastic, contoured to fit palm of hand

#5 Style, accepts blades #10-17. **KNF-221.00** \$5.25

#6 Style, accepts blades #20-26. **1605** \$5.00

#8 Style, accepts blades #60 & #70

Bard-Parker™. **371080**..... \$44.95

Scissors




Barraquer


Scissors, Spring, ambidextrous, stainless steel

Barraquer (butterfly), 6 cm long, economy. **3737** \$34.50
10 mm angled-on-flat blades

36205/5-855



M3720



Castroviejo, straight, sharp tips, economy

7 cm (2¾"), 5 mm blades. **36205** \$45.00

Effective cutting edge is 3 mm; tip-to-tip is 3 mm.

8.5 cm (3¼"), 5 mm blades. **5-855** \$45.00

Effective cutting edge is 4 mm; tip-to-tip is 3 mm.

11.5 cm (4½"), 17 mm blades. **M3720** \$26.00

Effective cutting edge is 12 mm.



Uniband™

Lifetime guarantee against rust and defects.



Featherlite™

Uniband™, sharp tips, Laschal Surgical, **USA**

Straight, 12 mm Blades. **LA-4** \$32.00

Straight, 20 mm Blades. **LA-1** \$32.00

Straight, 20 mm Serrated Blades. **LA-2** \$32.00

Curved, 12 mm Blades. **LA-4C** \$32.00

Curved, 20 mm Blades. **LA-1C** \$32.00

Featherlite™, Castroviejo, sharp tips, Laschal Surgical, **USA**

Greater tactile sensitivity than Uniband™ scissors.

Straight, 12 mm Blades. **MPF-4** \$80.80

Curved, 12 mm blades. **MPF-4C** \$80.80

Scissors, Bonn, 15 mm blades, 9 cm (3½"), Germany

Straight, Roboz. **RS-5840** \$107.00

Curved, Roboz. **RS-5841** \$107.00



Straight

Curved

Scissors, Delicate, stainless steel

Straight, sharp tips

9 cm (3½"), 20 mm blades, Roboz. **RS-5880** \$102.00


10.5 cm (4"), 22 mm blades, Standard. **13-003** \$21.30

Curved, sharp tips

9 cm (3½"), 20 mm blades, Roboz. **RS-5881** \$106.00

10 cm (4"), 22 mm blades, Standard. **13-059** \$22.45

Continued



Scissors, Dissecting, 11.5 cm (4½"), stainless steel

Straight, sharp tips

25 mm Blades, Economy. 310-040	\$4.95
25 mm Blades, Standard. 13-002	\$19.00
24 mm Blades, Roboz. RS-5912	\$91.10
35 mm Blades, Economy. 310-045	\$4.80


Curved, sharp tips

25 mm Blades, Economy. 310-035	\$4.80
25 mm Blades, Standard. 13-010	\$19.40
24 mm Blades, Roboz. RS-5913	\$94.05




Angled (angled-on-side), stainless steel

11.5 cm (4½"), Economy. 310-066	\$6.50
12 cm (4¾"), 30 mm blades, Roboz. RS-6708	\$99.95



Scissors, Deaver, 14 cm (5½"), stainless steel

Straight, S/S (sharp/sharp) tips, Standard. 10-632	\$21.80
Straight, S/B (sharp/blunt) tips, Standard. 10-630	\$21.30




Scissors, Mayo, semi-blunt tips, stainless steel

Straight

14 cm (5½"), Economy. 310-086	\$9.90
14.5 cm (5¾"), Standard. 15-575	\$22.55
14.5 cm (5¾"), Roboz. RS-6870	\$87.10
17.5 cm (6¾"), Economy. 310-088	\$12.15
17.5 cm (6¾"), Standard. 15-579	\$26.95
17.5 cm (6¾"), Standard-TC. 16-854	\$58.50
23 cm (9"), Standard. 15-585	\$42.35

Curved

14 cm (5½"), Economy. 310-087	\$10.40
14.5 cm (5¾"), Standard. 15-576	\$22.95
14 cm (5½"), Roboz. RS-6871	\$87.10



Scissors, Metzenbaum, stainless steel

Straight

15 cm (5¾"), Economy. 12-225	\$14.00
-------------------------------------	---------

Curved

14.5 cm (5¾"), Economy. 12-230	\$14.00
14.5 cm (5¾"), Standard. 15-552	\$23.05
14 cm (5½"), German. 105-569	\$46.20


Blunt/Blunt Tips



Sharp/Blunt Tips



Sharp/Sharp Tips



Scissors, Operating, stainless steel

Straight, S/S (sharp/sharp) tips

11.5 cm (4½"), Standard. 10-456	\$17.85
14 cm (5½"), Economy. 12-055	\$8.55
14 cm (5½"), LR (large rings), Economy. 11-107SSLR	\$6.05
14 cm (5½"), Standard. 10-552	\$21.30
14.5 cm (5¾"), German. 101-550	\$36.50
15.5 cm (6"), Economy. 12-060	\$10.75
16.5 cm (6½"), Economy. 12-065	\$10.75
16.5 cm (6½"), Standard. 10-602	\$24.75
16.5 cm (6½"), German. 101-600	\$40.20

Straight, S/B (sharp/blunt) tips

11.5 cm (4½"), Economy. 12-015	\$8.55
12 cm (4¾"), Standard. 10-453	\$17.85
12 cm (4¾"), German. 100-453	\$33.60
12 cm (4¾"), 30 mm blades, Roboz. RS-6700	\$99.95
14 cm (5½"), Economy. 12-025	\$8.55
14 cm (5½"), Standard. 10-550	\$21.30
14 cm (5½"), German. 100-552	\$36.50
16.5 cm (6½"), Economy. 310-111	\$11.00
16.5 cm (6½"), Standard. 10-600	\$24.75
16.5 cm (6½"), 60 mm blades, Roboz. RS-6824	\$76.25

Straight, B/B (blunt/blunt) tips

11.5 cm (4½"), Economy. 310-100	\$7.40
14 cm (5½"), Economy. 12-085	\$8.55
14 cm (5½"), Standard. 10-558	\$21.30
14 cm (5½"), German. 102-558	\$36.50
16.5 cm (6½"), Economy. 12-095	\$10.75
16.5 cm (6½"), Standard. 10-608	\$24.75
16.5 cm (6½"), 60 mm blades, Roboz. RS-6828	\$76.25

Curved, S/S (sharp/sharp) tips

14 cm (5½"), Economy. 12-155	\$8.55
14.5 cm (5¾"), Standard. 10-566	\$13.25
16.5 cm (6½"), Economy. 12-165	\$10.75
16.5 cm (6½"), Standard. 10-616	\$25.25

Scissors, Operating (cont.)

Curved, S/B (sharp/blunt) tips

11.5 cm (4½"), Economy. 310-101	\$7.65
14 cm (5½"), Economy. 12-125	\$8.55
14.5 cm (5¾"), Standard. 10-560	\$15.40 ↓
14 cm (5½"), German. 100-560	\$37.90
16.5 cm (6½"), Economy. 12-135	\$10.75
16.5 cm (6½"), Standard. 10-610	\$25.25
16.5 cm (6½"), 60 mm blades, Roboz. RS-6849	\$81.20


Curved, B/B (blunt/blunt) tips

11.5 cm (4½"), Economy. 310-102	\$7.85
14 cm (5½"), Economy. 12-185	\$8.55
14 cm (5½"), Standard. 10-568	\$21.80
16.5 cm (6½"), Economy. 12-195	\$10.75
16.5 cm (6½"), Standard. 10-618	\$25.25



Scissors, Stevens Tenotomy, 11.5 cm (4½"), stainless steel

Straight, blunt tips, Standard. **13-520** \$24.75




Scissors, Strabismus, stainless steel

Straight

10 cm (4"), Economy. 26-105	\$10.60
11.5 cm (4½"), Standard. 13-524	\$23.65

Curved

11 cm (4¼"), Economy. 26-110	\$10.60
11.5 cm (4½"), Standard. 13-526	\$24.15



Scissors, Bandage, stainless steel, economy

Knowles, straight, sharp probe tips, 14 cm (5½"). **1070** \$6.25

Lister, angled, blunt probe tips, 14 cm (5½"). **1050** \$5.00

Can be used as enterotomy scissors.



Cut through everything from wire to bone with these durable scissors. Bottom blade is serrated and has a footplate for safety. Blades are stainless steel & handles are plastic. Autoclavable to 120° C.

Utility Shears, 18 cm (7¼"), economy. **310-341** \$5.70

Scrapers

Remove paraffin, ice, etc., safely from your laboratory equipment with these plastic scrapers.



Scrapers, plastic

5 cm (2") Wide Blade. 311-245	\$3.70
10 cm (4") Wide Blade. 311-250	\$4.70
15 cm (6") Wide Blade. 311-255	\$5.70

Spatulas & Scoops



Spatulas, Hardwood Handle, stainless steel blade

10 cm (4") Blade. 14-365B	\$5.40
15 cm (6") Blade. 14-365D	\$7.05
20 cm (8") Blade. 20-497	\$6.05



Spatulas, Heavy-Duty, stainless steel

Round End/Square End, both flat

30.5 cm (12"). 320-425	\$16.10
46 cm (18"). 320-430	\$21.50



Spatulas, Micro, stainless steel

Square End/Spoon End

19 cm (7½"). 320-411	\$4.50
23 cm (9"). 320-414	\$4.65

Round End/Tapered End, 18 cm (7") long. **320-370** \$5.00

Round end: 8 x 51 mm; Tapered end: 6 x 51 mm

V-Shaped End/Round End, 16 cm (6¼") long. **HS15906** \$7.45



V-Shaped (Reagent Digger), 16.5 cm (6½"). 600LSCSS6	\$2.60
Handle, wooden, 13 cm long, each. 600LSHND1	\$3.00

Teasing Needles, Probes & Pin Vises

97-1302



Teasing Needles

Angular, wooden handle, 14.5 cm (5 3/4"), Pk/12. **97-1302** \$4.95
Straight, plastic handle, 14.5 cm (5 3/4"), Pk/12. **97-1313** \$4.65



Stainless Steel, straight, 14 cm. **M28** \$5.00

Blunt tip is angled and tapered. Hexagonal handle is 6 mm wide.



Probes, Light Mall, 15.5 cm (6")

Nickel Plated Steel. **15-800** \$2.90
(5800)
Stainless Steel. **15-800SS** \$3.80



Probes, Bowman Lacrimal, double-ended

Economy, sterling silver, 12.7 cm (5") long
Size 0000-000 (~0.4 – 0.5 mm). **24-706** \$11.65
Size 00-0 (~0.6 – 0.7 mm). **24-708** \$11.65
Size 1-2 (~0.9 – 1.0 mm). **24-710** \$11.65
Size 3-4 (~1.1 – 1.2 mm). **24-712** \$11.65




320-216A

Probe/Seekers, stainless steel

Double-Ended, 15.5 cm. **320-216A** \$4.90
One end tapered to a point; one end curved with a rounded tip.
Single End, 12.5 cm (5"). **320-217A** \$4.75
Curved, fine point tip with hexagonal handle.



Probe with Bulb Tips, 25 cm (10"), stainless, German. **541-807** \$12.00 ↓



Probe with Eye & Bulb Tip, flexible, 15.5 cm (6"). **23-496** \$2.25

PIN-220.00



PIN-212.00




Hold needles, probes, wires, drill bits, etc.

Pin Vise, Swivel-Head, 0-3.3 mm opening. **PIN-220.00** \$8.35
Two double-ended collets are provided; one in screw chuck, one inside handle.

Pin Vise, Mini, 0-0.8 mm opening. **PIN-212.00** \$5.20
blackened tool steel


Trocar



Trocar, 4-in-1 set, 13 cm, stainless steel. **1596** \$17.50

Accessories/Maintenance/Storage

105CB1218



1054319

Cutting Boards, white high-density polyethylene

10 x 6 x 1/2" (25.4 x 15.2 x 1.3 cm). **1054319** \$6.95
18 x 12 x 1/2" (45.7 x 30.5 x 1.3 cm). **105CB1218** \$16.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1 1/2 x 3" opening. **DEM-345.00** \$49.95



1257371



36900

Tools not included.

Dissecting Cases

6-Loop with snap fastener, black leatherette. **1257370** \$1.30
12-Loop with snap fastener, black leatherette. **1257371** \$1.40
13-Loop with zipper, brown leatherette. **36900** \$7.80
Features red velveteen interior, protector flap, and ID pouch.

Dissecting Kits

*Inquire about customizing a kit for your surgical procedure.
Kit components are discounted 10%.*

Basic kit. 38501.....\$5.30

Contents: Dissecting Case, 6-loop (1257370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

Mouse kit. 38551.....\$17.85

Contents: Dissecting Case, 6-loop (1257370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

Dissecting Pans, stainless steel (without wax)

20⁷/₈ x 12¹³/₁₆ x 2¹/₂". **AM20029**.....\$21.90
25³/₄ x 17³/₄ x 2¹/₄". **AM68357**.....\$64.90



Insect Pins, 38 mm (1½") long, nylon head, Czech Republic

Black Enameled, Pk/100.....\$7.00

#000, 0.25 mm dia. 355425	#1, 0.40 mm dia. 355440
#00, 0.30 mm dia. 355430	#2, 0.45 mm dia. 355445
#0, 0.35 mm dia. 355435	#3, 0.50 mm dia. 355450

Stainless Steel, Pk/100.....\$8.75

#000, 0.25 mm dia. 355525	#1, 0.40 mm dia. 355540
#00, 0.30 mm dia. 355530	#2, 0.45 mm dia. 355545
#0, 0.35 mm dia. 355535	#3, 0.50 mm dia. 355550




Instrument Cup, silicone, 3½" opening, 2½" deep. **150-1005**.....\$57.00



Instrument Storage Boxes, 10 x 4 x 2", 18-8 stainless steel

Small, 7 x 3 x 1½". **35-406**.....\$26.60
Large, 10 x 4 x 2". **35-408**.....\$35.00

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. **4461**.....\$0.95

White. **4460**.....\$0.55

Stainless Steel, 12 mm wide, **USA. 605-065**.....\$6.25



Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

Sterilizer, Hot Bead, Germinator™ 500, USA. 5-1450.....\$426.00



Glass Beads (replacement), 1.5 mm, Bg/12 oz. **5-1455**.....\$30.00



Safely incinerate microorganisms in 5-7 seconds at 815°C (1500°F) within this ceramic, stainless steel-encased chamber. Compact and flameless operation allows this unit to be used in anaerobic and aerobic chambers. Ideal for sterilizing loops and needles, glass tube/pipette mouths, and various glass and metal instruments. 1-year warranty.

Sterilizer, Infrared Heat, BactiZapper™. B1000.....\$476.00

6100A



Sterilizing Trays, plastic, autoclavable, PST, **USA**

B/L/M (base, lid & mat) included.

6 x 2½ x ¾" (15.2 x 6.4 x 1.9 cm). **2600A**.....\$23.00
(500252)

★ 6½ x 4 x ¾" (16.2 x 10.2 x 1.9 cm). **4600A**.....\$33.00

10 x 6 x ¾" (25.4 x 15.2 x 1.9 cm). **6100A**.....\$67.00
(500254)



Instrument Trays

Polypropylene, deep w/flat cover, autoclavable

8½ x 6 x 2¾" (22 x 15 x 7 cm), Pk/2. **60081736**.....\$25.15

Continued



15F



35-309

Stainless Steel

Shallow (Mayo), **USA**

15 $\frac{1}{8}$ x 10 $\frac{5}{8}$ x $\frac{3}{4}$ " (38.4 x 27 x 1.9 cm). **15F** \$51.65

Deep w/ Dome Cover

8 $\frac{3}{4}$ x 3 $\frac{1}{4}$ x 1 $\frac{1}{2}$ " (22.2 x 8.4 x 3.8 cm). **35-309** \$39.90

Deep – Tray Only, **USA**

8 $\frac{7}{8}$ x 5 x 2" (22.7 x 12.9 x 5.1 cm). **35-314NP** \$41.40



This natural Arkansas stone outlasts and produces a better finish than a man-made stone. Mounted in a wooden box. Dimensions: 4 x 2 x $\frac{1}{2}$ "

Sharpening Stone, fine grit (hard). **ABR-800.10** \$31.95

Includes a plastic cleaning basket.

Requires only 42 kHz of power.

Dimensions: 6 x 3 $\frac{1}{2}$ x 2"(D).



Ultrasonic Cleaner, 600 mL, mechanical timer. **CLN-520.20** \$45.95

Cleaning Solution, 1 Quart. **CLN-850.10** \$9.95

Biodegradable concentrate mixes 1:30.