

GENERAL LAB SUPPLIES & EQUIPMENT

Applicators




Applicators, 6" long

Cotton Tipped, wooden shaft

Nonsterile, 10 bg/bx, Bx/1000. 95-8702	\$15.30
Sterile Pk/1, Bx/100. 44024-106-1S	\$7.85
Sterile Pk/2, Bx/200. 44024-106-2S	\$10.50
Plain, wood, Bg/1000. 95-8701	\$9.35

Aprons



Protect yourself with this disposable 0.8-mil thick apron.

Aprons, disposable, polyethylene, 46"(L) x 28"(W)

0.8 mil Thick (0.02 mm), Bx/100. 105DA2846	\$16.65
1.75 mil Thick (0.044 mm), Bx/100. 59010485	\$38.00

Autoclaves



Autoclaves, benchtop, Benchmark (1-year warranty)

8 L, BioClave™ Mini. B4000-M	\$4,689.00
Chamber Dimensions: 6.7"(ID) x 12.5"(D)	
18 L, BioClave™ 18. B4000-18	\$6,414.00
Chamber Dimensions: 9.8"(ID) x 13.75"(D)	
28 L, BioClave™ 28. B4000-28-120	\$11,811.00
Chamber Dimensions: 12.6"(ID) x 16.4"(D)	

Bags



Collect biohazardous waste in this space-saving, autoclavable bag. Printed with a biohazard symbol and sterilization indicator patch. Measures 20.3 x 30.5 cm (8 x 12") and is 1.57 mil (0.04 mm) thick polypropylene.

34510044H



Bags, Biohazard, Benchtop, 8 x 12", Bx/500. 34510320 \$110.00

Biohazard Bag Holder, Benchtop. **34510044H** \$27.85

Epoxy-coated wire holder measures 9.1 x 12.8 x 21.4 cm.



Bags display a biohazard symbol and precautionary procedures. Sterilization indicator patch darkens when autoclaved. Polypropylene bags are bright red. Durable, extra heavy-duty bags provide maximum protection against dangerous tears and leaks.

Bags, Biohazard, Standard, 1.5 mil (0.038 mm) thick

14 x 19" (36 x 48 cm), Bx/200. 34510321	\$109.00
19 x 23" (48 x 58 cm), Bx/200. 131641923	\$157.00
25 x 35" (64 x 89 cm), Bx/200. 131642535	\$339.00

Bags, Biohazard, Heavy-Duty, 2.2 mil (0.055 mm) thick

12 x 18"(L), 9 L, Bx/200. 600UNU40105	\$70.25
12 x 24"(L), 12 L, Bx/200. 600UNU40106	\$79.75
18 x 22"(L), 16 L, Bx/200. 600UNU40103	\$97.50
25 x 35"(L), 37 L, Bx/200. 600UNU40104	\$164.00

Balances

Toploading

Benchmark
W3200-1200



Model	W3200-500	W3200-1200
Capacity (g)	500	1200
Readability (g)	0.001	0.01
Modes	g, oz, ct	g, oz, ct
Pan	11.5 cm (4.5")	16 cm (6.3")
Draft Shield	Yes, 3-door glass	No
Price	\$1,377.00	\$1,128.00

Warranty: 2 years

Analytical

Benchmark



Model	★ W3100-120	W3100A-120
Calibration	External (100 g weight provided)	Internal (Quick-Cal™ software)
Capacity (g)	120	120
Readability (g)	0.0001	0.0001
Modes	g, oz, lb, ct	g, oz, lb, ct
Pan	9 cm (3.5")	9 cm (3.5")
Price	\$1,997.00	\$2,119.00

Warranty: 2 years

Baths

BSH1002



Aluminum
Blocks



Dry Baths, digital, Benchmark (2-year warranty)

1-Block Capacity. BSH1001	\$549.00
2-Block Capacity. BSH1002	\$679.00
4-Block Capacity. BSH1004	\$1,029.00

Blocks sold separately – please specify configuration

Blocks, aluminum

Holds (48) 0.2 mL Tubes. BSW02	\$164.00
Holds (24) 0.5 mL Tubes. BSW05	\$164.00
Holds (24) 1.5 mL Tubes. BSW1500	\$164.00
Holds (24) 1.5/2.0 mL Tubes. BSW1520	\$164.00
Holds (20) 12/13 mL Tubes. BSW13	\$164.00
Holds (12) 15 mL Tubes. BSW15	\$164.00
Holds (5) 50 mL Tubes. BSW50	\$164.00



40 x 1.5/2.0 mL



15 x 50 mL



41 x 15 mL

Racks are only included with model B2000-4.

Water Baths, digital, MyBath™, Benchmark (2-year warranty)

4 L, 3 tube racks. B2000-4	\$1,119.00
8 L. B2000-8	\$1,194.00
12 L. B2000-12	\$1,409.00

Beakers

Polypropylene
Beakers



Beakers

Glass (borosilicate), low form, double scale

Economy

10 mL, no graduations, Pk/12 600BG10	\$35.50
50 mL, Pk/12. 600BG50	\$41.10
100 mL, Pk/12. 600BG100	\$43.15
150 mL, Pk/12. 600BG150	\$46.15
250 mL, Pk/12. 600BG250	\$49.85
400 mL, Pk/12. 600BG400	\$62.75
600 mL, Pk/6. 600BG600	\$45.80
800 mL, Pk/6. 600BG800	\$77.75
1000 mL, Pk/6. 600BG1000	\$74.00
2000 mL, each. 600BG2000	\$23.20
4000 mL, each. 600BG4000	\$69.75

Heavy-Duty

50 mL, Pk/12. 600BG100350	\$44.55
100 mL, Pk/12. 600BG1003100	\$75.50
150 mL, Pk/12. 600BG1003150	\$88.50
250 mL, Pk/12. 600BG1003250	\$102.00
400 mL, Pk/12. 600BG1003400	\$116.00
600 mL, Pk/6. 600BG1003600	\$88.75
1000 mL, each. 600BG100310000	\$20.65

Polypropylene, premium, low form, single scale

100 mL, 25 mL graduations, each. 262110000	\$4.05
250 mL, 50 mL graduations, each. 262130000	\$6.75
600 mL, 100 mL graduations, each. 262150000	\$10.50
1000 mL, 200 mL graduations, each. 262160000	\$10.60
2000 mL, 500 mL graduations, each. 262170000	\$32.15
4000 mL, 1000 mL graduations, each. 262190000	\$53.60

Bench Protector



Protect your bench surface using this 40 mil absorbent liner backed with 0.5 mil polyethylene film.

Bench Liner, Roll, 50' x 20", white, Labmat™. **246750000**..... \$95.45

Bottles




Deliver controlled drops from standard 0.030" tip opening with these LDPE bottles. Sealer cap snaps on easily.

Bottles, Dispensing, squeeze bottles

LDPE (low-density polyethylene), Yorker cap

30 mL (1 oz), Bg/12. 116370001	\$33.15
185 mL (6 oz), Bg/12. 116370006	\$63.85

A long nozzle, textured sides, and slim profile make this chemically inert, HDPE container ideal for dispensing immersion oil.



HDPE (high-density polyethylene)

20 mL, each. 18145	\$8.50
---------------------------	--------

Amber glass protects light-sensitive liquids. Screw-cap dropper features a black rubber bulb and glass pipet.



Boston Round

Bottles, Dropping, Amber Glass (flint)

30 mL (1 oz), Bx/12. 3104053	\$19.30
Bottle Dropper Assembly (replacement)	
20/400, 7 x 76 mm Pipet, Pk/12. 3104053D	\$11.10
60 mL (2 oz), Bx/12. 3104054	\$27.00
Bottle Dropper Assembly (replacement)	
20/400, 7 x 89 mm Pipet, Pk/12. 3104054D	\$11.15

French Square



Bottles, Dropping, Clear Glass (flint)

15 mL (1/2 oz), Bx/12. 3104062	\$33.80
Bottle Dropper Assembly (replacement)	
20/400, 7 x 60 mm Pipet, Pk/12. 3104074	\$11.55
30 mL (1 oz), Bx/12. 3104063	\$39.20
Bottle Dropper Assembly (replacement)	
24/400, 7 x 75 mm Pipet, Pk/12. 3104075	\$14.75

Continued

Borosilicate glass bottle;
blue PP cap (GL45)



Bottles, Media, graduated, autoclavable

100 mL, Cs/10. 600BM0100	\$96.85
250 mL, Cs/10. 600BM0250	\$110.00
500 mL, Cs/10. 600BM0500	\$125.00
1000 mL, Cs/10. 600BM1000	\$178.00

Replacement Items

Bottle Caps, GL45, blue, Pk/10. B3000-CAP	\$47.60
Bottle Drip-Free Rings, GL45, Pk/10. B3000-RIN	\$25.65



Globally Harmonized System
(GHS) wash bottles are vented
to prevent solvent drips.
Polyethylene bottles (LDPE)
with polypropylene (PP) caps.

Bottles, Wash, GHS Labeled, 500 mL (empty)

Acetone, Pk/6. 600370051	\$41.30
Bleach, Pk/6. 600370059	\$41.30
Deionized Water, Pk/6. 600370064	\$41.30
Ethanol, Pk/6. 600370052	\$41.30
Methanol, Pk/6. 600370053	\$41.30
Water, Pk/6. 600370061	\$41.30

Bottles, Wash, Safety Labeled, 500 mL (empty)

4-Color diamond, LDPE bottle/PP cap	
Soap, Pk/4. 117160014	\$41.30
Universal (label your own bottle), Pk/4. 117160009	\$41.30



Bottles, Wash, Integral Spout, LDPE bottle/PP cap

250 mL (8 oz), Cs/10. 601633	\$76.50
500 mL (16 oz), Cs/10. 601634	\$87.45



Bottles, Narrow-Mouth, HDPE bottle/polypropylene cap

125 mL, Pk/12. 106200015	\$25.45
250 mL, Pk/12. 106200016	\$36.70
500 mL, 28 mm cap, Pk/12. 106200007	\$42.95
1000 mL, Pk/6. 106200008	\$33.60



Bottles, Wide-Mouth

HDPE bottle/polypropylene cap

125 mL, 38 mm cap, Pk/12. 60033407	\$25.45
250 mL, Pk/12. 106250006	\$36.70
500 mL, 53 mm cap, Pk/12. 60033409	\$42.75
1000 mL, Pk/6. 60033410	\$35.95

Polypropylene bottle and cap, autoclavable

125 mL, 38 mm cap, Pk/12. 60033307	\$25.15
250 mL, 43 mm cap, Pk/12. 60033308	\$36.70
500 mL, 53 mm cap, Pk/12. 106320007	\$42.45
1000 mL, 63 mm cap, Pk/6. 60033310	\$37.40



Bottle Top Dispensers

0.5-5 mL. 600BTDR2	\$312.00
5-60 mL. 600BTDR5	\$372.00

Boxes



Box, Matchbox Style (Pill Box), cardboard, white

2 1/8 x 1 1/8 x 5/8", Pk/72. 20047	\$35.45
---	---------

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. **8342**

Each cardboard box can store approximately 540 paraffin blocks.



246530002

246530001



132040001


Boxes, Glass Disposal, cardboard, polyethylene bag included

Bench Model, 8 x 8 x 10", Pk/6. 246530002	\$122.00
Floor Model, 12 x 12 x 27", Pk/6. 246530001	\$163.00
Glass Disposal Box Cover, reusable, PP. 132040001	\$53.65

*Replacement item for the floor model. Nonrefundable.
Polypropylene cover is autoclavable.*

Continued

Hinged
Polystyrene



Boxes, Plastic, (see following chart)

Plastic Boxes				
Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
H-1	1 x 1 x ¼"	-	Reticles, Pointers	\$0.40
H-1/11B	1 x 1 x ¼"	1	Reticles, Pointers	\$0.50
H-54	2 9/16 x 2 9/16 x ½"	-	Filters	\$1.20
H-54/2916G	2 9/16 x 2 9/16 x ½"	1	Filters	\$1.35
H-28/334B2	3 ¾ x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40
H-10	2 7/8 x 1 3/16 x 3/8"	-	Small Parts	\$0.60
H-21/512B	5 ½ x 1 ¾ x ½"	1	Forceps & Scissors	\$1.80

Brushes



Brushes

Bench, horse hair, 10" long. **BRS-969.10** \$9.95
 Bristle head: 4¾"(L) x 1½"(W); bristle length: 1"
 Flat Bristles (paint style)
 Camel Hair, 1" wide, 6½" total length. **BR665B** (LB-37) \$6.65
 Horsehair, heavy-duty bristles, metal handle. **33761880** \$0.50 ↓
 Round Bristles, camel hair, 8½" long. **BR661** (LB-30) \$11.80
 Swivel-Head, black nylon bristles. **BR565** (LB-14) \$5.00
 Bristle head: 4¼"(L) x 1½"(W); total length: 12"

Burner



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing. Height is 15.6 cm (6⅞")



H-5550 gridded top, 3.8 cm (1½") diameter, produces a uniform blue cone of light. Height is 21.6 cm (8½")

Burners, Bunsen, natural gas, Humboldt Mfg. **USA**

Standard, flame stabilizer. **H-5870** \$32.65
 High Temperature, Meker-type. **H-5550** \$98.80
 Accessory
 Bunsen Burner Natural Gas Connector
 3' Long. **H-2087** \$39.15



Burner Lighter, replaceable single flint, round file. **SOL-290.00** ... \$4.40
 Flints (replacement), Pk/5. **SOL-290.01** \$2.95

Carboys



Carboys with Spigot, polypropylene, round

10 L (2.5 gal). **60034051** \$84.45
 20 L (5 gal). **60034052** \$117.00

Centrifuges

GCH-24



Centrifuge, Hematocrit, 24-place. GCH-24..... \$2,439.00
24-place capillary rotor and reading disk lid are included.



C1012



C1008-C

Centrifuges, Mini

5,500 RPM (2000 x g), 12-place, MyFuge™12. **C1012** \$448.00
(1-year warranty)

6,000 RPM (2000 x g), 8-place, MyFuge™Mini. **C1008-__** * \$333.00

*Specify clear (-C), green (-G), or purple (-P) for lid color.
(1-year warranty)

Clamps & Supports

Rectangular
Steel Base



Triangular
Cast-Iron Base



Support (Ring) Stands, finished metal base/plated support rod

Rectangular Steel Base

6 x 4" base, 5/16 x 18" rod. **3104SSS66** \$15.45

8 x 5" base, 3/8 x 20" rod. **3104SSS67** \$19.70

9 x 6" base, 1/2 x 24" rod. **3104SSS68** \$25.00

Rectangular Cast-Iron Base

Provides greater stability than steel base.

6 x 4" base, 5/16 x 18" rod. **3104SSC14** \$18.80


8 x 5" base, 3/8 x 20" rod. **3104SSC15** \$20.80

9 x 6" base, 1/2 x 24" rod. **3104SSC16** \$30.40

Triangular Cast-Iron Base

5" Legs, 1/2 x 24" rod. **1257G2324** \$25.00

6" Legs, 1/2 x 36" rod. **1257G23** \$30.65



600COBR3FR




600THCL01



600CLBURDD



97-4023



1257G49

Support Clamps

Burette, plastic coated jaws, 360° rotatable. **600COBR3FR** \$8.70

Burette, Double, deluxe. **600CLBURDD** \$22.95

Plastic covered spring jaws on a metal frame securely grip
burettes up to 30 mm OD. Fits rods up to 3/4" (19 mm) OD.

Hook Collar, nickel plated. **1257G49** \$6.45

Fits rods up to 1/2" (13 mm) in diameter.

Right-Angle, aluminum. **97-4023** \$8.60

Fits rods up to 5/8" in diameter.

Thermometer, aluminum, rotates 360°. **600THCL01** \$11.10

Fits rods up to 1/2" (13 mm) in diameter.



Support Rings, plated cast-iron

3" OD (2 1/2" ID). **600SRC103** \$9.55

4" OD (3 1/2" ID). **600SRC104** \$12.45

5" OD (4 1/2" ID). **600SRC105** \$13.25



Plain



Ceramic
Center

Wire Gauze Squares, galvanized iron wire w/turned edges

Plain

4 x 4". **31040004** \$1.45

5 x 5". **31040005** \$1.30

6 x 6". **31040006** \$1.55


Ceramic Center

4 x 4". **31040001** \$1.60

5 x 5". **31040002** \$1.70


6 x 6". **31040003** \$2.25

Continued



Wire Triangles, clay pipe, galvanized iron wire

1½" Pipe Length. 600CPT015	\$1.80
2" Pipe Length. 600CPT020	\$2.10
2½" Pipe Length. 600CPT025	\$2.30



Test Tube Clamp, brass wire, 13 cm (5"). **600TTCL02** \$2.50



Hoffman



Mohr

Tubing Clamps, nickel plated

Hoffman, swivel lower jaw, Pk/3. 97-1801	\$7.80
Mohr Pinchcock, each. 3104MC42	\$2.40

Counters



Counters

Colony, handheld electronic. 345120000	\$317.00
---	----------



Handheld Tally, 0-9999 range. 600A	\$9.45
---	--------

Crucibles



Crucibles, Tall Form with Lid, porcelain

5 mL Capacity, Pk/6. 600JCT005	\$26.55
10 mL Capacity, Pk/6. 600JCT010	\$30.35
15 mL Capacity, Pk/6. 600JCT015	\$34.05
30 mL Capacity, Pk/6. 600JCT030	\$42.45
50 mL Capacity, Pk/6. 600JCT050	\$48.90
100 mL Capacity, Pk/6. 600JCT100	\$68.75


Cylinders

Cylinders, graduated

Glass (borosilicate),

Class B (economy), Single Scale

10 mL, 0.2 mL grad., ea. 310CSS10	\$3.50
25 mL, 0.5 mL grad., ea. 310CSS25	\$4.00
50 mL, 1.0 mL grad., ea. 310CSS50	\$4.70
100 mL, 1.0 mL grad., ea. 310CSS100	\$5.90
250 mL, 2.0 mL grad., ea. 310CSS250	\$9.25
500 mL, 5.0 mL grad., ea. 310CSS500	\$14.35
1000 mL, 10.0 mL grad., ea. 310CSS1000	\$23.60




Class A, individually certified, double metric scale

50 mL, 1.0 mL grad., ea. 600CY302050	\$43.55
100 mL, 1.0 mL grad., ea. 600CY3020100	\$51.85
250 mL, 2.0 mL grad., ea. 600CY3020250	\$77.00
500 mL, 5.0 mL grad., ea. 600CY3020500	\$100.00
1000 mL, 10.0 mL grad., ea. 600CY30201000	\$190.00

Polypropylene, single scale

10 mL, 0.1 mL grad., ea. 284520000	\$9.70
50 mL, 1.0 mL grad., ea. 284540000	\$12.90
100 mL, 1.0 mL grad., ea. 601079	\$15.70
250 mL, 2.0 mL grad., ea. 284560000	\$17.30
500 mL, 5.0 mL grad., ea. 284570000	\$23.80
1000 mL, 10.0 mL grad., ea. 601082	\$35.60
2000 mL, 20.0 mL grad., ea. 601094	\$65.90

Desiccators



Non-Vacuum



Vacuum

Glazed porcelain plate has 5 mm perforations with a 25 mm center hole.

Desiccators

Non-Vacuum, borosilicate glass

8" (20 cm) ID. 600DSGL210	\$155.00
--	----------

Vacuum, polycarbonate top/polypropylene bottom

6" (15 cm), 13.8 cm plate diameter. 60055204	\$81.15
8" (20 cm), 19 cm plate diameter. 60055205	\$102.00
10" (25 cm), 23.8 cm plate diameter. 60055206	\$130.00

Drying Racks



189330023



Drying Racks, Benchtop, Bel-Art

Single-sided. **188050000**.....\$174.00

Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 43 x 17 x 48 cm(H).

Double-sided

2-Tier. **189330022**.....\$433.00

3-Tier. **189330023**.....\$589.00

Drying Rack, Benchtop/Wallmount, Heathrow

Single-sided. **HS23243A**.....\$148.00

Dry up to 90 pieces of labware with this self-standing or wall-mounted rack. Constructed of vinyl-coated wire hairpin loops with a polystyrene tray and rack. Dimensions: 46.2 x 18.2 x 52.5 cm(H).

Electrophoresis



Electrophoresis System, MyGel™ Mini. **E1101**.....\$649.00

Eyewear

Absorbs
99.9% UV



Glasses, Safety, integral side shield. **43574002**.....\$6.00



Flexible frame and adjustable elastic strap provide a tight seal.

Goggles, Safety, vented (chemical splash). **97-2001**.....\$4.40

Filter & Blotting Papers



Ahlstrom - Munksjo
Filter Paper

Qualitative, Standard Grades (for gravity filtration)

Slow flow (4 mL/min), 1.5 µm retention, Grade 610

(Whatman 5; S&S 602)

9 cm, Pk/100. **6100-0900**.....\$12.90

12.5 cm, Pk/100. **6100-1250**.....\$16.10

15 cm, Pk/100. **6100-1500**.....\$19.40

20.5 cm, Pk/100. **6100-2050**.....\$30.65

24 cm, Pk/100. **6100-2400**.....\$42.70

Medium flow (30 mL/min), 2 µm retention, Grade 642

(Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)

9 cm, Pk/100. **6420-0900**.....\$10.05

15 cm, Pk/100. **6420-1500**.....\$18.05

Medium flow (55 mL/min), 2.5 µm retention, Grade 601

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

4.25 cm, Pk/100. **6010-0425**.....\$7.70

9 cm, Pk/100. **6010-0900**.....\$10.35

12.5 cm, Pk/100. **6010-1250**.....\$14.85

15 cm, Pk/100. **6010-1500**.....\$18.10

18.5 cm, Pk/100. **6010-1850**.....\$26.25

20.5 cm, Pk/100. **6010-2050**.....\$28.30

33 cm, Pk/100. **6010-3300**.....\$76.40

Medium flow (60 mL/min), 6 µm retention, Grade 613

(Whatman Student; S&S 595, 0903)

9 cm, Pk/100. **6130-0900**.....\$8.80

15 cm, Pk/100. **6130-1500**.....\$15.35

20.5 cm, Pk/100. **6130-2050**.....\$23.00

33 cm, Pk/100. **6130-3300**.....\$51.15

Medium flow (200 mL/min), 10 µm retention, Grade 631

(Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)

20.5 cm, Pk/100. **6310-2050**.....\$28.50

33 cm, Pk/100. **6310-3300**.....\$77.20

Fast flow (235 mL/min), 25 µm retention, creped, Grade 615

9 cm, Pk/100. **6150-0900**.....\$9.10

12.5 cm, Pk/100. **6150-1250**.....\$13.20

15 cm, Pk/100. **6150-1500**.....\$15.85

20.5 cm, Pk/100. **6150-2050**.....\$23.70

33 cm, Pk/100. **6150-3300**.....\$52.70

Very Fast flow (360 mL/min), 35 µm retention, creped, Grade 617

9 cm, Pk/50. **6170-0900**.....\$8.30

12.5 cm, Pk/50. **6170-1250**.....\$11.05

15 cm, Pk/50. **6170-1500**.....\$12.85

20.5 cm, Pk/50. **6170-2050**.....\$19.60

33 cm, Pk/50. **6170-3300**.....\$42.00

Qualitative, Pleated (Fluted)

Slow flow (30 mL/min), 2 µm retention, Grade 542

(Whatman 2V)

12.5 cm, Pk/100. **5420-1250**.....\$82.45

18.5 cm, Pk/100. **5420-1850**.....\$93.30

Qualitative, Wet-Strengthened Grades (for vacuum filtration)

Medium flow (65 mL/min), 3 µm retention, **Grade 909**

(Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)

32 cm, Pk/100. **9090-3200**\$83.20

Fast flow (245 mL/min), 27 µm retention, creped, **Grade 962**

(Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)

11 cm, Pk/100. **9620-1100**\$12.35

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

Slow flow (4 mL/min), 1.5 µm retention, **Grade 94**

(Whatman 42; S&S 589/3)

9 cm, Pk/100. **0940-0900**\$34.85

15 cm, Pk/100. **0940-1500**\$60.00

Medium flow (30 mL/min), 2 µm retention, **Grade 74**

(Whatman 40; S&S 589/2, 589/4)

9 cm, Pk/100. **0740-0900**\$34.85

15 cm, Pk/100. **0740-1500**\$60.00

Fast flow (220 mL/min), 15 µm retention, **Grade 55**

(Whatman 54, 541; S&S 1505)

9 cm, Pk/100. **0550-0900**\$40.00

15 cm, Pk/100. **0550-1500**\$63.25

Blotting/Chromatography/Electrophoresis Papers

0.19 mm thick, 125 mm/30 min absorption, **Grade 601**

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

46 x 57 cm, Pk/100. **6018-4657**\$157.00

0.34 mm thick, 140 mm/30 min absorption, **Grade 238**

(Whatman 3MM, S&S 598, 231^{cr})

20 x 20 cm, Pk/100. **2388-2020**\$80.00

46 x 57 cm, Pk/100. **2388-4657**\$305.00

0.42 mm thick, 150 mm/30 min absorption, **Grade 237**

(Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703;

Albet 2040B, FP591; Munktel 3S/H, FN7a)

46 x 57 cm, Pk/100. **2378-4657**\$352.00

0.95 mm thick, 180 mm/30 min absorption, **Grade 243**

(Whatman 17CHR; S&S 2668)

11 x 14 cm, Pk/100. **2438-1114**\$31.75

20 x 20 cm, Pk/100. **2438-2020**\$66.15

2.48 mm thick, 240 mm/30 min absorption, **Grade 320**

(S&S 2727)

11 x 14 cm, Pk/50. **3208-1114**\$63.45

20 x 20 cm, Pk/50. **3208-2020**\$93.05



Smooth tips &
polypropylene grips

Forceps, Filter, 11.5 cm. 300-025\$4.40

Freezer Racks



Freezer Racks, stainless steel, USA

175V-13-2L

Restocking fee of 20% on all returned freezer racks listed below.

No returns on custom racks.

Horizontal (fits upright freezers)

Horizontal, Holds 2" Boxes

	1-9	≥10
9-Place (3L x 3H). 175HS-9-2	\$59.00	53.00
Dimensions: 16.31(L) x 5.5(W) x 6.5" (H)		
12-Place (3L x 4H). 175HS-12-2	\$78.00	71.00
Dimensions: 16.31(L) x 5.5(W) x 8.75" (H)		
15-Place (3L x 5H). 175HS-15-2	\$98.00	89.00
Dimensions: 16.31(L) x 5.5(W) x 11" (H)		
16-Place (4L x 4H). 175HL-16-2	\$105.00	95.00
Dimensions: 22(L) x 5.5(W) x 8.75" (H)		
20-Place (4L x 5H). 175HL-20-2	\$132.00	119.00
Dimensions: 22(L) x 5.5(W) x 11" (H)		

Horizontal, Holds 3" Boxes

9-Place (3L x 3H). 175HS-9-3	\$61.00	55.00
Dimensions: 16.31(L) x 5.5(W) x 9.5"(H)		
12-Place (4L x 3H). 175HL-12-3	\$80.00	73.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		



Horizontal Drawer (fits upright freezers)

Horizontal Drawer, Holds 2" Boxes

	1-9	≥10
9-Place (3L x 3H). 175HD-332	\$85.00	77.00
Dimensions: 16.5(L) x 5.5(W) x 7.125"(H)		
12-Place (3L x 4H). 175HD-342	\$111.00	101.00
Dimensions: 16.5(L) x 5.5(W) x 9.5" (H)		
15-Place (3L x 5H). 175HD-352	\$141.00	128.00
Dimensions: 16.5(L) x 5.5(W) x 11.812" (H)		
16-Place (4L x 4H). 175HD-442	\$151.00	137.00
Dimensions: 22(L) x 5.5(W) x 9.5" (H)		
20-Place (5L x 4H). 175HD-542	\$188.00	171.00
Dimensions: 27.25(L) x 5.5(W) x 9.5" (H)		

Horizontal Drawer, Holds 5" Boxes

	1-9	≥10
6-Place (3L x 2H). 175RD-325	\$74.00	67.00
Dimensions: 18.125(L) x 6(W) x 10.75" (H)		
8-Place (4L x 2H). 175RD-425	\$99.00	90.00
Dimensions: 24(L) x 6(W) x 10.75" (H)		

Horizontal Drawer, Holds 15 mL tubes

160-Place (80L x 2H). 175HD-32ML-15	\$99.00	90.00
Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)		
220-Place (110L x 2H). 175HD-42ML-15 ...	\$129.00	118.00
Dimensions: 22(L) x 5.5(W) x 9.5" (H)		

Horizontal Drawer, Holds 50 mL tubes

84-Place (42L x 2H). 175HD-42ML-50	\$129.00	118.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		

Continued

Horizontal Drawer-Tray (fits upright freezers)

Drawer-Tray, Holds 2" Boxes, aluminum trays	1-9	≥10
12-Place (3L x 4H). 175HDT-342	\$96.00	86.00
Dimensions: 16.63(L) x 5.5(W) x 8.5" (H)		
18-Place (3L x 6H). 175HDT-362	\$143.00	129.00
Dimensions: 16.63(L) x 5.5(W) x 12.75" (H)		
24-Place (4L x 6H). 175HDT-462	\$189.00	172.00
Dimensions: 22(L) x 5.5(W) x 12.75" (H)		


Vertical (fits chest freezers)

Vertical, Holds 2" Boxes

13-Place. 175V-13-2L	\$88.00	79.00
Dimensions: 5.5(L) x 5.625(W) x 28.562" (H)		

Tube Storage Boxes (cryogenic) – see in GENERAL LAB SUPPLIES & EQUIPMENT section

Funnels



Funnels

Disposable, polypropylene, short stems, Simport		
★ 55 mm Top ID, takes 11 cm filter, Cs/100. F490-3	\$40.70	
★ 65 mm Top ID, takes 12.5 cm filter, Cs/100. F490-4	\$40.70	
Glass (borosilicate), long stem (150 mm)		
50 mm Top ID, Pk/6. 600GF614050	\$49.95	
65 mm Top ID, Pk/6. 600GF614065	\$60.00	
75 mm Top ID, Pk/6. 600GF614075	\$62.00	
Heavy-Duty, Molded-In Handle, polypropylene		
250 mm Top ID, 2.5 L capacity, each. 60057232	\$26.55	
Micro, polypropylene, Bel-Art		
24 mm Top ID, 2.6 mL capacity, Pk/12. 1146850024	\$42.30	
28 mm Top ID, 5.2 mL capacity, Pk/12. 146850028	\$34.65	
43 mm Top ID, 15.7 mL capacity, Pk/12. 146850043	\$72.95	
Powder, polypropylene		
Parallel Stem		
78 mm Top ID, Pk/12. 60057222	\$16.85	
100 mm Top ID, Pk/6. 146600100	\$16.35	
145 mm Top ID, Pk/6. 60057224	\$21.50	

Funnels (cont.)

Standard Stem, Premium, polypropylene, Bel-Art	
Exclusive filter-step internal design provides speedier filtration; external ribs prevent air-lock.	
55 mm Top ID, takes 9 cm filter, Pk/12. 146760055	\$47.75
65 mm Top ID, takes 11 cm filter, Pk/12. 146760065	\$49.60
75 mm Top ID, takes 12.5 cm filter, Pk/6. 146760075	\$25.65
100 mm Top ID, takes 18.5 cm filter, Pk/4. 146760100	\$32.60
Standard Stem, Economy, polypropylene	
50 mm Top ID, each. 310FUNPP50	\$0.80
62 mm Top ID, each. 310FUNPP62	\$0.95
75 mm Top ID, each. 310FUNPP75	\$1.15
100 mm Top ID, each. 310FUNPP100	\$1.65

Gauze Pads

These 4 x 4" nonsterile gauze sponges are 12-ply.



Gauze Pads , 10 pk/200, Cs/2000. 44122000	\$88.45
---	---------

Germicidal Wipes

Sani-Cloth Plus



Germicidal Wipes , Tub/160 440Q89072	\$10.00
--	---------

Gloves & Hand Protection



Nitrile Gloves

Gloves

Purple Nitrile, powder-free, 5.5 mil, textured fingers, Halyard	
Nitrile gloves are hypoallergenic and provide greater puncture and chemical resistance than latex or vinyl gloves.	
Small, bx/100, Cs/10 Bx. 55081	\$139.00
Medium, bx/100, Cs/10 Bx. 55082	\$139.00
Large, bx/100, Cs/10 Bx. 55083	\$139.00
X-Large, bx/90, Cs/10 Bx. 55084	\$139.00
Purple Nitrile-XTRA, 12" Long, powder-free, 5.9 mil, Halyard	
Offers additional protection against chemotherapy drugs and common chemicals.	
Small, bx/50, Cs/10 Bx. 50601	\$213.00
Medium, bx/50, Cs/10 Bx. 50602	\$213.00
Large, bx/50, Cs/10 Bx. 50603	\$213.00
X-Large, bx/50, Cs/10 Bx. 50604	\$213.00


Continued



Safely grip hot or cold labware with this silicone rubber hand protector partially lined with mini suction cups. Withstands temperatures from 260°C (500°F) to -57°C (-70°F).

Hand Protector, red. 600HNDPTC \$20.40

Homogenizers



Rapidly homogenize, emulsify, suspend and disrupt 0.1 mL to 250 mL samples. Includes a 5 x 50 mm flat-bottom probe & 7 x 50 mm saw tooth probe for use with micro tubes.

Homogenizer, Handheld. D1000..... \$1,269.00



Benchmark BeadBlaster™ 24

Homogenize a variety of sample types, from soft tissue to bone, in 35 seconds or less.

Homogenizer, Benchtop, Microtube. D2400..... \$8,404.00

- Holds (24) 2.0 mL tubes, or (6) 5 mL tubes w/optional tube holder
- Speed adjusts from 4.00 to 7.00 m/s in 0.05 m/s increments
- Stores 50 programs in memory, each with a maximum of 10 cycles
- 2-year warranty

Homogenizer Prefilled Tubes

2.0 mL screw-capped tubes are prefilled with your choice of the following beads:

Glass, acid washed

0.1 mm, Pk/50. **D1031-01**.....\$135.00

0.5 mm, Pk/50. **D1031-05**.....\$135.00

Zirconium Beads, acid washed/heat treated/nuclease-free

1.5 mm, Pk/50. **D1032-15**.....\$189.00

3.0 mm, Pk/50. **D1032-30**.....\$189.00

Stainless Steel Beads, acid washed

2.8 mm, Pk/50. **D1033-28**.....\$189.00

Hot Plates, Hot Plate Stirrer & Magnetic Stirrers

310HOTPLT



H4000-H

Hot Plates, analog

6" Diameter Cast Iron Burner, 1000W. **310HOTPLT** \$32.40

7.5 x 7.5" Ceramic Platform, 80-380°C. **H4000-H** \$312.00

H4000-HS



Hot Plate Stirrer, 80-380°C, Benchmark. H4000-HS \$449.00



S1005

H4000-S



Magnetic Stirrers, analog, Benchmark

3.4 x 3.4" Platform, 100-1500 rpm, MiniMag™. **S1005**..... \$148.00
(1-year warranty)

7.1 x 7.1" Ceramic Platform, 60-1500 rpm. **H4000-S**..... \$312.00
(2-year warranty)

Incubators



Incubators, Mini Digital, MyTemp™, Benchmark (2-year warranty)

Heat Only, 1° above ambient to 60°C. **H2200-H**\$609.00

Heat/Cool, 15° below ambient to 60°C. **H2200-HC**\$917.00

Heat/Cool/CO₂, 10° below ambient to 55°C. **H2300-HC2** \$4,169.00

- 20 L capacity (0.8 cu ft)
- Digital temperature control with LED display
- Accepts bottles and flasks up to 2 L
- Two adjustable/removable shelves
- 9.3"(D) x 10.3"(W) x 12.8"(H) interior;
- 14.5"(D) x 13.2"(W) x 18.7"(H) exterior



Incubators, Dual-Convection, SureTemp™ (2-year warranty)

40 L. **H2505-40**\$2,139.00

13.7"(D) x 13.7"(W) x 13.7"(H) interior

70 L. **H2505-70**\$2,699.00

13.7"(D) x 17.7"(W) x 17.7"(H) interior

130 L. **H2505-130**\$3,449.00

17.7"(D) x 21.7"(W) x 21.7"(H) interior

- Precisely control temperature up to 75°C
- Record temperatures with SureCheck™ data logging
- Select either gravity or mechanical convection mode



Incubator, Shaking, Incu-Shaker™ Mini. **H1001-M**\$3,199.00

9.5" x 11.5" platform w/ rubber mat

- 5° above ambient to 70°C
- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Domed lid accepts flasks up to 1 L
- 16.2"(D) x 11"(W) x 13"(H) exterior

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

Steel Platform, Magic Clamp™, 9.5 x 11". **H1000-MR**\$326.00

Clamps, magnetic

250 mL Flask. **H1000-MR-250**\$49.85

500 mL Flask. **H1000-MR-500**\$49.85

1000 mL Flask. **H1000-MR-1000**\$76.25

Mortar and Pestle Sets

Mortar and Pestle Sets

Glass, **USA**

59 mL (2 oz) Capacity. **826-2**\$16.35

118 mL (4 oz) Capacity. **826-3**\$17.20

237 mL (8 oz) Capacity. **826-4**\$20.40

Porcelain, outside wall and rim glazed

30 mL Capacity, 63 mm top OD. **LPC-154**\$5.00

70 mL Capacity. **600JMD070**\$13.35

150 mL Capacity. **600JMD150**\$17.60

320 mL Capacity, 135 mm top OD. **LPC-156**\$11.50



Parafilm



Parafilm, thermoplastic sealing film, **USA**

2" x 250' (5 cm x 76 m) Roll. **HS234526A**\$34.00

4" x 125' (10 cm x 38 m) Roll. **HS234526B**\$34.90

Petri Dishes



Petri Dishes, polystyrene, Kord-Valmark

60 x 15 mm, stackable

Cs/500, sterile sleeve/20. **2901**\$118.00

100 x 15 mm

Slippable (rimless lid for automation)

Standard, Cs/500, sterile sleeve/25. **2917**\$126.00



Stackable

Standard, Cs/500, sterile sleeve/25. **2900**\$126.00

100 x 25 mm, stackable

Cs/340, sterile sleeve/17. **2906**\$137.00

pH Paper



Measure pH with 1 pH unit accuracy from 1-14.

Package contains 100 paper strips, 2 x 1/4", in a plastic vial with color chart.

pH Paper, pH 1-14, Pk/100. **600PHP100**\$4.00

Pipetting



Dropper Bulbs

Amber Latex

2 mL Capacity, 52 mm long, Pk/12. **20621B** \$7.35

2 mL Capacity, 52 mm long, Pk/72. **HS20622B** \$42.55

Natural Rubber

1 mL Capacity, 30 mm long, black, Pk/12. **200-1B** \$4.00


2 mL Capacity, 32 mm long, red, Pk/12. **200-2B** \$6.40

Natural Rubber w/Flange, blue

1 mL Capacity, 29 mm long, Pk/24. **200-1C** \$12.75



Natural Rubber w/Flange



glass w/rubber bulb; 75 mm (3")

Medicine (Eye) Dropper, Pk/12. 97-4301 \$4.00

3-Valve



Pipette Pump



Pipette Pump II

Pipet Fillers, manual

3-Valve, rubber, black. **207B** \$17.60

Wheel-Type, Pipette Pump™

Fits Pipets Up To 2 mL, blue. **378970000** \$38.75

Fits Pipets Up To 10 mL, green. **378980000** \$41.50

Fits Pipets Up To 25 mL, red. **378990000** \$55.10

Wheel-Type, Pipette Pump™ II, quick-release

Fits Pipets Up To 2 mL, blue. **379111002** \$48.35

Fits Pipets Up To 10 mL, green. **379111010** \$48.35

Fits Pipets Up To 25 mL, red. **379111025** \$63.85



Package includes: Diamond SeroLogic portable pipette controller; two rechargeable batteries; battery charger; three PTFE membrane filters – (1) 0.20 µm & (2) 0.45 µm (one installed); desktop stand and wall mounting bracket.

Pipet Filler, Motorized, rechargeable

Diamond SeroLogic (1-year warranty). **3380** \$369.00



Dispense fluids accurately and safely with these color-coded and plugged pipets. Negative graduations provide extra volume (not available on 25 mL). RNase-/DNase-/Pyrogen-free, and noncytotoxic. Calibrated "To Deliver" (TD).

Pipets, Serological, polystyrene, UniPlast™

Nonsterile

5 mL, blue band, pk/25, Cs/375. **1755** \$123.00

10 mL, orange band, pk/25, Cs/250. **1775** \$88.50

25 mL, red band, pk/10, Cs/100. **1784** \$82.50

Sterile Pk/1 (individually paper/plastic wrapped)

1 mL, yellow band, Cs/500. **1700** \$123.00

2 mL, green band, Cs/500. **1720** \$132.00

5 mL, blue band, reverse grad., Cs/250. **1740** \$98.60

10 mL, orange band, reverse grad., Cs/200. **1760** \$89.50

25 mL, red band, reverse grad., Cs/100. **1780** \$94.75


50 mL, black band, reverse grad., Cs/100. **1820** \$209.00

Sterile Bg/25


1 mL, yellow band, Cs/1000. **1710** \$202.00

10 mL, orange band, Cs/250. **1770** \$100.00

25 mL, red band, Cs/100. **1782** \$85.75



5 mL Capacity,
1 mL graduated



7 mL Capacity,
3 mL graduated

Pipets, Transfer, LDPE (low-density polyethylene)

Graduated, Nonsterile, 15 cm long

5 mL capacity, 1 mL graduated, 25 drops/mL



Bx/500. **PTP-02** \$23.10

Cs/10 Bx. **PTP-02** \$200.00

5 mL capacity, thick-walled

Cs/10 Bx. **188BP1200** \$164.00

7 mL capacity, 3 mL graduated, 14 drops/mL

Bx/500. **PTP-01** \$23.10

Cs/10 Bx. **PTP-01** \$200.00

Graduated, Sterile Pk/1 (individually wrapped), 15 cm long

5 mL capacity, 1 mL graduated

Bx/400. **P200-521S** \$49.85

Graduated, Sterile Pk/20, 15 cm long

5 mL capacity, 1 mL graduated

Bx/500. **HS206371B** \$49.00

NEW



These mini, air displacement pipettors provide an economical and accurate method for pipetting fixed volumes. Pipettors feature a two-step plunger and are fully autoclavable.

Pipettors, Single Channel

Fixed Volume, pocket-size, Diamond® Jr.

5 µL, red. 3322-5	\$32.95
10 µL, pink. 3322-10	\$32.95
20 µL, green. 3322-20	\$32.95
50 µL, purple. 3322-50	\$32.95
100 µL, orange. 3322-100	\$32.95

Single channel, adjustable volume pipettors are fully autoclavable.



Adjustable Volume, Diamond® Apex™ (3-year warranty)

0.5-10 µL, 0.02 µL increments, red. 3351-10	\$219.00
5-50 µL, 0.1 µL increments, yellow. 3351-50	\$219.00
10-100 µL, 0.1 µL increments, blue. 3351-100	\$219.00
100-1000 µL, 1 µL increments, orange. 3351-1000	\$219.00

Pitchers



Polypropylene



Stainless Steel

Pitchers

Polypropylene, graduated (mL), tall form, tapered

1000 mL, 50 mL graduations, each. 289910000	\$7.20
2000 mL, 100 mL graduations, each. 289920000	\$8.35
3000 mL, 100 mL graduations, each. 289930000	\$13.80

Stainless Steel, graduated (mL/oz), tapered, Polar Ware

1000 mL/32 oz, each. T1063	\$73.20
---	---------

Rockers, Shakers & Mixers



BR1000



B3D1020



B3D1320

Rockers, Fixed Speed, Benchmark (2-year warranty)

2-D (seesaw), 12 x 12" platform w/flat mat

BenchBlotter™. **BR1000** \$557.00

- Fixed speed of 12 rocks/minute
- Platform angle is $\pm 8^\circ$
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)

3-D (nutating)

Small, 10.5 x 7.5" platform w/dimpled mat

MiniMixer™. **B3D1020** \$459.00

Large, 12 x 12" platform w/dimpled mat

BioMixer™. **B3D1320** \$566.00

- Fixed speed of 24 rocks/minute
- Platform angle is $\pm 20^\circ$
- Cold room & incubator safe (4-45°C)
- Load capacity is 1.8 kg (4 lbs)

Continued



Rockers, Variable Speed

Benchmark (2-year warranty)

2-D (seesaw), 14 x 12" platform w/flat mat

BenchRocker™ 2D. **BR2000**\$849.00

Stacking Platform. **BR2000-STACK**\$167.00

14 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP**\$60.00

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 3.6 kg (8 lbs)

3-D (nutating), 12 x 12" platform w/dimpled and flat mats

BenchRocker™ 3D. **B3D2300**\$1,057.00


Stacking Platform. **BR1000-STACK**\$154.00

12 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP**\$60.00

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)



55S



55D



Rockers, Variable Speed

Reliable Scientific (5-year warranty), **USA**

2-D (seesaw)

Standard, 11 x 14" platform w/rubber mat

Single Platform. **55S1114**\$603.00

Double Platform, 3¼" apart. **55D1114**\$665.00

Large, 12 x 16" platform w/rubber mat

Single Platform. **55S1216**\$647.00

Double Platform, 3¼" apart. **55D1216**\$704.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



Rocker, Tube, 12 x 4.5" grooved mat, Benchmark. **M2100**\$416.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is ±30°
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs)
- 2-year warranty

Orbital Shakers



BT30

Aerate and mix with these low to high capacity orbital shakers.



BT302

Orbital Shakers, Benchmark (2-year warranty)

Low Speed, 14 x 12" platform w/flat rubber mat

Orbi-Blotter™. **BT30**\$957.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

Small, 9.5 x 8" platform w/flat rubber mat

Orbi-Shaker™ JR. **BT302**\$1,177.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min



BT3001

Large, 13 x 12" platform w/flat rubber mat

Orbi-Shaker™. **BT3001**\$1,769.00

- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L, or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR** ..\$354.00


Clamps, magnetic

250 mL Flask. **H1000-MR-250**\$49.85

500 mL Flask. **H1000-MR-500**\$49.85

1000 mL Flask. **H1000-MR-1000**\$76.25

Rotating Mixers



Tube Rotator, programmable, Benchmark (2-year warranty)

Roto-Bot Mini. **R4040** \$766.00
 Roto-Bot. **R4045** \$1,115.00

- Speed adjusts from 1-80 rpm
- Tube capacity (starter set)
 - R4040: 12 x 1.5 mL, 2 x 15 mL & 2 x 50 mL
 - R4045: 30 x 1.5 mL, 8 x 15 mL & 2 x 50 mL
- Timer is programmable from 1 sec – 99 hrs, 59 min / Cont.


Vortexers



Vortex Mixers, Benchmark (2-year warranty)

Vornado™, mini vortexer, fixed 2800 rpm. **BV101** * \$233.00
 *Specify Gray, Blue or Purple cup head
 BenchMixer™, single cup head, 200-3200 rpm. **BV1000** \$302.00
 Flat Head Attachment, 3". **BV1000-FLAT** \$46.55

Rulers, Meter Stick & Tape Measure



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible
 Transparent. **4461** \$0.95
 White. **4460** \$0.55
 Stainless Steel, 12 mm wide, 0.5 mm thick, **USA. 605-065** \$6.70

Meter Stick, lacquered wood, plain ends. **605-050** \$4.90

One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

Tape Measure, 152.4 cm/60", reinforced fiberglass. **11036** \$3.70

One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.

Sharps Containers



Sharps Containers, red polypropylene

5.4 Quart, horizontal entry lid, Prevent™, each. **4402269** \$9.90
 Dimensions: 12" x 4.75" x 11"(H)
 12 Gallon, hinged vertical entry lid, each. **440321391** \$32.20
 Dimensions: 18.25" x 12.75" x 18.75"(H)

Specimen Containers

Screw-Cap, 90 mL



Snap-Cap



Specimen Containers

Screw-Cap, polypropylene (PP)
 90 mL, sterile, Cs/400. **840175** \$150.00
 480 mL (16 oz), Cs/24. **59014510** \$71.95
 83 mm (3 1/4") opening
 960 mL (32 oz), Cs/12. **59014511** \$59.05
 111 mm (4 3/8") opening
 Snap-Cap, white HDPE bottom/white LDPE lid
 480 mL (16 oz), Cs/100. **271016W** \$95.75
 9.9 cm opening x 11.1 cm high
 1000 mL (34 oz), Cs/100. **271032W** \$144.00
 11.7 cm opening x 13.7 cm high

Spot Plate



Spot Plate, glazed white porcelain, 21 mm diameter wells

12-Well, 115 x 95 x 7 mm, each. **3104187701** \$4.50

Stirring Bars & Rods



Stirring Bars, Teflon®

Centered Sphere, Saturn Spinbar®

40 mm (1½") Long. 371160002	\$31.30
50 mm (2") Long. 371160003	\$60.25

Micro, Spinbar®

White

5 x 2 mm. 371190005	\$6.55
7 x 2 mm. 371190007	\$6.55
10 x 3 mm. 371190010	\$7.30
12.7 x 3 mm. 371190127	\$7.50
15 x 1.5 mm. 371190015	\$8.10

Red

7 x 2 mm. 371210013	\$6.80
8 x 1.5 mm. 371210001	\$6.05

Octagonal, molded-on pivot ring

White

12.7 x 8 mm (½ x 5/16"). 371101256	\$6.15
15.9 x 8 mm (5/8 x 5/16"). 371100058	\$6.15
25 x 8 mm (371100138). 600MSB25	\$7.65
40 x 8 mm (371101128). 600MSB40	\$11.70
48.5 x 10 mm (371100238). 600MSB48	\$15.15

Red, 12.7 x 3.2 mm (½ x 1/8"), no pivot ring. **371090034**.....

Stirring Bar Assortment, octagonal, Pk/6. **371600000**.....

Includes two of each size: 25.4 x 8 mm (1 x 5/16"),
38.1 x 8 mm (1½ x 5/16"), and 50.8 x 8 mm (2 x 5/16").

Polygonal, molded-on pivot ring, Spinbar®

12 x 4.5 mm (7/16 x 3/16"). 371220012	\$6.50
--	--------

Triangular, vanes w/cross-mounted magnets, Spinvane®

Fits 1 mL Vials. 371340000	\$12.85
Fits 3/5 mL Vials. 371350000	\$12.85



Stirring Bar Retriever, magnetic

45.7 cm (18") Long. 377730000	\$25.85
--	---------



Stirring Rods, glass, annealed ends

150 x 5 mm (6 x 3/16"), Pk/12. 600GSR006	\$5.15
Rubber Policeman, fits 5 mm rod, Pk/12. 201BB	\$6.65
254 x 6 mm (10 x 1/4"), Pk/12. 600GSR010	\$7.30

Surgical Drape

Non-fenestrated towels
are sterile, latex-free
and poly-lined.



Surgical Drape, 18 x 26", Bx/50. 44025-515	\$35.30
---	---------

Tape

Kraft paper with
rubber resin adhesive.



Tape, Label, Write-On™, Bel-Art

1" Core, 15 yd (540") roll

¾" Wide

Green, Pk/4. 134820075	\$34.35
White, Pk/4. 134800075	\$32.50

3" Core, 40 yd roll


½" Wide

White, each. 134630005	\$13.55
White, Pk/4. 134830075	\$71.75

¾" Wide

Green, each. 134633075	\$17.95
White, Pk/4. 134830075	\$71.75
Yellow, each. 134632075	\$17.95

134730000




Label Tape Dispensers

Multi-Roll. 134730000	\$68.80
Single-Roll. 134610000	\$152.00

Thermal Cycler



Versatile multi-format block holds (96) 0.2 mL tubes/strips, (39) 0.5 mL tubes, or (1) nonskirted 96-well PCR plate. In-situ adapter sold separately. Touchscreen control and height adjustable heated lid. 2-year warranty.



0.2 mL Tubes

96-Well Plate

0.5 mL Tubes

In-Situ Adapter

Thermal Cycler, multi-format block, Benchmark. **T5000-96** .. \$5,555.00
In-Situ Adapter, holds 4 slides. **T5000-A-IS**\$597.00

Trays



Trays

Cafeteria Style, polypropylene, 12 x 16¼", NSF, **USA**
Orange. **105CT1216O** \$5.95
Red. **105CT1216R**..... \$5.95
Black. **105CT1216B** \$5.95



Contain spills or transport lab products with this rimmed, durable, large white tray. Holds 5 gallons (18.9 L). *Not autoclavable.*

HDPE (high-density polyethylene), **USA**. **1053506**..... \$43.75
Dimensions: 26 x 18 x 3½" (66 x 45.7 x 8.9 cm)



Polypropylene, autoclavable

12 x 7¾ x 5⅞", Bel-Art, each. **162600000**.....\$32.70
15 x 14 x 5", Pk/6. **60081721**\$162.00
20 x 17 x 5", Pk/6. **60081722**\$171.00

Tubes & Vials

15 mL



Conical Bottom

50 mL

Tubes, Centrifuge, Conical bottom, sterile, white writing patch

RNase-/DNase-/Pyrogen-free

15 mL, PP (polypropylene) tube/HDPE cap, natural

Rack/25, 20 racks, Cs/500. **BP1501** \$124.00

Maximum RCF is 12,000 x g. Frosted writing patch.

15 mL, PS (polystyrene)

Rack/50, 10 racks, Cs/500. **6274** \$269.00

Maximum RCF is 3,000 x g.

50 mL, natural, polypropylene (PP) tube/HDPE cap

Bag/25, 20 bags, Cs/500. **BP5002**..... \$136.00

Maximum RCF is 12,000 x g. Frosted writing patch.

Bag/25, 20 bags, Cs/500. **6297** \$214.00

Maximum RCF is 20,000 x g

Rack/25, 20 racks, Cs/500. **BP5001** \$143.00

Maximum RCF is 12,000 x g. Frosted writing patch.

Self-Standing



Tubes, Centrifuge, Self-Standing (skirted)

30 mL

Nonsterile, PP tube/PP cap, Cs/500. **109220** \$199.00

Dimensions: 23(OD) x 90 mm

Sterile, PP tube/HDPE cap

Cs/500, 20 bags/25. **C572-1** \$158.00

Molded-in graduations from 5 to 30 mL in 5 mL increments. Dimensions: 25.3(OD) x 111 mm

Plastic culture tubes with dual-position caps.



Tubes, Culture,

12 x 75 mm (5 mL)	
Polystyrene, Sterile, dual-position caps	
Cs/500, 20 sterile bags/25. 110428	\$98.70
Polystyrene, Nonsterile, Cs/1000. 110410	\$64.00
Polypropylene, Sterile, dual-position caps	
Cs/500, 20 sterile bags/25. 110438	\$98.90
13 x 100 mm (10 mL)	
Glass (Type 1, Class A borosilicate), rimless	
Cs/1000. 1510	\$106.00
17 x 100 mm (16 mL)	
Polypropylene, Sterile, dual-position caps	
Cs/500, 20 sterile bags/25. 110178	\$159.00
16 x 125 mm (19 mL)	
Glass (Type 1, Class A borosilicate), rimless	
Cs/1000. 1515	\$175.00
16 x 150 mm (23 mL)	
 Glass (Type 1, Class A borosilicate), rimless	
Cs/1000. 1517	\$180.00



Culture Tube Caps, polypropylene

Up to 16 mm Tubes

Blue, Cs/1000. 118152B	\$110.00
Green, Cs/1000. 118152G	\$110.00
Red, Cs/1000. 118152R	\$110.00
Yellow, Cs/1000. 118152Y	\$110.00



15 mL

Black



RNase/DNase-free

Tubes, Microcentrifuge, Economy, polypropylene

0.5 mL, 0.1 mL graduations, natural, Pk/1000. HS4422	\$39.10
1.5 mL, 0.25 mL graduations, Heathrow	
Natural, max RCF is 16,400 x g	
Pk/500. HS4323	\$19.00
Pk/1000. HS4325	\$35.95
Black, max RCF is 18,000 x g, Pk/500. HS4323K	\$21.15
5 mL, natural, max RCF is 25,000 x g, Benchmark	
Nonsterile, 4 bags/50, Pk/200. C1005T5	\$68.45
Sterile, 10 bags/20, Pk/200. C1005T5ST	\$75.10



0.6 mL



1.5 mL



2.0 mL

Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable. Free of RNase/DNase.

Tubes, Microcentrifuge, Premium, PP, G-Tube™, Bio Plas, USA

Amber color protects light sensitive specimens.

0.6 mL, Pk/500	\$28.20
Amber. 4048 Lavender. 4047 Natural. 4040	
Mixed (blue, green, orange, pink, & yellow). 4049	\$34.00
1.5 mL, 0.5 mL graduations, frosted area on side, Pk/500	\$28.20
Amber. 4038 Blue. 4033 Lavender. 4037 Natural. 4030	
Mixed (blue, green, orange, pink, & yellow). 4039	\$34.00
2.0 mL, 0.5 mL graduations, frosted area on side	
Natural, Pk/500. 4050	\$30.30



Tubes, PCR, 0.2 mL thin-wall, polypropylene, Simport

Individual, flat caps, natural, Pk/1000. T325-2VN	\$51.20
--	---------



Vials, Cryogenic, round bottom, internal thread

O-Ring Seal, 4 mL, Cs/1000. T301-4	\$302.00
Washer Seal, 4 mL, Cs/1000. T311-4	\$302.00



Vinyl grip, USA

Vial Decapping Pliers

13 mm Crimp Seals. D-13	\$72.00
20 mm Crimp Seals. D-20	\$72.00

Tube Storage

Cryogenic Storage Boxes

SB2CC-64 Cryogenic
Tube Storage Box



Tube Storage Boxes, Cardboard, vented

2" Box Only, 5¼ x 5¼ x 2", each. SB2CC	\$1.95
2" Box w/dividers, white	
100-Cell (10 x 10), Pk/12. SB2CC-100	\$31.00
81-Cell (9 x 9), Pk/12. SB2CC-81	\$30.75
64-Cell (8 x 8), Pk/12. SB2CC-64	\$30.75
2" Mini box w/dividers, 3 x 3 x 2"	
25-Cell (5 x 5), Pk/12. MB2C-25	\$33.40
Holds 13 mm vials and 4.0-5.0 mL cryovials.	
3" Box Only, 5¼ x 5¼ x 3", each. SB3C	\$1.95
3" Box w/dividers	
100-Cell (10 x 10), Pk/12. SB3C-100	\$31.00
81-Cell (9 x 9), Pk/12. SB3C-81	\$30.75
64-Cell (8 x 8), Pk/12. SB3C-64	\$30.75
Dividers only, Fit 2 & 3" Boxes	
100-Cell (10 x 10), Pk/12. CD-100	\$12.90
Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.	
81-Cell (9 x 9), Pk/12. CD-81	\$12.60
Holds 13 mm vials and 4.0-5.0 mL cryovials.	
64-Cell (8 x 8), Pk/12. CD-64	\$12.60
Holds 14 mm vials.	
49-Cell (7 x 7), Pk/12. CD-49	\$12.90
Holds 14 mm vials.	
25-Cell (5 x 5), 25 mm cell, Pk/12. CD-25	\$12.60
16-Cell (4 x 4), 31 mm cell, Pk/12. CD-16	\$12.60
5" Box Only, 5¾ x 5¾ x 5", each. SB4C	\$3.60
5" Box w/dividers	
64-Cell (8 x 8), Pk/12. SB4C-64	\$47.40
36-Cell (6 x 6), 15 mL tubes, Pk/12. SB4C-36	\$47.40
16-Cell (4 x 4), 50 mL tubes, Pk/12. SB4C-16	\$47.40
Dividers only, Fit 5" Box	
64-Cell (8 x 8), Pk/12. CD4-64	\$12.60
36-Cell (6 x 6), Pk/12. CD4-36	\$12.60
16-Cell (4 x 4), Pk/12. CD4-16	\$12.60



Holds 1.0-2.0 mL tubes and cryovials.
Features a locking hinged lid, and base
with molded alphanumeric coordinates.
Dimensions: 12.8 x 12.8 x 5.4 cm.

Tube Storage Boxes, Polypropylene


81-Well (9x9), 2 blue/2 yellow boxes, Pk/4. 60066501	\$36.50
---	---------

Tube Racks




Tube Racks, Grip, polypropylene, No-Wire™, Bel-Art

18-20 mm Tubes, 40-place, each. 187490002	\$36.95
--	---------




Tube Racks, Interlocking, polypropylene, Heathrow

4-Way, mixed colors, Pk/5. HS29022G	\$47.40
Holds (4) 50 mL tubes, (12) 15 mL tubes, (32) 1.5 mL tubes, or (32) 0.5 mL tubes.	
4-Way Rotatable, Rota-Rack™, four connected modules	
Large Tubes, 22.4 x 11 x 11 cm, each. HS29040A	\$49.30
Each module holds (2) 50 mL tubes, (8) 15 mL tubes, (6) 20 mm tubes, or (10) 10-13 mm tubes. Mixed colors.	



Tube Racks, Mega Rack, polypropylene, fully-assembled

10-13 mm Tubes	
216-Place, 33.7 x 26.7 x 9.1 cm(H)	
Blue, each. 34524322A	\$51.60
White, each. 34524322E	\$51.60
432-Place, 67.4 x 26.7 x 9.1 cm(H)	
Blue, each. 345120113	\$90.25
White, each. 345120112	\$90.25



Tube Racks, Microtube, polypropylene, Heathrow

80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm	
Holds 1.5-2.0 mL tubes in conical wells.	
Blue, Pk/5. HS29025B	\$28.50
Pink, Pk/5. HS29025D	\$28.50
Mixed, Pk/5. HS29025G	\$28.50


Available in blue,
green, orange,
purple & yellow.



Tube Rack, PCR, polypropylene, mixed colors, Heathrow
Reversible, Pk/5. **HS2344A** \$62.30
Holds (168) 0.2, (40) 0.5 and (12) 1.5 mL tubes.



Half-Size



Full-Size

Tube Racks, Teardrop, OneRack™, POM, Heathrow
Half-Size
16 mm tubes, 36-place, 12.7 x 12.7 x 7 cm(H)
Blue, Pk/8. **HS27502B** \$124.00
Full-Size
13 mm tubes, 72-place, 20.2 x 10.4 x 5.9 cm(H)
Blue, Pk/8. **HS27511B** \$164.00
Magenta, Pk/8. **HS27511D** \$164.00
16 mm tubes, 72-place, 25 x 12.7 x 7 cm(H)
Blue, Pk/8. **HS27512B** \$164.00
Magenta, Pk/8. **HS27512D** \$164.00



Polypropylene cube and lid has
two snap-in, gridded inserts to
hold either (16) 15 mL tubes or
(9) 50 mL tubes.

Tube Racks, Transformer Cube, Pk/5. **HS23271A** \$64.50



HS232461

Tube Racks, Wire, epoxy-coated
10-13 mm Tubes
108-place (9 x 12), each. **345120759** \$34.90
13-16 mm Tubes
108-Place (9 x 12), each. **345120763** \$34.90
15-16 mm Tubes
180-Place (9 x 20), each. **187860770** \$63.15
Dimensions: 38.6 x 18.3 x 6.4 cm(H)
25-30 mm Tubes (50 mL)
8-Place (2 x 4), each. **HS232461** **NEW** \$36.00

Weighing Dishes



Hexagonal



Standard



Pour

Weigh Boats, polystyrene, Heathrow

Standard

★ Small, 4.6 x 4.6 x 0.8 cm(D), Bx/500. **HS1420A** \$25.30
★ Medium, 8.5 x 8.5 x 2.4 cm(D), Bx/500. **HS1420B** \$41.15
★ Large, 14 x 14 x 2.2 cm(D), Bx/500. **HS1420C** \$67.65

Hexagonal, anti-static

Small, 3.8 cm top dia., Bx/500. **HS14251A** \$18.60
Medium, 6.4 cm top dia., Bx/500. **HS14251B** \$32.10
Large, 12 x 14 cm top OD, Bx/500. **HS14251C** \$84.20
(25433-104)
X-Large, 14 x 16 cm top OD, Bx/500. **HS14251D** \$107.00
(25433-106)

Pour (Canoes)

Small, 5.8 x 4.3 x 1.3 cm(D), Bx/500. **HS1419A** \$42.15
Medium, 13.2 x 8.3 x 2.6 cm(D), Bx/500. **HS1419B** \$82.70
Large, 18.3 x 10.8 x 2.6 cm(D), Bx/500. **HS1419C** \$120.00

Wipes



Extra low-lint, nonabrasive,
highly absorbent, 1-ply wipes
are perfect for microscope
maintenance.

Kimwipes®
4½ x 8½" Wipes. **34155** bx/280 cs/60 bx
\$4.00 214.00



WypAll® L40

Lab Wipes, ¼ fold, low lint
WypAll® L40, limited use. **43005701** pk/56 cs/18 pk
Measures 12 x 12.5" \$7.45 117.00