



EDUCATION & NOVELTY ITEMS

Books



Adult

- How to Use a Light Microscope* by Frank M. Carryl \$12.00
Great, straightforward, undecorated primer.
- Optimizing Light Microscopy* by Barbara Foster \$37.00
- Photography Through the Microscope*
by Eastman Kodak \$24.95

Youth

- The Ultimate Guide to Your Microscope* \$9.95
by S. Levine & L. Johnstone



Mounted on a telescoping aluminum handle that extends from 18" to 32", this conical-shaped, 36" long nylon net has a 14" opening.

- Insect Net. 12571402** \$18.45
Dimensions: 17.7 x 15.7 x 12"


Sizes:

#5 (4000 μ), #10 (2000 μ),
#35 (500 μ), #60 (250 μ),
#120 (125 μ), #230 (63 μ)
Sieves measure 6.5 x 2.5" (H).



- Sieve Set, 6-pc, plastic frame & wire screen. 12530706** \$118.00
Includes bottom pan, lid, and teacher's guide.

Demonstrate Koehler illumination with this unique tool. Visualize the function of the condenser by tracking light rays as they travel through this "smoked" acrylic block.



- Smoked Acrylic Block, 40 x 35 x 20 mm. 85502** \$29.00
Refractive index is 1.48 at 546 nm.

Novelty Items

- Micro Bible. 89585** \$15.00
View with a 10X objective all 1,245 pages of the King James version of the Holy Bible on this approximately 1" square microfilm.



- Ruler, 1 Decimeter, stainless steel, Germany. 605-060** \$6.00
Ruler is 1 dm long, 1 cm wide, and 1 mm thick. Its volume is 1 cm³ and its weight in grams is the specific gravity. One edge of one side is graduated in centimeters with 1 mm divisions.

Reduce waste & save money with this heavy-duty wringer by completely extracting contents from even the most rigid tube. Features 3" aluminum rollers with steel handle.



- Tube-Wringer®, USA. 295401** \$38.95